



Europäisches
Patentamt

European
Patent Office

Office européen
des brevets

REC'D 23 JUL 1997
WIPO PCT

Bescheinigung

Certificate

Attestation

Die angehefteten Unterla-
gen stimmen mit der
ursprünglich eingereichten
Fassung der auf dem näch-
sten Blatt bezeichneten
europäischen Patentanmel-
dung überein.

The attached documents
are exact copies of the
European patent application
described on the following
page, as originally filed.

Les documents fixés à
cette attestation sont
conformes à la version
initialement déposée de
la demande de brevet
européen spécifiée à la
page suivante.

Patentanmeldung Nr. Patent application No. Demande de brevet n°

96810421.6

PRIORITY DOCUMENT

Der Präsident des Europäischen Patentamts;
Im Auftrag

For the President of the European Patent Office

Le Président de l'Office européen des brevets
p.o.

J-N. Haon
tel: (070)3403149

DEN HAAG, DEN
THE HAGUE,
LA HAYE, LE

06/97



Europäisches
Patentamt

European
Patent Office

Office européen
des brevets

**Blatt 2 der Bescheinigung
Sheet 2 of the certificate
Page 2 de l'attestation**

Anmeldung Nr.:
Application no.: 96810421.6
Demande n°:

Anmeldetag:
Date of filing: 24/06/96
Date de dépôt:

Anmelder:
Applicant(s):
Demandeur(s):
ROTKREUZSTIFTUNG ZENTRALLABORATORIUM BLUTSPENDEDIENST SRK
CH-3000 Bern 22
SWITZERLAND

Bezeichnung der Erfindung:
Title of the invention:
Titre de l'invention:

Polypeptides capable of forming antigen binding structures with specificity for the Rhesus D antigens, the DNA encoding them and the process for their preparation and use

In Anspruch genommene Priorität(en) / Priority(ies) claimed / Priorité(s) revendiquée(s)

Staat:
State:
Pays:

Tag:
Date:
Date:

Aktenzeichen:
File no.
Numéro de dépôt:

Internationale Patentklassifikation:
International Patent classification:
Classification internationale des brevets:

C12N15/13, C12N15/62, C07K16/34, A61K39/395, G01N33/80

Am Anmeldetag benannte Vertragsstaaten.
Contracting states designated at date of filing: AT/BE/CH/DE/DK/ES/FI/FR/GB/GR/IE/IT/LI/LU/MC/NL/PT/SE
Etats contractants désignés lors du dépôt:

Bemerkungen:
Remarks:
Remarques:

Polypeptides capable of forming antigen binding structures with specificity for the Rhesus D antigens, the DNA encoding them and the process for their preparation and use

5 This invention relates to polypeptides forming antigen binding structures with specificity for Rhesus D antigens and especially to Fab molecules with specificity for the Rhesus D antigen. The invention also relates to their application to provide pharmacological and diagnostic compositions. The above Fab fragments when genetically engineered to be part of complete antibodies are useful for the prophylaxis of hemolytic
10 disease of the newborn (HDN). This invention provides the novel DNA and amino acid sequences of the above polypeptides.

Thus, the antibodies can be used for the protection of Rhesus negative women before or immediately after the birth of a Rhesus positive child to prevent HDN in subsequent pregnancies.

15 The invention also includes the application of the said Fab molecules either alone or in combination with Fc constant regions as complete antibodies for the purposes of treating other illnesses which might benefit from anti-Rhesus D immunoglobulin e.g. treatment of idiopathic thrombocytopenic purpura (ITP).

20 In addition anti-Rhesus D immunoglobulin can be used after mistransfusions of Rhesus positive blood to Rhesus negative recipients in order to prevent sensitization to the Rhesus D antigen. Further the invention relates to the application of these Fab fragments and antibodies as diagnostic reagents.

25 HDN is the general designation for hemolytic anemia of fetuses and newborn babies caused by antibodies of the mother. These antibodies are directed against antigens on the surface of the fetal erythrocytes. These antigens can belong to the Rhesus, ABO or other blood group systems.

The Rhesus blood group system includes 5 major antigens: D, C, c, E and e (Issitt, P.D., Med. Lab. Sci. 45:395, 1988). The D antigen is the most important of these antigens as it is highly immunogenic eliciting anti-Rhesus D antibodies during Rhesus incompatible pregnancies and following transfusion of Rhesus incompatible blood. The D antigen is found in approximately 85% of Caucasians in Europe and those individuals are said to be Rhesus positive. Individuals lacking the D antigen are called Rhesus negative. The expression of the D antigen can vary due to either low antigen density, hereafter known as weak D or D^u, or due to partial antigenicity, hereafter known as partial D antigens.

The Rhesus D antigen, a membrane protein of the erythrocyte, has recently been cloned and its primary structure described (Le Van Kim, C., et al., PNAS 89:10925, 1992). Modeling studies suggest that the Rhesus D antigen has 12 transmembrane domains with only very short connecting regions extending outside the cell membrane or protruding into the cytoplasm.

The partial D phenotypes were first identified in people who carried D antigen on their red cells but who had an alloanti-D in their sera (Rose, R. R. and Sanger, R., Blood groups in man, Blackwell Scientific, Oxford, U.K. 1975; Tippett, P. et al., Vox Sanguinis. 70:123, 1996). This can be explained by regarding the D antigen as a mosaic structure with at least 9 different epitopes (epD1 to epD9). Thus in some D variant people the red cells lack part of this mosaic and antibodies are made to the missing D epitopes. Rhesus positive individuals that make antibodies against partial D antigens have been classified into six main different categories (D^{II} to D^{VII}) each having a different abnormality in the D antigen. More recently it has been shown that these D categories gave different patterns of reaction when tested against panels of human monoclonal anti-D antibodies (Tippett, P., et al., Vox Sanguinis. 70:123, 1996). The different reaction patterns identified the 9 epitopes and so define the different partial D categories. The number of epitopes present on the D antigen varies from one partial D category to another with the D^{VI} category expressing the least, epD3, 4 and 9. The D^{VI} category is clinically important as a D^{VI} woman can be immunized strongly enough to cause hemolytic disease of the newborn.

The prophylactic efficacy of anti-RhD IgG for prevention of hemolytic disease of the newborn is well established and has been in routine use for many years. As a result this severe disease has become a rarity. Nevertheless the underlying cause of the disease, i.e. RhD incompatibility
5 between a RhD negative mother carrying a RhD positive child still remains and thus requires a continual supply of therapeutic anti-RhD IgG.

In recent years the assurance of a continual supply of anti-RhD IgG has become an increasing problem. The pool of available hyperimmune serum from alloimmunized multiparous Rhesus negative women has
10 drastically decreased due to the success of prophylactic anti-RhD. Thus the current methods of production require repeated immunization of an increasingly reluctant pool of donors for the production of high titer antiserum (Selinger, M., Br. J. Obstet. Gynaecol. 98:509, 1991). There are also associated risk factors and technical problems such as the use of Rhesus
15 positive red blood cells for repeated immunization carrying the risk of transmission of viral diseases like hepatitis B, AIDS and other as yet unknown viruses (Hughes-Jones, N.C., Br. J. Haematol. 70:263, 1988). Therefore an alternative method for production of anti-RhD antibodies is required.

In the past few years various alternative sources of hyperimmune
20 serum have been tried but all are associated with disadvantages. Epstein Barr Virus (EBV) transformation of lymphocytes creating B lymphoblastoid cell lines that secrete specific antibody including against the Rhesus D antigen (Crawford et al., Lancet. 386:Feb.19th, 1983) are unstable and require extensive cloning. Also due to the low transformation efficiencies (1-3% of B
25 cells) only a restricted range of antibody specificities can be obtained from the potential repertoire. Additionally it seems that mice do not respond to the Rhesus D antigen and thus no murine monoclonal antibodies are available which could be used for producing chimaeric or humanised antibodies. Until recently the only other alternative was production of human antibodies by the
30 hybridoma technique which was also restricted by the lack of a suitable human myeloma cell fusion partner (Kozbor, D. and Roder, J.C., Immunol. Today. 4:72, 1983).

It is thus the object of the present invention to provide Fab fragments having a reactivity against the Rhesus D antigen as well as complete antibodies comprising the Fab fragments which are free from the above mentioned drawbacks.

5 In the last few years the technique of repertoire cloning and the construction of phage display libraries has opened up new possibilities to produce human antibodies of defined specificity (Williamson, R.A. et al., PNAS 90:4141, 1993). These methods were thus applied to the preparation of polypeptides capable of forming antigen binding structures with specificity for
10 Rhesus D antigens, especially of Fab fragments having an activity against Rhesus D and partial D antigens.

The generation of human antibodies by repertoire cloning as described in recent years (Barbas III, C.F. and Lerner, R.A., Companion Methods Enzymol. 2:119, 1991) is based on isolating mRNA from peripheral
15 B cells. This method offers the tools to isolate natural antibodies, autoantibodies or antibodies generated during the course of an immune response (Zebedee, S.L., et al., PNAS 89:3175, 1992; Vogel, M. et al., Eur.J. Immunol. 24:1200, 1994). This method relies on constructing a recombinant antibody library from a particular donor starting from the mRNA coding for
20 immunoglobulin (Ig) molecules. As only the peripheral blood lymphocytes (PBL) can be isolated from a donor the chances of finding specific antibody producing B cells in the periphery are increased if an individual is boosted with the desired antigen shortly before harvesting the PBL (Persson, M.A.A., et al., PNAS 88:2432, 1991). The total RNA is then isolated and the mRNA of
25 the Ig repertoire can be cloned using Ig specific primers in the polymerase chain reaction (PCR) followed by the co-expression of heavy and light chains of the Ig molecule on the surface of a filamentous phage particle thereby forming an „organism“ that in analogy to a B cell can bind to an antigen. In the literature this method is also known as the combinatorial approach as it
30 allows the independent combining of heavy and light chains to form a functional Fab antibody fragment attached to one of the tail proteins, called pIII, of a filamentous phage. Phages carrying the Fab molecules (hereafter known as Phab particles) are selected for the desired antigen specificity, by a

process known as bio-panning. The antigen can be applied to a solid support, specific Phab bind to the antigen whilst non specific Phab are washed away and finally the specific Phab are eluted from the solid support. The specific Phab are then amplified in bacteria, allowed to re-bind to the antigen on the solid support and the whole process of bio-panning is repeated.

The successive rounds of panning and amplification of selected Phab in bacteria result in an enrichment of specific Phab that can be seen from a rise in titer of colony forming units (cfu) plated out after each round of panning. Our previous experience and published data indicate that specific phage can usually be detected after 4 to 6 panning rounds (Vogel, M. et al., Eur.J. Immunol. 24:1200, 1994). In the above cited related art there is, however, no hint that the indicated steps can be used for a successful preparation of Fab fragments of anti-Rh D antibodies.

In the appended figures 1a to 17b; DNA sequences coding for variable regions (V regions) of anti Rh D Fab fragments and the corresponding polypeptide sequences are disclosed.

Fig. 18 shows the pComb3 expression system used according to the present invention.

Figs. 19 and 20 show the separate preparation of genes of the heavy and light chains of the complete antibody according to the description in example 6.

Subjects of the present invention are polypeptides capable of forming antigen binding structures with specificity for Rhesus D antigens according to the definition of claim 1. The table in claim 1 refers to the appended figures. The identification number for each sequence is given. The locations of the Rhesus D specific CDR1 (complementarity determining region 1), CDR2 and CDR3 regions are indicated in the figures and according to base pair number in the table of claim 1. Preferred polypeptides according to the invention are anti-Rhesus D antibodies which include the variable regions of the heavy and light chains according to the sequences given in

Figs. 1a -17b. The Figs. 1a, 2a, ... 17a are related to the variable regions of the heavy chain and the Figs. 1b, 2b, ... 17b are related to the variable regions of the light chain.

Further subjects of the present invention are the DNA sequences
5 coding for antigen binding polypeptides according to the definition of claim 5. Preferred DNA sequences are those coding for variable regions of Fab fragments of anti-Rh D antibodies according to the Figs. 1a -17b. The Figs. 1a, 2a, ... 17a are related to the heavy chain and the Figs. 1b, 2b, ... 17b are related to the light chain.

10 A further subject of the present invention is the process for preparing recombinant Fab polypeptides according to the definition in claim 9.

Further subjects of the present invention are anti-Rh D antibodies according to the definition of claim 10, preferably anti-Rh D immunoglobulin molecules comprising the heavy and light chain variable regions according to
15 the Figs. 1a to 17b combined with known heavy and light chain constant regions.

Further subjects of the present invention are pharmaceutical and diagnostic compositions comprising polypeptides, anti-Rh D antibodies or Fab fragments according to the invention.

20 The total re-amplified Phab population obtained after each panning can be tested for specificity using various methods such as ELISA and immunodot assays. It is also defined by the nature of the antigen e.g. anti-Rhesus D Phabs are detected by indirect haemagglutination using a rabbit anti-phage antibody or equivalent Coombs reagent as the cross linking
25 antibody. Once a total Phab population has been identified as positive for the desired antigen, individual Phab clones are isolated and the DNA coding for the desired Fab molecules is sequenced. Individual Fab can then be produced by use of the pComb3 expression system which is illustrated in Fig. 18. In this system the gIII gene, coding for the tail protein pIII, is cut out from
30 the phagemid vector pComb3. This allows production of soluble Fab in the

bacterial periplasm. Such individual Fab fragments can then be tested for antigen specificity.

The phage display approach has also been used as a means of rescuing monoclonal antibodies from unstable hybridoma cell lines. This has
5 been reported for anti-Rhesus D antibodies (Siegel, D.L. and Silberstein, L.E., Blood. 83:2334, 1994; Dziegiel, M. et al., J. Immunol. Methods. 182:7, 1995). A phage display library constructed from non-immunized donors has also been used to select Fv fragments (i.e. variable regions of heavy and light chains, V_H and V_L) specific for human blood group antigens which included
10 one Fv fragment reacting against the Rhesus D antigen (Marks, J.D. et al., Biotechnology. 11:1145, 1993).

Important considerations when constructing combinatorial libraries are the source of cells used for RNA extraction and the nature of the antigen used for panning. Therefore, this invention uses a hyperimmune donor who
15 was boosted i.v. with Rhesus D⁺ red blood cells (rbc). The PBL of the donor were harvested at +5 and +18 days after the i.v. boost and were used to construct 2 combinatorial libraries hereafter known as library D1 (LD1) and library D2 (LD2) respectively. Double immunofluorescence analysis of the harvested PBL, using the markers CD20 and CD38 for pan B cells and
20 lymphoblastoid cells respectively, showed a higher than normal percentage of lymphoblastoid B cells, of plasma cell morphology. The high number of plasma cells found in the peripheral blood is most unusual as normally there are less than 1% in the periphery and probably indicates that the donor had a high percentage of circulating B cells with specificity for the Rhesus D
25 antigen.

After construction of the library, the selection of Phabs specific for the Rhesus D antigen was achieved by bio-panning on fresh whole rbc of phenotype R1R1 (CDe/CDe) i.e. the reference cells used for Rhesus D
typing. This was necessary since the Rhesus D antigen, an integral
30 membrane protein of 417 amino acids (Le Van Kim, C. et al, PNAS 89:10925, 1992), loses its immunogenicity during purification (Paradis, G. et al, J. Immunol. 137:240, 1986) and therefore a chemically purified D antigen

cannot be bound to a solid phase for selection of immunoreactive Phabs as for other antigen specificities previously selected in this system (Vogel, M. et al., Eur.J. Immunol. 24:1200, 1994). Modelling studies have suggested that only very short connecting regions of the Rhesus D antigen extend outside
 5 the cell membrane or protrude into the cytoplasm (Chérif-Zahar, B. et al, PNAS 87:6243, 1990). Thus the parts of the RhD antigen visible to antibodies are relatively restricted and may be under conformational constraint. This aspect of the Rhesus D antigen becomes even more important when considering selection of Phabs with reactivity against the partial D
 10 phenotypes which essentially lack certain defined epitopes of the D membrane protein (Mouro, I. et al, Blood. 83:1129, 1994).

Furthermore, since whole rbc do not only express the D antigen, a series of negative absorptions had to be performed on Rhesus D negative rbc in order to absorb out those Phabs reacting with the other antigenic proteins
 15 found on the rbc.

This panning procedure performed on Phabs coming from both LD1 and LD2 libraries resulted in the isolation of 7 different Fab producing clones from library LD1 and 9 different Fab producing clones from library LD2.

20 The nomenclature and the figures where the sequences are listed are given in table 1.

Table 1

| LIBRARY LD1 Clone No. | V _H - Sequence Figure | V _L - Sequence Figure | LIBRARY LD2 Clone No. | V _H - Sequence Figure | V _L - Sequence Figure |
|-----------------------------|--|--|-----------------------------|--|--|
| LD1-28 | 1a | 1b | LD2-1 | 8a | 8b |
| LD1-40 | 2a | 2b | LD2-4 | 9a | 9b |
| LD1-52 | 3a | 3b | LD2-5 | 10a | 10b |
| LD1-84 | 4a | 4b | LD2-10 | 11a | 11b |
| LD1-98 | 5a | 5b | LD2-11 | 12a | 12b |
| LD1-110 | 6a | 6b | LD2-14 | 13a | 13b |
| LD1-117 | 7a | 7b | LD2-17 | 14a | 14b |
| | | | LD2-18 | 15a | 15b |
| | | | LD2-20 | 16a | 16b |

The above Fab clones show exclusive reactivity against the Rhesus D antigen, 3 of 5 D^u rbc tested and agglutinating reactivity against the Partial D phenotypes as follows: Rh33, DIII, DIVa, DIVb, DVII.

5 However using the above mentioned R1R1 rbc for panning of the Phabs no clones were isolated which reacted against the Partial DVI phenotype. As the serum of the original hyperimmune donor tested at the time of construction of the recombinant library, was known to react against the DVI phenotype the recombinant library should also contain the anti-DVI
10 specificity.

In order to select for the DVI reactivity the panning conditions were changed in that different cells were used. A special donor whose rbc had been typed and were known to express the Partial DVI phenotype was used as the source of cells for re-panning the LD1 and LD2 librairies. This second
15 series of pannings was essentially performed in the same way as the first series except for the substitution of DVI rbc for R1R1 rbc and the addition of bromelase treatment to the DVI rbc. The DVI phenotype expresses the least number of Rhesus D epitopes and it is therefore difficult to make antibodies against it. It has been reported that only 15% of unselected polyclonal anti-D
20 and 35% of selected anti-D made by Rhesus D negative subjects reacted with DVI+ cells (Mouro, I. et al, Blood. 83:1129, 1994). Bromelase treatment which removes N- acetylneuraminic acid (sialic acid) from the rbc membrane, was performed in order to render the Rhesus DVI epitopes more accessible during the panning with the pre-absorbed Phabs.

25 This second series of pannings on the LD1 and LD2 librairies resulted in 1 Fab producing clone reacting with DVI+ rbc.

The nomenclature is given below:

| LIBRARY LD1 | V _H -Sequence figure | V _L -Sequence figure |
|--------------------|------------------------------------|------------------------------------|
| Clone No: LD1-6-17 | 17a | 17b |

Thus a total of 17 different anti-Rhesus D Fab clones have been isolated. The DNA from these clones has been isolated and sequenced using Fluorescent Cycle Sequencing on an ABI 373A Sequencing System. The
5 nucleotide and corresponding amino acid sequences of the said Fab clones form the basis of this invention.

The DNA sequences obtained and Fab fragments are useful for the preparation of complete antibodies having an activity against the Rhesus D antigen.

10 The examples which follow explain the invention in detail, without any restriction of the scope of the invention.

Example 1 describes the construction of 2 combinatorial libraries; especially the aforementioned LD1 and LD2 libraries.

15 Example 2 describes a series of pannings using R1R1 rbc on the said LD1 and LD2 libraries in detail.

Example 3 describes a series of pannings using bromelase treated DVI+ rbc on the said LD1 and LD2 libraries.

20 Example 4 describes an indirect haemagglutination assay using a rabbit anti-phage antibody, as an equivalent Coombs reagent, to monitor the enrichment and specificity of Rhesus D specific Phabs after panning.

Example 5 describes the preparation and purification of Fab antibody fragments for application as diagnostic reagents.

Example 6 describes the preparation of complete anti-Rhesus D immunoglobulins using the sequences of the present invention.

Example 1

Construction of the recombinant LD1 and LD2 libraries

a) Source of the lymphocytes

A male adult who was a member of the volunteer pool of
 5 hyperimmune Rhesus D donors was given an i.v. boost of 2 ml of packed rbc
 from a known male donor of blood group O RhD⁺. The PBL were harvested at
 +5 and +18 days after the boost and the mononuclear cells (MNC) isolated by
 Ficoll gradient centrifugation (Lymphoprep, Pharmacia, Milwaukee, WI). The
 +5 day MNC were used directly for RNA preparation using a phenol-
 10 chloroform guanidinium isothiocyanate procedure (Chomczynski, P. and
 Sacchi, N., Anal. Biochem. 162:156, 1987). The +18 day MNC were first
 cultured for 3 days in RPMI-1640 medium (Seromed, Basel) containing 10³
 U/ml of IL-2 (Sandoz Research Center, Vienna, Austria) and 10 µg/ml of
 pokeweed mitogen (PWM; Sigma L9379, Buchs, Switzerland) before
 15 extracting RNA.

Immunofluorescence analysis of donor lymphocytes +5 days after rbc i.v. boost

| Cell Surface Antigen | % Positive cells | Cell Surface Antigen | % Positive cells |
|----------------------|------------------|----------------------|------------------|
| CD20 | 15 | CD8 | 12 |
| CD38 | 20 | CD25 | 7.6 |
| CD20/38 | 15 | CD57 | 12.5 |
| CD3 | 47 | CD14 | 6 |
| CD4 | 34 | HLA-DR | 18 |

b) Construction of Library

20 Two separate libraries were constructed called LD1 and LD2 (as
 detailed above) corresponding to the cells harvested at +5 days and +18 days
 (finally +21 days including the +3 days PWM stimulation) after the i.v. boost

respectively. Total RNA was then prepared from these cells using a phenol-chloroform guanidinium isothiocyanate method. From this RNA, 10 µg were used to make cDNA using an oligo(dT) primer (400 ng) and reverse transcribed with M-MuLV reverse transcriptase according to the conditions specified by the supplier (Boehringer Mannheim Germany). PCR amplification was performed as described in Vogel, M. et al., E.J. of Immunol. 24:1200, 1994. Briefly, 100 µl PCR reaction contained Perkin-Elmer buffer with 10 mM MgCl₂, 5 µl cDNA, 150 ng of each appropriate 5' and 3' primer, all four dNTP at 200 µM each and 2 U/ml Taq Polymerase (Perkin Elmer, NJ). The PCR amplification of the heavy and light chains of the Fab molecule was performed separately with a set of primers from Stratacyte (details given below). For the heavy chain six upstream primers were used that hybridize to each of the six families of the V_H genes whereas one kappa and one lambda chain primer were used for the light chain. The downstream primers were designed to match the hinge region of the constant domains γ1 and γ3 for the heavy chain. For the light chain the downstream primers were matched to the 3' end of kappa and lambda constant domains. The heavy and light chain PCR products were pooled separately, gel purified and cut with Xho1/Spe1 and Sac1/ Xba1 restriction enzymes (Boehringer Mannheim), respectively. After digestion the PCR products were extracted once with phenol : chloroform : isoamylalcohol and purified by gel excision. The insertion of the Xho1/Spe1 digested Fd fragment and subsequent ligation of the Sac1/Xba1 digested light chain into the pComb3 vector, the transformation into XL1-Blue cells, and the production of phages were performed as described by (Barbas III, C.F. and Lerner, R.A., Companion Methods Enzymol. 2:119, 1991).

After transformation of the XL1-Blue E.coli cells samples were withdrawn and titrated on plates to determine the library size. These results indicated expression libraries of 7.5×10^6 and 7.7×10^6 cfu (colony forming units) for LD1 and LD2 respectively.

30. c) PCR Primers

VHI 5'-CAC TCC CAG GTG CAG CTG CTC GAG TCT GG-3'
 VHII 5'-GTG CTG TCC CAG GTC AAC TTA CTC GAG TCT GG-3'

VHIII 5'-GTC CAG GTG GAG GTG CAG CTG CTC GAG TCT GG-3'
 VHIV 5'-GTC CTG TCC CAG GTG CAG CTG CTC GAG TCG GG-3'
 VHV 5'-GTC TGT GCC GAG GTG CAG CTG CTC GAG TCT GG-3'
 VHVI 5'-GTC CTG TCA CAG GTA CAG CTG CTC GAG TCA GG-3'
 5 CHI(gI) 5'-AGC ATC ACT AGT ACA AGA TTT GGG CTC-3'
 VL(k) 5'-GT GCG AGA TGT GAG CTC GTG ATG ACC CAG TCT CCA-3'
 CL(k) 5'-T CCT TCT AGA TTA CTA ACA CTC TCC CCT GTT GAA GCT
 CTT TGT GAC GGG CGA ACT C-3'
 VL(l) 5'C TGC ACA GGG TCC TGG GCC GAG CTC GTG GTG ACT CA-3'
 10 CL(l) 5'G CAT TCT AGA CTA TTA TGA ACA TTC TGT AGG GGC-3'

d) Vectors and bacterial strains

The pComb3 vector used for cloning of the Fd and the light chain was obtained from the Scripps Research Institute La Jolla, CA; (Barbas III, C.F. and Lerner, R.A., Companion Methods Enzymol. 2:119, 1991). The
 15 *Escherichia coli* strain XL1-Blue used for transformation of the pComb3 vector and the VCSM13 helper phage were purchased from Stratacyte (La Jolla, CA).

Example 2

Selection of Rhesus D Phabs from LD1 and LD2 librairies.on R1R1 rbc

a) Absorption and Bio-Panning

A series of three negative absorptions on rbc group O Rh negative were performed for each panning round before positive selection on rbc group O Rh positive (R1R1). Fresh rbc were collected in ACD (acid citrate dextrose) anticoagulant and washed 3 times in 0.9% NaCl. The rbc were
 25 counted in Hayems solution and adjusted to 40×10^6 /ml. Absorption : 1 ml of phage preparation in PBS/3%BSA was added to rbc group O Rh negative pellet (16×10^6 rbc) in 12 ml tubes (Greiner 187261, Reinach, Switzerland) and incubated at RT for 30 min. with careful shaking. All tubes were pre-blocked in PBS/3% BSA for a minimum of 1hr at RT. The rbc were pelleted by
 30 centrifuging for 5 min. 300 x g at 4°C. The resulting phage supernatant was

carefully harvested and the process repeated twice more. After the final absorption the phage supernatant was added to the rbc group O Rh positive pellet (16×10^6 rbc) and again incubated at RT for 30 min. with gentle shaking. Then the rbc were washed at least 5 times in 10 ml ice cold PBS, centrifuged 5 min. $300 \times g$ at 4°C , followed by elution with $200 \mu\text{l}$ of 76 mM citric acid pH 2.8 for 6 min. at R.T. and neutralisation with $200 \mu\text{l}$ 1M Tris. The rbc were centrifuged $300 \times g$, 5 min. at 4°C and the resulting supernatant containing the eluted phages was carefully removed and stored with carrier protein (0.3% BSA) at 4°C ready for re-amplification.

10 **Selection of Rhesus D+ Phabs from the LD1 and LD2 librairies on R1R1 rbc**

| Panning Round No.a) | No. of eluted Rhesus D specific phages | |
|------------------------|--|-------------------|
| | Library D1 cfu | Library D2 cfu |
| 1 | 8×10^6 | 4.6×10^7 |
| 2 | 6×10^7 | 1.4×10^7 |
| 3 | 1×10^8 | 7.9×10^7 |
| 4 | 3×10^8 | 1.3×10^8 |
| 5 | 3×10^8 | 1×10^8 |
| 6 | nd | 2.8×10^8 |

15 a) For each round 10^{12} Phabs were incubated in tubes with rbc Group O Rhesus negative (absorption phase) followed by elution from rbc Group O Rhesus positive (R1R1)

nd = not done

cfu = colony forming units

Example 3

Selection of Rhesus D Phabs from LD1 and LD2 librairies on DVI+ rbc

a) Absorption on rbc group O Rh negative

A series of three negative absorptions on rbc group O Rh negative
 5 were performed for each panning round before positive selection on rbc
 group O Rh DVI positive. Fresh rbc were collected in ACD anticoagulant and
 washed 3 times in 0.9% NaCl. The rbc were counted in Hayems solution and
 adjusted to 40×10^6 /ml. Absorption : 1 ml of phage preparation in PBS/3%BSA
 was added to rbc group O Rh negative pellet (16×10^6) in 12 ml tubes (Greiner
 10 187261, Reinach, Switzerland) and incubated at RT for 30 min. with careful
 shaking. All tubes were pre-blocked in PBS/3% BSA for a minimum of 1hr at
 RT. The rbc were pelleted by centrifuging for 5 min. 300 x g at 4°C. The
 resulting phage supernatant was carefully harvested and the process
 repeated twice more. Treatment of the rbc group O Rhesus DVI+ with the
 15 enzyme bromelase was performed at this stage in order to enhance
 accessibility of the antigens.

b) Treatment of rbc Rhesus DVI+ with bromelase

Bromelase 30 (Baxter, Düringen, Switzerland) was used to treat
 rbc Rhesus DVI+ in the same proportions as used in a routine
 20 haemagglutination assay, i.e. 10 µl bromelase per 2×10^6 rbc. Thus 80 µl of
 bromelase was added to 16×10^6 DVI+ rbc and incubated at 37°C for 30 min.
 followed by washing 3 times in 0.9% NaCl, re-counting in Hayems solution
 and adjusting to 40×10^6 /ml in PBS/3% BSA ready for Phab panning.

c) Bio-Panning on bromelase treated Rhesus DVI+ rbc

25 After the final absorption the phage supernatant was added to the
 enzyme treated rbc group O Rh DVI+ pellet (16×10^6) and again incubated at
 RT for 30 min. with gentle shaking. Then the rbc were washed at least 5 times
 in 10 ml ice cold PBS, centrifuged 5 min. 300 x g at 4°C, followed by elution
 with 200 µl of 76 mM citric acid pH 2.8 for 6 min. at R.T. and neutralisation
 30 with 200 µl 1M Tris. The rbc were centrifuged 300 x g, 5 min. at 4°C and the

resulting supernatant containing the eluted phages was carefully removed and stored with carrier protein (0.3% BSA) at 4°C ready for re-amplification.

Selection of Rhesus D Phabs from LD1 and LD2 librairies on Rhesus DVI+ red blood cells

5

| Panning Round No. ^{a)} | No. of eluted Rhesus D specific phages | |
|------------------------------------|--|--------------------|
| | Library LD1 cfu | Library LD2 cfu |
| 1 | 5×10^7 | 4.6×10^7 |
| 2 | 1.8×10^7 | 1.4×10^7 |
| 3 | 4×10^8 | 7.9×10^7 |
| 4 | 6.8×10^8 | 1.3×10^8 |
| 5 | 5.8×10^8 | 1×10^8 |

a) For each round 10^{12} Phabs were incubated in tubes with rbc Group O Rhesus negative (absorption phase) followed by elution from rbc Group O Rhesus DVI+

10 **Example 4**

Monitoring of the panning rounds and determination of the specificity of the enriched Phabs using a rabbit anti-phage antibody

Indirect haemagglutination assay

Freshly collected rbc of different ABO and Rhesus blood groups
 15 were washed 3 times in 0.9% NaCl and adjusted to a 3-5% solution ($45-50 \times 10^7$ /ml) in either 0.9% NaCl or PBS/3% BSA. For each test condition 50 μ l rbc and 100 μ l test (precipitated and amplified phage or control antibodies) were incubated together in glass blood grouping tubes (Baxter, Düringen, Switzerland) for 30 min. at 37°C. The rbc were washed 3 times in 0.9% NaCl
 20 and then incubated with 2 drops of Coombs reagent (Baxter, Düringen, Switzerland) for positive controls or with 100 μ l of 1/1000 diluted rabbit anti-phage antibodies (made by immunising rabbits with phage VCSM13

preparation, followed by purification on an Affi-Gel Blue column and absorption on *E. coli* to remove *E. coli*-specific antibodies). The tubes were incubated for 20 min at 37°C, centrifuged 1min at 125xg and rbc examined for agglutination by careful shaking and using a magnifier viewer.

- 5 When purified Fab were tested for agglutination an affinity purified anti-Fab antibody (The Binding Site, Birmingham, U.K.) was used instead of the rabbit anti-phage antibody.

Monitoring of Phabs from LD1 and LD2 librairies by indirect haemagglutination after Panning on R1R1 rbc

| Phab sample Panning Round | Library LD1 tested on rbc O Rh D+ (a) | Library LD2 |
|------------------------------|--|-------------|
| <i>No. 4</i> | | |
| undiluted | + | + |
| 1/4 | + | +/- |
| 1/20 | - | - |
| <i>No. 5</i> | | |
| undiluted | ++ | + |
| 1/4 | ++ | + |
| 1/20 | - | - |
| <i>No. 6</i> | | |
| undiluted | nd | +++ |
| 1/4 | nd | ++ |
| 1/20 | nd | nd |
| <i>Helper phage (b)</i> | | |
| undiluted, 1/4, 1/20 | - | - |

10

a) Indirect haemagglutination was performed in glass tubes using 50 µl rbc (40×10^7 /ml) and 100 µl Phabs starting at 4×10^{11} /ml. After 30 min. at 37°C the rbc were washed 3 times and further incubated for 20 min. at 37°C with a 1/1000 dilution of rabbit anti-phage antibody.

- 15 b) The M13 helper phage was used as a negative control and showed no non-specific agglutination due to the phage particle alone.

Agglutination was scored by visual assessment from +++ (strong agglutination) descending to - (no agglutination). nd = not done

Reaction Pattern of Fab Rhesus D Clones, resulting from R1R1 rbc panning, against Partial D Variants

| | Partial D Variants | | | | | |
|---|--------------------|------|------|------|-----|------|
| (a) Reaction Pattern of Clones from Library LD2 | Rh33 | DIII | DIVa | DIVb | DVI | DVII |
| <i>Pattern 1 n=2 (b)</i> | +++ | nd | +++ | +++ | - | +++ |
| <i>Pattern 2 n=1</i> | - | +++ | - | + | - | +++ |
| <i>Pattern 3 n=1</i> | - | +++ | +++ | + | - | +++ |
| <i>Pattern 4 n=2</i> | - | nd | +++ | +++ | - | +++ |
| <i>Pattern 5 n=1</i> | +++ | +++ | +++ | +++ | - | +++ |
| <i>Pattern 6 n=1</i> | - | +++ | - | - | - | +++ |
| <i>Pattern 7 n=1</i> | - | +++ | +++ | - | - | +++ |

- 5 a) soluble Fab preparations were made of each clone (as detailed in example 5) followed by indirect haemagglutination on the above panel of D variants
 b) n represents number of clones in each reaction pattern

Pattern 1: clones: LD2-1,14
 Pattern 2: clone: LD2-4
 Pattern 3: clone: LD2-17
 10 Pattern 4: clones: LD2-5,10
 Pattern 5: clone: LD2-18
 Pattern 6: clone: LD2-11
 Pattern 7: clone: LD2-20
 nd = not done

15 **Reaction Pattern of Rhesus D Clones against Partial DVI Variant**

Only the LD1-6-17 clone resulting from panning against DVI+rbc reacted with DVI+ and R1R1 rbc in the indirect haemagglutination assay using rabbit anti phage antibody as the Coombs equivalent reagent and in the presence of Enlist (Baxter).

Example 5

Preparation and purification of Fab antibody fragments for application as diagnostic reagents

After the bio-panning procedures detailed in Examples 2 and 3 a
 5 phage population which showed specific agglutination on Rhesus D+ rbc was selected and used to prepare phagemid DNA. More precisely the Phabs selected on R1R1 rbc were used after the 5th and 6th rounds of bio-panning for LD1 and LD2 libraries respectively and after the 5th bio-panning on DVI+ rbc for isolation of the LD1-6-17 clone. In order to produce soluble Fab, the
 10 sequence gIII coding for the pIII tail protein of the phage particle must be deleted.

Phagemid DNA was prepared using a Nucleotrap kit (Machery-Nagel) and the gIII sequence was removed by digesting the so isolated phagemid DNA with Nhe1/Spe1 as described (Burton, D.R., et al., PNAS,
 15 1989). After transformation into XL1-Blue individual clones were selected (nomenclature given in table 1) and grown in LB (Luria Broth) containing 50 µg/ml carbenicillin at 37°C to an OD of 0.6 at 600 nm. Cultures were induced with 2 mM isopropyl β-D-thiogalactopyranoside (IPTG) (Biofinex, Praroman, Switzerland) and grown overnight at 37°C. The whole culture was spun at
 20 10,000xg for 30 min. at 4°C to pellet the bacteria. The bacterial pellet was treated with a lysozyme/DNase solution to liberate the Fab fragments inside the cells. As some Fab were released into the culture supernatant this was also harvested separately. These Fab preparations were then pooled and precipitated with 60% ammonium sulphate (Merck, Darmstadt, Germany) to
 25 concentrate the Fab followed by extensive dialysis in phosphate buffered saline (PBS) and ultracentrifugation at 200,000xg to pellet any insoluble complexes. The Fab preparations were then purified on a ceramic hydroxyapatite column (HTP Econo cartridge, BioRad, Glattbrugg, Switzerland) using a gradient elution of PBS (Buffer A) and PBS + 0.5M NaCl
 30 (Buffer B). The linear gradient was programmed to increase from 0-100% Buffer B in 40 min. The Fab was eluted as a single peak between 40-60% Buffer B. The positive fractions as identified by immunodot assay using an

anti-Fab peroxidase conjugate (The Binding Site, Birmingham, U.K.) were pooled, concentrated using polyethylene glycol and extensively dialysed against PBS. The positive fractions from the hydroxyapatite column for each clone were used in a classical indirect haemagglutination assay in glass
 5 tubes using either the standard Coombs reagent (Baxter Diagnostics AG Dade, anti-human serum) or an anti-Fab (The Binding Site, Birmingham, U.K.) as the cross linking reagent. These Fab of defined specificity on the Partial D variants as shown on page 18 can be used to type rbc of unknown Partial D phenotype.

10 Example 6

Construction of complete immunoglobulin genes

The LD2-14 heavy chain V gene (V_H gene) was amplified from the anti-Rhesus D-Fab-encoding plasmid LD2-14 with the polymerase chain reaction (PCR) using gene-specific primers. The 5'-primer had the sequence
 15 5'-AGGTGTCGACGCACAGGTGAACTGCTCGAG-3' whereas the 3'-primer was of the sequence: 5'-GAGGAGACGGTGACCGTGGT-3'. The PCR reaction was performed with Taq DNA Polymerase and the buffer solution from Boehringer Mannheim (Mannheim, Germany) at the conditions recommended by the manufacturer including 60 pmol of each primer and the
 20 four deoxynucleotides at a concentration of 0.25 μ M each. The reaction was run for 35 cycles with the following temperature steps: 60 s at 94°C (extended by 5 min. during the first cycle), 75 s at 55°C and 90 s at 72°C (extended by 5 min. during the last cycle). The PCR product was purified with the QIAquick kit from Qiagen, Basel, Switzerland, digested with the restriction enzymes *Sal*
 25 *I* and *Bst*E II, extracted with phenol, purified by preparative agarose gel electrophoresis, concentrated by precipitation with ethanol and solubilized in a minimal volume of LTE buffer (10 mM Tris-HCl pH 8.0, 0.1 mM EDTA).

Vector # 150 (Sandoz Pharma, Basel) which contained an irrelevant intact human genomic immunoglobulin V_H gene was cut with *Sal* I
 30 and *Bst*E II, treated with calf intestinal phosphatase, extracted twice with

phenol, once with chloroform/isoamylalcohol (24:1), concentrated by precipitation with ethanol and solubilized in a minimal volume of LTE. 100 ng of this vector and 15 ng of the digested and purified PCR product were ligated according to the recommendations of the manufacturer in a total volume of 20 μ l (Ligase and buffer concentrate from Boehringer Mannheim (Mannheim)).
 5 The ligase was inactivated for 10 min. at 65°C, the reaction mix diluted to 100 μ l with H₂O and 3 μ l of the diluted solution were electroporated with a GenePulser (BioRad, Gaithersburg) into 40 μ l of competent *E. coli* JA221 according to the recommendations of the manufacturer (0.1 cm cuvettes, 1.8
 10 kV, 200 Ω , 25 μ FD), diluted with SOC medium, incubated at 37°C without shaking for 1 h and plated on LB plates containing ampicillin (50 μ g/ml). Plasmid-minipreps of the resulting colonies were checked with restriction digests for the presence of the appropriate insert.

With this procedure, the irrelevant resident V_H gene in vector # 150
 15 was replaced by the amplified anti-Rhesus D V_H sequence and yielded plasmid # 150-LD2-14. The structure of the resulting genomic immunoglobulin V_H gene construct was confirmed by sequencing, cut out by digestion with *EcoR* I and *Bam*H I and gel purified as described above for the PCR product. Expression vector # 10 (Sandoz Pharma, Basel) containing the human
 20 genomic immunoglobulin C γ 1 gene segment was also digested with *EcoR* I and *Bam*H I, treated with calf intestinal phosphatase, extracted with phenol, concentrated by precipitation with ethanol as described above, ligated with the *EcoR* I/*Bam*H I-V_H gene segments previously obtained from plasmids # 150-LD2-14 and transfected into *E. coli* JA221 as described above. This
 25 second cloning step yielded a complete anti-Rhesus D heavy chain immunoglobulin gene in plasmid # 10-LD2-14 (Figures 19 and 20).

The LD2-14 light chain V gene (V_L gene) was amplified from the same anti-Rhesus D-Fab plasmid LD2-14 by PCR using gene-specific primers. The 5'-primer had the sequence:
 30 5'-GGTACGCGTTGTGAGCTCGTGATGACCCAG-3'
 whereas the 3'-primer was of the sequence
 5'-TTTGATCTCAAGCTTGGTCCCAGGGCC-3'.

PCR reaction, product purification and subsequent cloning steps were performed as described for the V_H gene, except that the appropriate light chain vectors were used. Briefly, the V_L PCR product was digested with the restriction enzymes *Mlu*I and *Hind* III, extracted with phenol, purified from an agarose gel, ligated into vector # 151 (Sandoz Pharma, Basel) and propagated in *E. coli* JA221. This vector had been cut with *Mlu*I and *Hind* III, treated with calf intestinal phosphatase, extracted with phenol and concentrated by precipitation with ethanol. Plasmid # 151 contained an irrelevant intact human genomic immunoglobulin V_L gene. With this procedure, the irrelevant resident V_L gene was replaced by the amplified anti-Rhesus_D sequence and yielded plasmids # 151-LD2-14. The structure of the resulting genomic construct was confirmed by sequencing, cut out by digestion with *Eco*R I and *Xba* I and gel purified as described above. Subsequently, expression vector # 98 (Sandoz Pharma, Basel, Switzerland) containing the human genomic immunoglobulin C_k gene segment was digested with *Eco*R I and *Xba* I, treated with calf intestinal phosphatase, extracted with phenol, concentrated by precipitation with ethanol, ligated with the *Eco*R I/*Xba* I- V_L gene segment obtained from plasmid # 151-LD2-14 and transfected into *E. coli* JA221. This procedure yielded a complete anti-Rhesus D light chain immunoglobulin gene in plasmid # 98-LD2-14.

Restriction digests confirmed the structure of the expression constructs for LD2-14 heavy and light chain. The mouse myeloma cell line SP2/0-Ag 14 (ATCC CRL 1581) was cotransfected by electroporation with the plasmids # 10-LD2-14 and # 98-LD2-14. The electroporation was performed as follows: Exponentially growing cells were washed twice and suspended in phosphate buffered sucrose (272 mM sucrose, 1 mM $MgCl_2$, 7 mM NaH_2PO_4 , pH 7.4) at a density of 2×10^7 cells/ml. 0.8 ml of cells were added to a 0.4 cm cuvette, mixed with 15 μ g of linearized plasmids # 10-LD2-14 and # 98-LD2-14, held on ice for 15 min., electroporated with 290 Volts, 200 Ω , 25 μ FD, put back on ice for 15 min., transferred to a T75 cell culture flask with 20 ml of cold RPMI 1640 medium (10% heat inactivated fetal bovine serum, 50 μ M beta-mercaptoethanol), left for 2 h at room temperature and then incubated for 60 h at 37°C. After this period, the cells were transferred to 50 ml of medium containing 1 mg/ml G418 for selection. Two weeks later, all non-

transformed cells had died and the supernatants tested positive for the presence of human IgG/k by an enzyme linked immuno-sorbent assay (ELISA). The cultures were cloned by limiting dilution in microtiter plates and the supernatants quantitated by ELISA. The best producer clone was used for

5 production cultures. From 2 - 3 liters of culture supernatant the antibodies were purified by affinity chromatography on a Protein G Sepharose (Pharmacia, Uppsala) according to the recommendations of the manufacturer. Briefly, the supernatant was diluted with 1 volume of 10 mM phosphate buffer pH 7.0, 0.02% NaN_3 and pumped over a column with 2 ml of Protein G

10 Sepharose. The column was washed with 30 ml of 10 mM phosphate buffer pH 7.0, 0.02% NaN_3 and the bound antibodies were eluted with 0.1 M Glycine-HCl pH 2.0, 0.02% NaN_3 and immediately neutralised with solid NaHCO_3 . These antibodies (rD2-14) were positive for agglutination of Rhesus D positive and negative for Rhesus D negative red blood cells.

15 The same procedure was used to produce complete recombinant antibodies from the other clones from LD1 and LD2. Because of sequence diversity, for some clones different PCR primers were synthesized. For the heavy chain this was LD1-84, where the sequence of the 5'-primer was 5'-AGGTGTCGACGCACAGGTAACTGCTCGAG-3'. For the light chain the 5'-

20 primer was the same as for LD2-14 noted above (except for LD2-1, -10, where the sequence was 5'-GGTACGCGTTGTGA-GCTCGTGGTGACTCAG-3') whereas the 3'-primer was of the sequence 5'-TTTGATTTCAAGC-TTGGTCCCTTGGCC-3', except for some clones where the following 3'-primers were used:

25 LD1-52, -98: 5'-TTTGATCTCAAGCTTGGTCCCAGGGCC-3',

LD1-40: 5'-TTTGATCTCAAGCTTTGTCCCTTGGCC-3',

LD1-84: 5'-TTTGATGTCAAGCTTGGTCCCCCGGCC-3',

LD2-11: 5'-TTTGATCTGAAGCTTGGTCCCCTGCCC-3'.

LD2-1, -10: 5'-TAGGACGGTAAGCTTGGTCCCTCCGCC-3').

The cloning steps for the production of complete antibody heavy and light chain genes were the same as those described above for LD2-14. Expression in SP2/0-Ag14 cells yielded antibodies with the same anti-Rhesus D specificities as the Fab expressed in *E. coli*.

Claims

1. Polypeptides capable of forming antigen binding structures with specificity for Rhesus D antigens which include Rhesus D-specific CDR 1, CDR 2 and CDR 3 regions of the amino acid sequences V_H and V_L with the identification numbers according to the figures given in the table below:

| Identification No. | V_H | | | | V_L | | | |
|--------------------|----------|---------------------------|---------------------------|---------------------------|----------|---------------------------|---------------------------|---------------------------|
| | Figure | CDR 1 base pair No. | CDR 2 base pair No. | CDR 3 base pair No. | Figure | CDR 1 base pair No. | CDR 2 base pair No. | CDR 3 base pair No. |
| LD1-28 | Fig. 1a | 91-104 | 148-198 | 295-342 | Fig. 1b | 64-96 | 142-162 | 259-288 |
| LD1-40 | Fig. 2a | 91-105 | 148-198 | 295-342 | Fig. 2b | 64-96 | 142-162 | 259-288 |
| LD1-52 | Fig. 3a | 91-105 | 148-198 | 295-342 | Fig. 3b | 64-96 | 142-162 | 259-288 |
| LD1-84 | Fig. 4a | 91-105 | 148-198 | 295-342 | Fig. 4b | 64-96 | 142-162 | 259-285 |
| LD1-98 | Fig. 5a | 91-105 | 148-198 | 295-342 | Fig. 5b | 64-96 | 142-162 | 259-288 |
| LD1-110 | Fig. 6a | 91-105 | 148-198 | 295-342 | Fig. 6b | 64-96 | 142-162 | 259-285 |
| LD1-117 | Fig. 7a | 91-105 | 148-198 | 295-345 | Fig. 7b | 64-96 | 142-162 | 259-288 |
| LD2-1 | Fig. 8a | 91-105 | 148-198 | 295-342 | Fig. 8b | 61-99 | 145-165 | 262-294 |
| LD2-4 | Fig. 9a | 91-105 | 148-198 | 295-342 | Fig. 9b | 64-96 | 142-162 | 259-282 |
| LD2-5 | Fig. 10a | 91-105 | 148-198 | 295-342 | Fig. 10b | 64-96 | 142-162 | 259-288 |
| LD2-10 | Fig. 11a | 91-105 | 148-198 | 298-345 | Fig. 11b | 61-102 | 148-168 | 265-294 |
| LD2-11 | Fig. 12a | 91-105 | 148-198 | 295-342 | Fig. 12b | 64-96 | 142-162 | 259-285 |
| LD2-14 | Fig. 13a | 91-105 | 148-198 | 295-342 | Fig. 13b | 64-96 | 142-162 | 259-285 |
| LD2-17 | Fig. 14a | 91-105 | 148-198 | 295-342 | Fig. 14b | 64-96 | 142-162 | 259-285 |
| LD2-18 | Fig. 15a | 91-105 | 148-198 | 295-342 | Fig. 15b | 64-96 | 142-162 | 259-288 |
| LD2-20 | Fig. 16a | 91-105 | 148-198 | 295-342 | Fig. 16b | 64-96 | 142-162 | 259-285 |
| LD1-6-17 | Fig. 17a | 91-105 | 148-198 | 295-351 | Fig. 17b | 64-96 | 142-162 | 259-285 |

2. Polypeptides according to claim 1 which include regions with the amino acid sequences V_H and V_L and have identification numbers according to the figures given in the table of claim 1.

3. Polypeptides according to claim 1 or 2 characterised as antigen binding Fab fragments.

4. Polypeptides according to claim 1 or 2 comprising immunoglobulin heavy and light chains capable of forming complete anti-
5 Rhesus D antibodies.

5. DNA sequences coding for polypeptides capable of forming antigen binding structures with specificity for Rhesus D antigens which include regions with the Rhesus D-specific CDR 1, CDR 2 and CDR 3 segments of the DNA sequences V_H and V_L with the identification numbers
10 according to the figures given in the table below, and functional equivalents thereof:

| Identi- fication No. | V_H | | | | V_L | | | |
|----------------------------|----------|---------------------------|---------------------------|---------------------------|----------|---------------------------|---------------------------|---------------------------|
| | Figure | CDR 1 base pair No. | CDR 2 base pair No. | CDR 3 base pair No. | Figure | CDR 1 base pair No. | CDR 2 base pair No. | CDR 3 base pair No. |
| LD1-28 | Fig. 1a | 91-104 | 148-198 | 295-342 | Fig. 1b | 64-96 | 142-162 | 259-288 |
| LD1-40 | Fig. 2a | 91-105 | 148-198 | 295-342 | Fig. 2b | 64-96 | 142-162 | 259-288 |
| LD1-52 | Fig. 3a | 91-105 | 148-198 | 295-342 | Fig. 3b | 64-96 | 142-162 | 259-288 |
| LD1-84 | Fig. 4a | 91-105 | 148-198 | 295-342 | Fig. 4b | 64-96 | 142-162 | 259-285 |
| LD1-98 | Fig. 5a | 91-105 | 148-198 | 295-342 | Fig. 5b | 64-96 | 142-162 | 259-288 |
| LD1-110 | Fig. 6a | 91-105 | 148-198 | 295-342 | Fig. 6b | 64-96 | 142-162 | 259-285 |
| LD1-117 | Fig. 7a | 91-105 | 148-198 | 295-345 | Fig. 7b | 64-96 | 142-162 | 259-288 |
| LD2-1 | Fig. 8a | 91-105 | 148-198 | 295-342 | Fig. 8b | 61-99 | 145-165 | 262-294 |
| LD2-4 | Fig. 9a | 91-105 | 148-198 | 295-342 | Fig. 9b | 64-96 | 142-162 | 259-282 |
| LD2-5 | Fig. 10a | 91-105 | 148-198 | 295-342 | Fig. 10b | 64-96 | 142-162 | 259-288 |
| LD2-10 | Fig. 11a | 91-105 | 148-198 | 298-345 | Fig. 11b | 61-102 | 148-168 | 265-294 |
| LD2-11 | Fig. 12a | 91-105 | 148-198 | 295-342 | Fig. 12b | 64-96 | 142-162 | 259-285 |
| LD2-14 | Fig. 13a | 91-105 | 148-198 | 295-342 | Fig. 13b | 64-96 | 142-162 | 259-285 |
| LD2-17 | Fig. 14a | 91-105 | 148-198 | 295-342 | Fig. 14b | 64-96 | 142-162 | 259-285 |
| LD2-18 | Fig. 15a | 91-105 | 148-198 | 295-342 | Fig. 15b | 64-96 | 142-162 | 259-288 |
| LD2-20 | Fig. 16a | 91-105 | 148-198 | 295-342 | Fig. 16b | 64-96 | 142-162 | 259-285 |
| LD1-6-17 | Fig. 17a | 91-105 | 148-198 | 295-351 | Fig. 17b | 64-96 | 142-162 | 259-285 |

6. DNA sequences according to claim 5 which include regions with the DNA sequences V_H and V_L with the identification numbers according to the figures given in claim 5.

7. DNA sequences according to claim 5 or 6 coding for polypeptides capable of forming antigen binding Fab fragments.

8. DNA sequences according to claim 5 or 6 coding for polypeptides capable of forming complete anti-Rhesus D antibodies.

9. A process for preparing recombinant polypeptides capable of forming antigen binding structures, e.g. Fab fragments, with specificity for Rhesus D antigens which process comprises the following steps in sequential order:

- a) boosting of an individual capable of forming anti-Rhesus D antibodies with Rhesus D positive red blood cells,
- b) isolating mononuclear cells from the individual,
- 15 c) isolating total RNA from the mononuclear cells,
- d) preparing a cDNA by using an oligo(dT)primer and reverse transcribing of the mRNA with M-MuLV reverse transcriptase and amplifying the cDNA repertoire by a polymerase chain reaction using immunoglobulin gene family specific primers,
- 20 e) creating a phage display library by inserting the DNA coding for the heavy and light chain of the Fab polypeptide into a phagemid vector; the DNA for the heavy chain is inserted in frame to the gene coding for the phage protein pIII which allows the expression of a Fab pIII fusion protein on the surface of the phage,
- 25 f) transforming bacterial cells with the obtained recombinant plasmids, cultivating of the transformed bacterial cells and co-expression of the heavy and the light chain of a Fab on filamentous phage particles,
- g) amplifying the Fab-carrying phage in bacteria,
- h) selecting individual phage clones by several rounds of panning on
30 Rhesus positive red blood cells.

- i) isolating the plasmid DNA from the selected clones and cutting out the gIII gene,
- j) transforming bacterial cells with the obtained plasmid, cultivating of the transformed bacterial cells expressing the Fab, and isolating the Fab fragments.

5

10. Anti-Rhesus D antibodies having heavy and light chain variable regions comprising the Rhesus D-specific CDR 1, CDR 2 and CDR 3 sequences of the amino acid sequences V_H and V_L given the identification numbers as indicated in the table below:

| Identification No. | V_H | | | | V_L | | | |
|--------------------|----------|---------------------------|---------------------------|---------------------------|----------|---------------------------|---------------------------|---------------------------|
| | Figure | CDR 1 base pair No. | CDR 2 base pair No. | CDR 3 base pair No. | Figure | CDR 1 base pair No. | CDR 2 base pair No. | CDR 3 base pair No. |
| LD1-28 | Fig. 1a | 91-104 | 148-198 | 295-342 | Fig. 1b | 64-96 | 142-162 | 259-288 |
| LD1-40 | Fig. 2a | 91-105 | 148-198 | 295-342 | Fig. 2b | 64-96 | 142-162 | 259-288 |
| LD1-52 | Fig. 3a | 91-105 | 148-198 | 295-342 | Fig. 3b | 64-96 | 142-162 | 259-288 |
| LD1-84 | Fig. 4a | 91-105 | 148-198 | 295-342 | Fig. 4b | 64-96 | 142-162 | 259-285 |
| LD1-98 | Fig. 5a | 91-105 | 148-198 | 295-342 | Fig. 5b | 64-96 | 142-162 | 259-288 |
| LD1-110 | Fig. 6a | 91-105 | 148-198 | 295-342 | Fig. 6b | 64-96 | 142-162 | 259-285 |
| LD1-117 | Fig. 7a | 91-105 | 148-198 | 295-345 | Fig. 7b | 64-96 | 142-162 | 259-288 |
| LD2-1 | Fig. 8a | 91-105 | 148-198 | 295-342 | Fig. 8b | 61-99 | 145-165 | 262-294 |
| LD2-4 | Fig. 9a | 91-105 | 148-198 | 295-342 | Fig. 9b | 64-96 | 142-162 | 259-282 |
| LD2-5 | Fig. 10a | 91-105 | 148-198 | 295-342 | Fig. 10b | 64-96 | 142-162 | 259-288 |
| LD2-10 | Fig. 11a | 91-105 | 148-198 | 298-345 | Fig. 11b | 61-102 | 148-168 | 265-294 |
| LD2-11 | Fig. 12a | 91-105 | 148-198 | 295-342 | Fig. 12b | 64-96 | 142-162 | 259-285 |
| LD2-14 | Fig. 13a | 91-105 | 148-198 | 295-342 | Fig. 13b | 64-96 | 142-162 | 259-285 |
| LD2-17 | Fig. 14a | 91-105 | 148-198 | 295-342 | Fig. 14b | 64-96 | 142-162 | 259-285 |
| LD2-18 | Fig. 15a | 91-105 | 148-198 | 295-342 | Fig. 15b | 64-96 | 142-162 | 259-288 |
| LD2-20 | Fig. 16a | 91-105 | 148-198 | 295-342 | Fig. 16b | 64-96 | 142-162 | 259-285 |
| LD1-6-17 | Fig. 17a | 91-105 | 148-198 | 295-351 | Fig. 17b | 64-96 | 142-162 | 259-285 |

11. Anti-Rhesus D antibodies according to claim 10 with the amino acid sequences V_H and V_L and the identification numbers according to the figures, as indicated in the table of claim 10.

12. Anti-Rhesus D antibodies according to claim 10 or 11 wherein
 5 the immunoglobulin constant regions are of at least one of the defined isotypes IgG1, IgG2, IgG3 or IgG4.

13. A process for preparing complete anti-Rhesus D antibodies according to one of the claims 10 to 12, comprising in sequential order the steps of

10 a) amplifying separately the members of a pair of a heavy chain V gene segment and a light chain V gene segment containing Rhesus D-specific CDR 1, CDR 2 and CDR 3 regions as depicted in Figs. 1a - 17a and 1b - 17b, respectively, from an anti-Rhesus D-Fab-encoding plasmid by carrying out a polymerase chain reaction with specific
 15 primers,

b) preparing separately the genes of a complete anti-Rhesus D immunoglobulin heavy chain and a complete anti-Rhesus D immunoglobulin light chain in suitable plasmids containing the immunoglobulin constant region gene segments coding for either one
 20 of the human $\gamma 1$, $\gamma 2$, $\gamma 3$ and $\gamma 4$ heavy chains and for the human κ or λ light chain and transforming the obtained plasmids separately in suitable E. coli bacteria, and

c) cotransfecting the obtained plasmids into a suitable mouse myeloma cell line, cultivating of the cells, separating the non-transformed cells,
 25 cloning of the cultures, selecting the best producing clone, using it as a production culture and isolating the complete antibodies from the supernatant of the cell culture.

14. A pharmaceutical composition comprising at least one polypeptide according to the definition of claim 1 or 2 or at least one anti-
 30 Rhesus D antibody according to one of the claims 10 to 12 for the prophylaxis

of haemolytic disease of the newborn, for the treatment of idiopathic thrombocytopenic purpura and mistransfusions of Rhesus incompatible blood.

15. A diagnostic composition for Rhesus D typing comprising Fab fragments according to claim 3 or anti-Rhesus D antibodies according to one
5 of the claims 10 to 12.

Abstract

Polypeptides capable of forming antigen binding structures specific for Rhesus D antigens include the sequences indicated in the figs. 1a to 17 b. They are prepared in a process comprising the following steps in sequential order:

- boosting of an individual capable of forming anti-Rhesus D antibodies with Rhesus D positive red blood cells,
- isolating mononuclear cells from the individual,
- isolating total RNA from the mononuclear cells,
- preparing a cDNA by using an oligo(dT)primer and reverse transcribing of the mRNA with M-MuLV reverse transcriptase and amplifying the cDNA repertoire by a polymerase chain reaction using immunoglobulin gene family specific primers,
- creating a phage display library by inserting the DNA coding for the heavy and light chain of the Fab polypeptide into a phagemid vector; the DNA for the heavy chain is inserted in frame into the gene coding for the phage protein pIII which allows the expression of a Fab pIII fusion protein on the surface of the phage,
- transforming bacterial cells with the obtained recombinant plasmids, cultivation of the transformed bacterial cells and co-expression of the heavy and the light chain of a Fab on filamentous phage particles,
- amplifying the Fab-carrying phage in bacteria,
- selecting individual phage clones by several rounds of panning on Rhesus positive red blood cells,
- isolating the plasmid DNA from the selected clones and cutting out the gIII gene,
- transforming bacterial cells with the obtained plasmid, cultivating of the transformed bacterial cells expressing the Fab, and isolating the Fab fragments.

The obtained polypeptides, being Fab fragments, can be used directly as an active ingredient in pharmaceutical and diagnostic compositions. The Fab and their DNA sequences can also be used for the preparation of complete recombinant Anti-Rhesus D antibodies.

(1) GENERAL INFORMATION

(i) APPLICANT:

- (A) NAME: Rotkreuzstiftung Zentrallaboratorium Blutspendedienst SRK
- (B) STREET: Wankdorfstrasse 10
- (C) CITY: Bern 22
- (D) STATE OR PROVINCE:
- (E) COUNTRY: Switzerland
- (F) POSTAL CODE: CH-3000

(ii) TITLE OF INVENTION: Polypeptides capable of forming antigen binding structures with specificity for the Rhesus D antigens, the DNA encoding them and the process for their preparation and use

(iii) NUMBER OF SEQUENCES: 34

(iv) COMPUTER-READABLE FORM:

- (A) MEDIUM TYPE: 3.5" Floppy disk, 1.44 MB
- (B) COMPUTER: IBM compatible PC
- (C) OPERATING SYSTEM: IBM-DOS 6.3/Windows 3.1
- (D) SOFTWARE: MS Word for Windows 6.0 /saved as MS-DOS text

(v) CURRENT APPLICATION DATA: n.a.

- (A) APPLICATION NUMBER: n.a.

(2) INFORMATION FOR SEQ ID NO: 1 (LD1-28-VH)
(i) SEQUENCE CHARACTERISTICS:
(A) LENGTH: 375 Base Pairs
(B) TYPE: Nucleic acid with corresponding protein
(C) STRANDEDNESS: Single
(D) TOPOLOGY: Linear
(ii) MOLECULE TYPE: cDNA to mRNA
(iii) HYPOTHETICAL: No
(iv) ANTI-SENSE: No
(vi) ORIGINAL SOURCE: N-terminal fragment
(A) ORGANISM: Homo sapiens
(B) STRAIN: Not applicable
(C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
(D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
(E) HAPLOTYPE: Diploid
(F) TISSUE TYPE: Not applicable
(G) CELL TYPE: Peripheral lymphocyte B
(H) CELL LINE: Not applicable
(I) ORGANELLE: Not applicable
(vii) IMMEDIATE SOURCE: Escherichia coli
(A) LIBRARY: cDNA LIBRARY, LD1
(B) CLONE: LD1-28
(viii) POSITION IN GENOME:
(A) CHROMOSOME / SEGMENT: Chromosome 14
(B) MAP POSITION: q32.3
(C) UNITS: Chromosome band number
(xi) FEATURE:
(A) NAME / KEY: CDR1, CDR2, CDR3
(B) LOCATION: join base pair (91..105, 148..198, 295..342)
(C) IDENTIFICATION METHOD: By similarity with a known sequence
(D) OTHER INFORMATION:

(xi) SEQUENCE DESCRIPTION: SEQ ID NO. 1 (LD1-28-VH)

| | |
|---|-----|
| CAG GTG AAA CTG CTC GAG TCT GGG GGA GGC GTG GTC CAG CCG GGG GGG | 48 |
| Gln Val Lys Leu Leu Glu Ser Gly Gly Gly Val Val Gln Pro Gly Gly | |
| 5 10 15 | |
| TCC CTG AGA CTC TCC TGT GAA GCG TCT GGA TTC GCC CTC AGA AGT TCT | 96 |
| Ser Leu Arg Leu Ser Cys Glu Ala Ser Gly Phe Ala Leu Arg Ser Ser | |
| 20 25 30 | |
| GGC ATG CAC TGG GTC CGC CAG GCT CCT GGC AAG GGG ATG GAG TGG GTG | 144 |
| Gly Met His Typ Val Arg Gln Ala Pro Gly Lys Gly Met Glu Typ Val | |
| 35 40 45 | |
| GCA CTT ATA TGG TTT GAT GGA AGT ATC AGA TCG TAT GCA GAA TCC GCG | 192 |
| Ala Leu Ile Typ Phe Asp Gly Ser Ile Arg Ser Tyr Ala Glu Ser Ala | |
| 50 55 60 | |
| AAG GGC CGA TTC ACC ATC TCC AGA GAC ACT TCC AAG AAC ACC CTA TAT | 240 |
| Lys Gly Arg Phe Thr Ile Ser Arg Asp Thr Ser Lys Asn Thr Leu Tyr | |
| 65 70 75 80 | |
| CTC CAT ATG CGC AGT CTG AGT GCC GAC GAC ACG GAT GTG TTT TAC TGT | 288 |
| Leu His Met Arg Ser Leu Ser Ala Asp Asp Thr Asp Val Phe Tyr Cys | |
| 85 90 95 | |
| GCG AGA GAC AAG GCG GTT CGG GGA ATT AAC AGG TAC AAC TAT TAC ATG | 336 |
| Ala Arg Asp Lys Ala Val Arg Gly Ile Asn Arg Tyr Asn Tyr Tyr Met | |
| 100 105 110 | |
| GAC GTC TGG GTC AAA GGG ACC ACG GTC ACC GTC TCC TCA | 375 |
| Asp Val Typ Val Lys Gly Thr Thr Val Thr Val Ser Ser | |
| 115 120 125 | |

- (2) INFORMATION FOR SEQ ID NO:2 (LD1-28-VL)
- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 318 Base Pairs
 - (B) TYPE: Nucleic acid with corresponding protein
 - (C) STRANDEDNESS: Single
 - (D) TOPOLOGY: Linear
- (ii) MOLECULE TYPE: cDNA to mRNA
- (iii) HYPOTHETICAL: No
- (iv) ANTI-SENSE: No
- (vi) ORIGINAL SOURCE: N-terminal fragment
- (A) ORGANISM: Homo sapiens
 - (B) STRAIN: Not applicable
 - (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
 - (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
 - (E) HAPLOTYPE: Diploid
 - (F) TISSUE TYPE: Not applicable
 - (G) CELL TYPE: Peripheral lymphocyte B
 - (H) CELL LINE: Not applicable
 - (I) ORGANELLE: Not applicable
- (vii) IMMEDIATE SOURCE: Escherichia coli
- (A) LIBRARY: cDNA LIBRARY, LD1
 - (B) CLONE: LD1-28
- (viii) POSITION IN GENOME:
- (A) CHROMOSOME / SEGMENT: Chromosome 2
 - (B) MAP POSITION: p12
 - (C) UNITS: Chromosome band number
- (xi) FEATURE:
- (A) NAME / KEY: CDR1, CDR2, CDR3
 - (B) LOCATION: join base pair (64..96, 142..162, 259..288)
 - (C) IDENTIFICATION METHOD: By similarity with a known sequence
 - (D) OTHER INFORMATION:

- (2) INFORMATION FOR SEQ ID NO: 3 (LD1-40-VH)
- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 375 Base Pairs
 - (B) TYPE: Nucleic acid with corresponding protein
 - (C) STRANDEDNESS: Single
 - (D) TOPOLOGY: Linear
- (ii) MOLECULE TYPE: cDNA to mRNA
- (iii) HYPOTHETICAL: No
- (iv) ANTI-SENSE: No
- (vi) ORIGINAL SOURCE: N-terminal fragment
- (A) ORGANISM: Homo sapiens
 - (B) STRAIN: Not applicable
 - (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
 - (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
 - (E) HAPLOTYPE: Diploid
 - (F) TISSUE TYPE: Not applicable
 - (G) CELL TYPE: Peripheral lymphocyte B
 - (H) CELL LINE: Not applicable
 - (I) ORGANELLE: Not applicable
- (vii) IMMEDIATE SOURCE: Escherichia coli
- (A) LIBRARY: cDNA LIBRARY, LD1
 - (B) CLONE: LD1-40
- (viii) POSITION IN GENOME:
- (A) CHROMOSOME / SEGMENT: Chromosome 14
 - (B) MAP POSITION: q 32.3
 - (C) UNITS: Chromosome band number
- (xi) FEATURE:
- (A) NAME / KEY: CDR1, CDR2, CDR3
 - (B) LOCATION: join base pair (91..105, 148..198, 295..342)
 - (C) IDENTIFICATION METHOD: By similarity with a known sequence
 - (D) OTHER INFORMATION:

(xi) SEQUENCE DESCRIPTION: SEQ ID NO. 3 (LD1-40-VH)

| | |
|---|-----|
| CAG GTG AAA CTG CTC GAG TCT GGG GGA GGC GTG GTC CAG CCT GGG AGG | 48 |
| Gln Val Lys Leu Leu Glu Ser Gly Gly Gly Val Val Gln Pro Gly Arg | |
| 5 10 15 | |
| TCC CTG AGA CTC TCC TGT ATA GCG TCT GGA TTC ACC CTC AGG AAT TAT | 96 |
| Ser Leu Arg Leu Ser Cys Ile Ala Ser Gly Phe Thr Leu Arg Asn Tyr | |
| 20 25 30 | |
| GCC ATG CAC TGG GTC CGC CAG GCT CCA GGC AAG GGG CTG GAG TGG GTG | 144 |
| Ala Met His Typ Val Arg Gln Ala Pro Gly Lys Gly Leu Glu Typ Val | |
| 35 40 45 | |
| GCA GGT ATA TGG TTT GAT GGA AGT AAC AAA AAC TAT GCA GAC TCC GTG | 192 |
| Ala Gly Ile Typ Phe Asp Gly Ser Asn Lys Asn Tyr Ala Asp Ser Val | |
| 50 55 60 | |
| AAG GGC CGA TTC ACC ATC TCC AGA GAC AAT TCC AAG AAC ACG CTG TTT | 240 |
| Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ser Lys Asn Thr Leu Phe | |
| 65 70 75 80 | |
| CTG CAA CTG AAC AGC CTG AGA GAC GAG GAC ACG GCT GTG TAT TAT TGT | 288 |
| Leu Gln Leu Asn Ser Leu Arg Asp Glu Asp Thr Ala Val Tyr Tyr Cys | |
| 85 90 95 | |
| GCG AGA GAG CGA GCA GCA CGT GGT ATT TCT AGG TTC TAT TAC TAC ATG | 336 |
| Ala Arg Glu Arg Ala Ala Arg Gly Ile Ser Arg Phe Tyr Tyr Tyr Met | |
| 100 105 110 | |
| GAC GTC TGG GGC AAA GGG ACC ACG GTC ACC GTC TCC CCA | 375 |
| Asp Val Typ Gly Lys Gly Thr Thr Val Thr Val Ser Pro | |
| 115 120 125 | |

(2) INFORMATION FOR SEQ ID NO: 4 (LD1-40-VL)

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 318 Base Pairs
- (B) TYPE: Nucleic acid with corresponding protein
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

- (A) ORGANISM: Homo sapiens
- (B) STRAIN: Not applicable
- (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
- (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
- (E) HAPLOTYPE: Diploid
- (F) TISSUE TYPE: Not applicable
- (G) CELL TYPE: Peripheral lymphocyte B
- (H) CELL LINE: Not applicable
- (I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

- (A) LIBRARY: cDNA LIBRARY, LD1
- (B) CLONE: LD1-40

(viii) POSITION IN GENOME:

- (A) CHROMOSOME / SEGMENT: Chromosome 2
- (B) MAP POSITION: 2p12
- (C) UNITS: Chromosome band number

(xi) FEATURE:

- (A) NAME / KEY: CDR1, CDR2, CDR3
- (B) LOCATION: join base pair (64..96, 142..162, 259..288)
- (C) IDENTIFICATION METHOD: By similarity with a known sequence
- (D) OTHER INFORMATION:

(2) INFORMATION FOR SEQ ID NO: 5 (LD1-52-VH)

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 375 Base Pairs

(B) TYPE: Nucleic acid with corresponding protein

(C) STRANDEDNESS: Single

(D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

(A) ORGANISM: Homo sapiens

(B) STRAIN: Not applicable

(C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor

(D) DEVELOPMENT STAGE: Adult, Rearranged development pattern

(E) HAPLOTYPE: Diploid

(F) TISSUE TYPE: Not applicable

(G) CELL TYPE: Peripheral lymphocyte B

(H) CELL LINE: Not applicable

(I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

(A) LIBRARY: cDNA LIBRARY, LD1

(B) CLONE: LD1-52

(viii) POSITION IN GENOME:

(A) CHROMOSOME / SEGMENT: Chromosome14

(B) MAP POSITION: q32.3

(C) UNITS: Chromosome band number

(xi) FEATURE:

(A) NAME / KEY: CDR1, CDR2, CDR3

(B) LOCATION: join base pair (91..105, 148..198, 295..342)

(C) IDENTIFICATION METHOD: By similarity with a known sequence

(D) OTHER INFORMATION:

(xi) SEQUENCE DESCRIPTION: SEQ ID NO.5 (LD1-52-VH)

| | |
|---|-----|
| CAG GTG AAA CTG CTC GAG TCT GGG GGA GGC GTG GTC CAG CCG GGG GGG | 48 |
| Gln Val Lys Leu Leu Glu Ser Gly Gly Gly Val Val Gln Pro Gly Gly | |
| 5 10 15 | |
| TCC CTG AGA CTC TCC TGT GAA GCG TCT GGA TTC GCC CTC AGA AGT TCT | 96 |
| Ser Leu Arg Leu Ser Cys Glu Ala Ser Gly Phe Ala Leu Arg Ser Ser | |
| 20 25 30 | |
| GGA ATG CAC TGG GTC CGC CAG GCT CCT GGC AAG GGG CTG GAG TGG GTG | 144 |
| Gly Met His Typ Val Arg Gln Ala Pro Gly Lys Gly Leu Glu Typ Val | |
| 35 40 45 | |
| GCA CTT ATA TGG TTT GAT GGA AGT ATC AGA TCG TAT GCA GAA TCC GTG | 192 |
| Ala Leu Ile Typ Phe Asp Gly Ser Ile Arg Ser Tyr Ala Glu Ser Val | |
| 50 55 60 | |
| AAG GGC CGA TTC ACC ATC TCC AGA GAC ACT TCC AAG AAC ACC CTA TAT | 240 |
| Lys Gly Arg Phe Thr Ile Ser Arg Asp Thr Ser Lys Asn Thr Leu Tyr | |
| 65 70 75 80 | |
| CTC CAA ATG CGC AGT CTG AGT GCC GAC GAC ACG GCT GTG TAT TAC TGT | 288 |
| Leu Gln Met Arg Ser Leu Ser Ala Asp Asp Thr Ala Val Tyr Tyr Cys | |
| 85 90 95 | |
| GCG AGA GAC AAG GCG GTT CGG GGA ATT AGC AGG TAC AAC TAT TAC ATG | 336 |
| Ala Arg Asp Lys Ala Val Arg Gly Ile Ser Arg Tyr Asn Tyr Tyr Met | |
| 100 105 110 | |
| GAC GTC TGG GGC AAA GGG ACC ACG GTC ACC GTC TCC TCA | 375 |
| Asp Val Typ Gly Lys Gly Thr Thr Val Thr Val Ser Ser | |
| 115 120 125 | |

(2) INFORMATION FOR SEQ ID NO: 6 (LD1-52-VL)

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 318 Base Pairs
- (B) TYPE: Nucleic acid with corresponding protein
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

- (A) ORGANISM: Homo sapiens
- (B) STRAIN: Not applicable
- (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
- (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
- (E) HAPLOTYPE: Diploid
- (F) TISSUE TYPE: Not applicable
- (G) CELL TYPE: Peripheral lymphocyte B
- (H) CELL LINE: Not applicable
- (I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

- (A) LIBRARY: cDNA LIBRARY, LD1
- (B) CLONE: LD1-52

(viii) POSITION IN GENOME:

- (A) CHROMOSOME / SEGMENT: Chromosome 2
- (B) MAP POSITION: p12
- (C) UNITS: Chromosome band number

(xi) FEATURE:

- (A) NAME / KEY: CDR1, CDR2, CDR3
- (B) LOCATION: join base pair (64..96, 142..162, 259..288)
- (C) IDENTIFICATION METHOD: By similarity with a known sequence
- (D) OTHER INFORMATION:

(2) INFORMATION FOR SEQ ID NO:7 (LD1-84-VH)

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 375 Base Pairs

(B) TYPE: Nucleic acid with corresponding protein

(C) STRANDEDNESS: Single

(D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

(A) ORGANISM: Homo sapiens

(B) STRAIN: Not applicable

(C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor

(D) DEVELOPMENT STAGE: Adult, Rearranged development pattern

(E) HAPLOTYPE: Diploid

(F) TISSUE TYPE: Not applicable

(G) CELL TYPE: Peripheral lymphocyte B

(H) CELL LINE: Not applicable

(I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

(A) LIBRARY: cDNA LIBRARY, LD1

(B) CLONE: LD1-84

(viii) POSITION IN GENOME:

(A) CHROMOSOME / SEGMENT: Chromosome 14

(B) MAP POSITION: q32.3

(C) UNITS: Chromosome band number

(xi) FEATURE:

(A) NAME / KEY: CDR1, CDR2, CDR3

(B) LOCATION: join base pair (91..105, 148..198, 295..342)

(C) IDENTIFICATION METHOD: By similarity with a known sequence

(D) OTHER INFORMATION:

(xi) SEQUENCE DESCRIPTION: SEQ ID NO.7 (LD1-84-VH)

| | | | | | | | | | | | | | | | | |
|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| CAG | GTA | AAA | CTG | CTC | GAG | TCT | GGG | GGA | GGC | GTG | GTC | CAG | CCG | GGG | GGG | 48 |
| Gln | Val | Lys | Leu | Leu | Glu | Ser | Gly | Gly | Gly | Val | Val | Gln | Pro | Gly | Gly | |
| | | | 5 | | | | | | 10 | | | | | 15 | | |
| TCC | CTG | AGA | CTC | TCC | TGT | GAA | GCG | TCT | GGA | TTC | ACC | CTC | AGA | AGT | TCT | 96 |
| Ser | Leu | Arg | Leu | Ser | Cys | Glu | Ala | Ser | Gly | Phe | Thr | Leu | Arg | Ser | Ser | |
| | | | 20 | | | | | | 25 | | | | | 30 | | |
| GGC | ATG | CAC | TGG | GTC | CGC | CAG | GCT | CCT | GGC | AAG | GGG | CTG | GAG | TGG | GTG | 144 |
| Gly | Met | His | Typ | Val | Arg | Gln | Ala | Pro | Gly | Lys | Gly | Leu | Glu | Typ | Val | |
| | | 35 | | | | | 40 | | | | | 45 | | | | |
| ACA | CTT | ATA | TGG | TTT | GAT | GGA | AGT | ATC | AGA | TCG | TAT | GCA | GAA | TCC | GTG | 192 |
| Thr | Leu | Ile | Typ | Phe | Asp | Gly | Ser | Ile | Arg | Ser | Tyr | Ala | Glu | Ser | Val | |
| | | 50 | | | | 55 | | | | | | 60 | | | | |
| AAG | GGC | CGA | TTC | ACC | ATC | TCC | AGA | GAC | ACT | TCC | AAG | AAC | ACC | CTA | TAT | 240 |
| Lys | Gly | Arg | Phe | Thr | Ile | Ser | Arg | Asp | Thr | Ser | Lys | Asn | Thr | Leu | Tyr | |
| 65 | | | | | 70 | | | | | 75 | | | | 80 | | |
| CTC | CAA | ATG | CGC | AGT | CTG | AGT | GCC | GAC | GAC | ACG | GCT | GTG | TAT | TAC | TGT | 288 |
| Leu | Gln | Met | Arg | Ser | Leu | Ser | Ala | Asp | Asp | Thr | Ala | Val | Tyr | Tyr | Cys | |
| | | | | 85 | | | | | 90 | | | | | 95 | | |
| GCG | AGA | GAC | AAG | GCG | GTT | CGG | GGA | ATT | AGC | AGG | TAC | AAC | TAT | TAC | ATG | 336 |
| Ala | Arg | Asp | Lys | Ala | Val | Arg | Gly | Ile | Ser | Arg | Tyr | Asn | Tyr | Tyr | Met | |
| | | | 100 | | | | | 105 | | | | | 110 | | | |
| GAC | GTC | TGG | GGC | AAA | GGG | ACC | ACG | GTC | ACC | GTC | TCC | TCA | | | | 375 |
| Asp | Val | Typ | Gly | Lys | Gly | Thr | Thr | Val | Thr | Val | Ser | Ser | | | | |
| | | | 115 | | | | 120 | | | | | 125 | | | | |

- (2) INFORMATION FOR SEQ ID NO: 8 (LD1-84-VL)
- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 315 Base Pairs
 - (B) TYPE: Nucleic acid with corresponding protein
 - (C) STRANDEDNESS: Single
 - (D) TOPOLOGY: Linear
- (ii) MOLECULE TYPE: cDNA to mRNA
- (iii) HYPOTHETICAL: No
- (iv) ANTI-SENSE: No
- (vi) ORIGINAL SOURCE: N-terminal fragment
- (A) ORGANISM: Homo sapiens
 - (B) STRAIN: Not applicable
 - (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
 - (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
 - (E) HAPLOTYPE: Diploid
 - (F) TISSUE TYPE: Not applicable
 - (G) CELL TYPE: Peripheral lymphocyte B
 - (H) CELL LINE: Not applicable
 - (I) ORGANELLE: Not applicable
- (vii) IMMEDIATE SOURCE: Escherichia coli
- (A) LIBRARY: cDNA LIBRARY, LD1
 - (B) CLONE: LD1-84
- (viii) POSITION IN GENOME:
- (A) CHROMOSOME / SEGMENT: Chromosome 2
 - (B) MAP POSITION: p12
 - (C) UNITS: Chromosome band number
- (xi) FEATURE:
- (A) NAME / KEY: CDR1, CDR2, CDR3
 - (B) LOCATION: join base pair (64..96, 142..162, 259..285)
 - (C) IDENTIFICATION METHOD: By similarity with a known sequence
 - (D) OTHER INFORMATION:

(2) INFORMATION FOR SEQ ID NO: 9 (LD1-98-VH)

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 375 Base Pairs
- (B) TYPE: Nucleic acid with corresponding protein
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

- (A) ORGANISM: Homo sapiens
- (B) STRAIN: Not applicable
- (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
- (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
- (E) HAPLOTYPE: Diploid
- (F) TISSUE TYPE: Not applicable
- (G) CELL TYPE: Peripheral lymphocyte B
- (H) CELL LINE: Not applicable
- (I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

- (A) LIBRARY: cDNA LIBRARY, LD1
- (B) CLONE: LD1-98

(viii) POSITION IN GENOME:

- (A) CHROMOSOME / SEGMENT: Chromosome 14
- (B) MAP POSITION: q32.3
- (C) UNITS: Chromosome band number

(xi) FEATURE:

- (A) NAME / KEY: CDR1, CDR2, CDR3
- (B) LOCATION: join base pair (91..105, 148..198, 295..342)
- (C) IDENTIFICATION METHOD: By similarity with a known sequence
- (D) OTHER INFORMATION:

(xi) SEQUENCE DESCRIPTION: SEQ ID NO. 9 (LD1-98-VH)

| | |
|---|-----|
| CAG GTG AAA CTG CTC GAG TCT GGG GGA GGC GTG GTC CAG CCG GGG GGG | 48 |
| Gln Val Lys Leu Leu Glu Ser Gly Gly Gly Val Val Gln Pro Gly Gly | |
| 5 10 15 | |
| TCC CTG AGA CTC TCC TGT GAA GCG TCT GGA TTC GCC CTC AGA AGT TCT | 96 |
| Ser Leu Arg Leu Ser Cys Glu Ala Ser Gly Phe Ala Leu Arg Ser Ser | |
| 20 25 30 | |
| GAC ATA CAC TGG GTC CGC CAG GCT CCT GGC AAG GGG CTG GAG TGG GTG | 144 |
| Asp Ile His Typ Val Arg Gln Ala Pro Gly Lys Gly Leu Glu Typ Val | |
| 35 40 45 | |
| GCA CTT ATA TGG TTT GAT GGA AGT ATC AGA TCG TAT GCA GAA TCC GTG | 192 |
| Ala Leu Ile Typ Phe Asp Gly Ser Ile Arg Ser Tyr Ala Glu Ser Val | |
| 50 55 60 | |
| AAG GGC CGA TTC ACC ATC TCC AGA GAC ACT TCC AAG AAC ACC CTA TAT | 240 |
| Lys Gly Arg Phe Thr Ile Ser Arg Asp Thr Ser Lys Asn Thr Leu Tyr | |
| 65 70 75 80 | |
| CTC CAA ATG CGC AGT CTG AGT GCC GAC GAC ACG CGT GTG TAT TAC TGT | 288 |
| Leu Gln Met Arg Ser Leu Ser Ala Asp Asp Thr Ala Val Tyr Tyr Cys | |
| 85 90 95 | |
| GCG AGA GAC AAG GCG GTT CGG GGA ATT AGC AGG TAC AAC TAT TAC ATG | 336 |
| Ala Arg Asp Lys Ala Val Arg Gly Ile Ser Arg Tyr Asn Tyr Tyr Met | |
| 100 105 110 | |
| GAC GTC TGG GGC AAA GGG ACC ACG GTC ACC GTC TCC TCA | 375 |
| Asp Val Typ Gly Lys Gly Thr Thr Val Thr Val Ser Ser | |
| 115 120 125 | |

(2) INFORMATION FOR SEQ ID NO: 10 (LD1-98-VL)

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 318 Base Pairs
- (B) TYPE: Nucleic acid with corresponding protein
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

- (A) ORGANISM: Homo sapiens
- (B) STRAIN: Not applicable
- (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
- (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
- (E) HAPLOTYPE: Diploid
- (F) TISSUE TYPE: Not applicable
- (G) CELL TYPE: Peripheral lymphocyte B
- (H) CELL LINE: Not applicable
- (I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

- (A) LIBRARY: cDNA LIBRARY, LD1
- (B) CLONE: LD1-98

(viii) POSITION IN GENOME:

- (A) CHROMOSOME / SEGMENT: Chromosome 2
- (B) MAP POSITION: p12
- (C) UNITS: Chromosome band number

(xi) FEATURE:

- (A) NAME / KEY: CDR1, CDR2, CDR3
- (B) LOCATION: join base pair (64..96, 142..162, 259..288)
- (C) IDENTIFICATION METHOD: By similarity with a known sequence
- (D) OTHER INFORMATION:

(xi) SEQUENCE DESCRIPTION: SEQ ID NO. 10 (LD1-98-VL)

| | |
|---|-----|
| GTG ATG ACC CAG TCT CCA TCC TCC CTG TCT GCA TCT GTA GGA GAC AGA | 48 |
| Val Met Thr Gln Ser Pro Ser Ser Leu Ser Ala Ser Val Gly Asp Arg | |
| 5 10 15 | |
| GTC ACC ATC ACT TGC CGG GCA AGT CAG AAC ATT ATC CGC TAT TTA AAT | 96 |
| Val Thr Ile Thr Cys Arg Ala Ser Gln Asn Ile Ile Arg Tyr Leu Asn | |
| 20 25 30 | |
| TGG TAT CAG CAG AAG CCA GGG AAA GCC CCT AGG CTC CTG ATC TAT GGT | 144 |
| Typ Tyr Gln Gln Lys Pro Gly Lys Ala Pro Arg Leu Leu Ile Tyr Gly | |
| 35 40 45 | |
| GCG TCC ACT TTG CAA AGT GGG GTC CCA TCA AGG TTC AGT GGC AGT GGA | 192 |
| Ala Ser Thr Leu Gln Ser Gly Val Pro Ser Arg Phe Ser Gly Ser Gly | |
| 50 55 60 | |
| TCT GGG ACA GAT TTC ACT CTC ACC ATC AGT AGT CTG CAA CCT GAA GAT | 240 |
| Ser Gly Thr Asp Phe Thr Leu Thr Ile Ser Ser Leu Gln Pro Glu Asp | |
| 65 70 75 80 | |
| TTT GCA ACT TAC TAC TGT CAA CAG AGT TAC CGT ACC CCT CCA TTC ACT | 288 |
| Phe Ala Thr Tyr Tyr Cys Gln Gln Ser Tyr Arg Thr Pro Pro Phe Thr | |
| 85 90 95 | |
| TTC GGC CCT GGG ACC AAA GTG GAG ATC AAA | 318 |
| Phe Gly Pro Gly Thr Lys Val Glu Ile Lys | |
| 100 105 | |

(2) INFORMATION FOR SEQ ID NO: 11 (LD1-110-VH)

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 375 Base Pairs
- (B) TYPE: Nucleic acid with corresponding protein
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

- (A) ORGANISM: Homo sapiens
- (B) STRAIN: Not applicable
- (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
- (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
- (E) HAPLOTYPE: Diploid
- (F) TISSUE TYPE: Not applicable
- (G) CELL TYPE: Peripheral lymphocyte B
- (H) CELL LINE: Not applicable
- (I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

- (A) LIBRARY: cDNA LIBRARY, LD1
- (B) CLONE: LD1-110

(viii) POSITION IN GENOME:

- (A) CHROMOSOME / SEGMENT: Chromosome 14
- (B) MAP POSITION: q32.3
- (C) UNITS: Chromosome band number

(xi) FEATURE:

- (A) NAME / KEY: CDR1, CDR2, CDR3
- (B) LOCATION: join base pair (91..105, 148..198, 295..342)
- (C) IDENTIFICATION METHOD: By similarity with a known sequence
- (D) OTHER INFORMATION:

(xi) SEQUENCE DESCRIPTION: SEQ ID NO. 11 (LD1-110-VH)

| | |
|---|-----|
| CAG GTG AAA CTG CTC GAG TCT GGG GGA GGC GTG GTC CAG CCT GGG AGG | 48 |
| Gln Val Lys Leu Leu Glu Ser Gly Gly Gly Val Val Gln Pro Gly Arg | |
| 5 10 15 | |
| TCC CTG AGA CTC TCC TGT ATA GCG TCT GGA TTC ACC CTC AGG AAT TAT | 96 |
| Ser Leu Arg Leu Ser Cys Ile Ala Ser Gly Phe Thr Leu Arg Asn Tyr | |
| 20 25 30 | |
| GCC ATG CAC TGG GTC CGC CAG GCT CCA GGA AAG GGG CTG GAG TGG GTG | 144 |
| Ala Met His Typ Val Arg Gln Ala Pro Gly Lys Gly Leu Glu Typ Val | |
| 35 40 45 | |
| GCA GGT ATA TGG TTT GAT GGA AGC AAC AAA AAC TAT GCA GAC TCC GTG | 192 |
| Ala Gly Ile Typ Phe Asp Gly Ser Asn Lys Asn Tyr Ala Asp Ser Val | |
| 50 55 60 | |
| AAG GGC CGA TTC ACC ATC TCC AGA GAC AAC TCC AAG AAC ACT CTG TTT | 240 |
| Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ser Lys Asn Thr Leu Phe | |
| 65 70 75 80 | |
| CTG CAC ATG AAC AGC CTG AGA GCC GAG GAC ACG GCT ACA TAT TAC TGT | 288 |
| Leu His Met Asn Ser Leu Arg Ala Glu Asp Thr Ala Thr Tyr Tyr Cys | |
| 85 90 95 | |
| GCG AGA GAG AGG GCG ATT CGG GGA ATC AGT AGA TAC AAT TAC TAC ATG | 336 |
| Ala Arg Glu Arg Ala Ile Arg Gly Ile Ser Arg Tyr Asn Tyr Tyr Met | |
| 100 105 110 | |
| GAC GTC TGG GGC AAG GGG ACC ACG GTC ACC GTC TCC TCA | 375 |
| Asp Val Typ Gly Lys Gly Thr Thr Val Thr Val Ser Ser | |
| 115 120 125 | |

- (2) INFORMATION FOR SEQ ID NO: 12 (LD1-110-VL)
- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 315 Base Pairs
 - (B) TYPE: Nucleic acid with corresponding protein
 - (C) STRANDEDNESS: Single
 - (D) TOPOLOGY: Linear
- (ii) MOLECULE TYPE: cDNA to mRNA
- (iii) HYPOTHETICAL: No
- (iv) ANTI-SENSE: No
- (vi) ORIGINAL SOURCE: N-terminal fragment
- (A) ORGANISM: Homo sapiens
 - (B) STRAIN: Not applicable
 - (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
 - (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
 - (E) HAPLOTYPE: Diploid
 - (F) TISSUE TYPE: Not applicable
 - (G) CELL TYPE: Peripheral lymphocyte B
 - (H) CELL LINE: Not applicable
 - (I) ORGANELLE: Not applicable
- (vii) IMMEDIATE SOURCE: Escherichia coli
- (A) LIBRARY: cDNA LIBRARY, LD1
 - (B) CLONE: LD1-110
- (viii) POSITION IN GENOME:
- (A) CHROMOSOME / SEGMENT: Chromosome 2
 - (B) MAP POSITION: p12
 - (C) UNITS: Chromosome band number
- (xi) FEATURE:
- (A) NAME / KEY: CDR1, CDR2, CDR3
 - (B) LOCATION: join base pair (64..96, 142..162, 259..285)
 - (C) IDENTIFICATION METHOD: By similarity with a known sequence
 - (D) OTHER INFORMATION:

| | | | | | | | | | | |
|-----|-----|-----|-----|-----|-----|-----|-----|-----|--|-----|
| GGC | CAA | GGG | ACC | AAG | GTG | GAA | ATC | AAA | | 315 |
| Gly | Gln | Gly | Thr | Lys | Val | Glu | Ile | Lys | | |
| | | | 100 | | | | | 105 | | |

- (2) INFORMATION FOR SEQ ID NO: 13 (LD1-117-VH)
- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 378 Base Pairs
 - (B) TYPE: Nucleic acid with corresponding protein
 - (C) STRANDEDNESS: Single
 - (D) TOPOLOGY: Linear
- (ii) MOLECULE TYPE: cDNA to mRNA
- (iii) HYPOTHETICAL: No
- (iv) ANTI-SENSE: No
- (vi) ORIGINAL SOURCE: N-terminal fragment
- (A) ORGANISM: Homo sapiens
 - (B) STRAIN: Not applicable
 - (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
 - (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
 - (E) HAPLOTYPE: Diploid
 - (F) TISSUE TYPE: Not applicable
 - (G) CELL TYPE: Peripheral lymphocyte B
 - (H) CELL LINE: Not applicable
 - (I) ORGANELLE: Not applicable
- (vii) IMMEDIATE SOURCE: Escherichia coli
- (A) LIBRARY: cDNA LIBRARY, LD1
 - (B) CLONE: LD1-117
- (viii) POSITION IN GENOME:
- (A) CHROMOSOME / SEGMENT: Chromosome 14
 - (B) MAP POSITION: q32.3
 - (C) UNITS: Chromosome band number
- (xi) FEATURE:
- (A) NAME / KEY: CDR1, CDR2, CDR3
 - (B) LOCATION: join base pair (91..105, 148..198, 295..345)
 - (C) IDENTIFICATION METHOD: By similarity with a known sequence
 - (D) OTHER INFORMATION:

(xi) SEQUENCE DESCRIPTION: SEQ ID NO. 13 (LD1-117-VH)

| | |
|---|-----|
| CAG GTG AAA CTG CTC GAG TCA GGA GGA GGC GTG GTC CAG CCT GGG AAG | 48 |
| Gln Val Lys Leu Leu Glu Ser Gly Gly Gly Val Val Gln Pro Gly Lys | |
| 5 10 15 | |
| TCC CTG AGA CTT TCC TGT GCA GCG TCT GGA TTC AGT TTC AAT AGC CAT | 96 |
| Ser Leu Arg Leu Ser Cys Ala Ala Ser Gly Phe Ser Phe Asn Ser His | |
| 20 25 30 | |
| GGN ATG CAC TGG GTC CGC CAG GCT CCA GGC AAG GGG CTG GAG TGG GTG | 144 |
| Gly Met His Typ Val Arg Gln Ala Pro Gly Lys Gly Leu Glu Typ Val | |
| 35 40 45 | |
| GCA TTT ATT TGG TTT GAT GGC AGT AAT AAA TAC TAT GCA GAC TCC GTG | 192 |
| Ala Phe Ile Typ Phe Asp Gly Ser Asn Lys Tyr Tyr Ala Asp Ser Val | |
| 50 55 60 | |
| AAG GGC CGT TTC ACC ATC ACC AGA GAC AAC TCC AAG AAC ACG CTG TAT | 240 |
| Lys Gly Arg Phe Thr Ile Thr Arg Asp Asn Ser Lys Asn Thr Leu Tyr | |
| 65 70 75 80 | |
| CTN CAA ATG AAC AGC CTG AGA GCC GAG GAC ACG GCT GTC TAT TAC TGT | 288 |
| Leu Gln Met Asn Ser Leu Arg Ala Glu Asp Thr Ala Val Tyr Tyr Cys | |
| 85 90 95 | |
| GCG AGA GAG ACC TCA GTA AGG CTA GGG TAT AGC CGC TAC AAT TAC TAC | 336 |
| Ala Arg Glu Thr Ser Val Arg Leu Gly Tyr Ser Arg Tyr Asn Tyr Tyr | |
| 100 105 110 | |
| ATG GAC GTC TGG GCC AAA GGG ACC ACG GTC ACC ATC TCG TCA | 378 |
| Met Asp Val Typ Ala Lys Gly Thr Thr Val Thr Ile Ser Ser | |
| 115 120 125 | |

(2) INFORMATION FOR SEQ ID NO: 14 (LD1-117-VL)

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 318 Base Pairs
- (B) TYPE: Nucleic acid with corresponding protein
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

- (A) ORGANISM: Homo sapiens
- (B) STRAIN: Not applicable
- (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
- (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
- (E) HAPLOTYPE: Diploid
- (F) TISSUE TYPE: Not applicable
- (G) CELL TYPE: Peripheral lymphocyte B
- (H) CELL LINE: Not applicable
- (I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

- (A) LIBRARY: cDNA LIBRARY, LD1
- (B) CLONE: LD1-117

(viii) POSITION IN GENOME:

- (A) CHROMOSOME / SEGMENT: Chromosome 2
- (B) MAP POSITION: p12
- (C) UNITS: Chromosome band number

(xi) FEATURE:

- (A) NAME / KEY: CDR1, CDR2, CDR3
- (B) LOCATION: join base pair (64..96, 142..162, 259..288)
- (C) IDENTIFICATION METHOD: By similarity with a known sequence
- (D) OTHER INFORMATION:

- (2) INFORMATION FOR SEQ ID NO: 15 (LD2-1-VH)
- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 375 Base Pairs
 - (B) TYPE: Nucleic acid with corresponding protein
 - (C) STRANDEDNESS: Single
 - (D) TOPOLOGY: Linear
- (ii) MOLECULE TYPE: cDNA to mRNA
- (iii) HYPOTHETICAL: No
- (iv) ANTI-SENSE: No
- (vi) ORIGINAL SOURCE: N-terminal fragment
- (A) ORGANISM: Homo sapiens
 - (B) STRAIN: Not applicable
 - (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
 - (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
 - (E) HAPLOTYPE: Diploid
 - (F) TISSUE TYPE: Not applicable
 - (G) CELL TYPE: Peripheral lymphocyte B
 - (H) CELL LINE: Not applicable
 - (I) ORGANELLE: Not applicable
- (vii) IMMEDIATE SOURCE: Escherichia coli
- (A) LIBRARY: cDNA LIBRARY, LD2
 - (B) CLONE: LD2-1
- (viii) POSITION IN GENOME:
- (A) CHROMOSOME / SEGMENT: Chromosome 14
 - (B) MAP POSITION: q32.3
 - (C) UNITS: Chromosome band number
- (xi) FEATURE:
- (A) NAME / KEY: CDR1, CDR2, CDR3
 - (B) LOCATION: join base pair (91..105, 148..198, 295..342)
 - (C) IDENTIFICATION METHOD: By similarity with a known sequence
 - (D) OTHER INFORMATION:

(xi) SEQUENCE DESCRIPTION: SEQ ID NO. 15 (LD2-1-VH)

| | |
|---|-----|
| CAG GTG AAA CTG CTC GAG TCT GGG GGA GGC GTG GTC CAG CCG GGG GGG | 48 |
| Gln Val Lys Leu Leu Glu Ser Gly Gly Val Val Gln Pro Gly Gly | |
| 5 10 15 | |
| TCC CTG AGA CTC TCC TGT GTA GCG TCT GGA TTC ACC CTC AGG AGT TAT | 96 |
| Ser Leu Arg Leu Ser Cys Val Ala Ser Gly Phe Thr Leu Arg Ser Tyr | |
| 20 25 30 | |
| GGC ATG CAC TGG GTC CGC CAG GCT CCA GGA AAG GGC CTG GAG TGG GTG | 144 |
| Gly Met His Typ Val Arg Gln Ala Pro Gly Lys Gly Leu Glu Typ Val | |
| 35 40 45 | |
| GCT TTT ATA TGG TTT GAT GGA AGT AAT AAA GGA TAT GTA GAC TCC GTG | 192 |
| Ala Phe Ile Typ Phe Asp Gly Ser Asn Lys Gly Tyr Val Asp Ser Val | |
| 50 55 60 | |
| AAG GGC CGA TTC ACC ATC TCC CGA GAC AAT TCC AAG AAC ATG GTC TAT | 240 |
| Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ser Lys Asn Met Val Tyr | |
| 65 70 75 80 | |
| CTC CAA ATG AAC AGC CTG AGA GCC GAT GAC ACG GCT GTA TAT TAT TGT | 288 |
| Leu Gln Met Asn Ser Leu Arg Ala Asp Asp Thr Ala Val Tyr Tyr Cys | |
| 85 90 95 | |
| GCG AGA GAG AAG GCG CTT CGG GGA ATC AGC AGA TAC AAC TAT TAC CTG | 336 |
| Ala Arg Glu Lys Ala Leu Arg Gly Ile Ser Arg Tyr Asn Tyr Tyr Leu | |
| 100 105 110 | |
| GAC GTC TGG GGC AAG GGG ACC ACG GTC ACC GTC TCC TCA | 375 |
| Asp Val Typ Gly Lys Gly Thr Thr Val Thr Val Ser Ser | |
| 115 120 125 | |

(2) INFORMATION FOR SEQ ID NO: 16 (LD2-1-VL)

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 333 Base Pairs

(B) TYPE: Nucleic acid with corresponding protein

(C) STRANDEDNESS: Single

(D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

(A) ORGANISM: Homo sapiens

(B) STRAIN: Not applicable

(C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor

(D) DEVELOPMENT STAGE: Adult, Rearranged development pattern

(E) HAPLOTYPE: Diploid

(F) TISSUE TYPE: Not applicable

(G) CELL TYPE: Peripheral lymphocyte B

(H) CELL LINE: Not applicable

(I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

(A) LIBRARY: cDNA LIBRARY, LD2

(B) CLONE: LD2-1

(viii) POSITION IN GENOME:

(A) CHROMOSOME / SEGMENT: Chromosome 22

(B) MAP POSITION: q11.2

(C) UNITS: Chromosome band number

(xi) FEATURE:

(A) NAME / KEY: CDR1, CDR2, CDR3

(B) LOCATION: join base pair (61..99, 145..165 , 262..294)

(C) IDENTIFICATION METHOD: By similarity with a known sequence

(D) OTHER INFORMATION:

(xi) SEQUENCE DESCRIPTION: SEQ ID NO. 16 (LD2-1-VL)

| | | | | | | | | | | | | | | | | |
|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| GTG | GTG | ACT | CAG | CCA | CCC | TCA | GCG | TCT | GGG | ACC | CCC | GGA | CAG | AGG | GTC | 48 |
| Val | Val | Thr | Gln | Pro | Pro | Ser | Ala | Ser | Gly | Thr | Pro | Gly | Gln | Arg | Val | |
| | | | | 5 | | | | | 10 | | | | | 15 | | |
| ACC | ATC | TCT | TGT | TCT | GGA | AGC | AAC | TCC | ATC | CTT | GGA | AGT | AAG | TAT | GTA | 96 |
| Thr | Ile | Ser | Cys | Ser | Gly | Ser | Asn | Ser | Ile | Leu | Gly | Ser | Lys | Tyr | Val | |
| | | | 20 | | | | | 25 | | | | | 30 | | | |
| TAC | TGG | TAC | CAG | AAA | CTC | CCA | GGA | ACG | GCC | CCC | AAA | CTC | CTC | ATC | TAT | 144 |
| Tyr | Typ | Tyr | Gln | Lys | Leu | Pro | Gly | Thr | Ala | Pro | Lys | Leu | Leu | Ile | Tyr | |
| | | | 35 | | | | 40 | | | | | 45 | | | | |
| AAG | AAT | GAT | CAG | CGG | CCC | TCA | GGG | GTC | TCT | GAC | CGA | TTC | TCT | GGC | TCC | 192 |
| Lys | Asn | Asp | Gln | Arg | Pro | Ser | Gly | Val | Ser | Asp | Arg | Phe | Ser | Gly | Ser | |
| | 50 | | | | | 55 | | | | | 60 | | | | | |
| AAG | TCT | GGC | ACC | TCG | GCC | TCC | CTG | GCC | ATC | AGT | GGG | CTC | CGG | TCC | GAG | 240 |
| Lys | Ser | Gly | Thr | Ser | Ala | Ser | Leu | Ala | Ile | Ser | Gly | Leu | Arg | Ser | Glu | |
| 65 | | | | | 70 | | | | 75 | | | | | 80 | | |
| GAT | GAG | GCT | GAC | TAT | TAC | TGT | GCA | CCA | TGG | GAT | GCC | AAC | CTG | GGT | GGC | 288 |
| Asp | Glu | Ala | Asp | Tyr | Tyr | Cys | Ala | Pro | Typ | Asp | Ala | Asn | Leu | Gly | Gly | |
| | | | | 85 | | | | | 90 | | | | | 95 | | |
| CCG | GTG | TTC | GGC | GGA | GGG | ACC | AAG | CTG | ACC | GTC | CTA | AGT | CAG | CCC | | 333 |
| Pro | Val | Phe | Gly | Gly | Gly | Thr | Lys | Leu | Thr | Val | Leu | Ser | Gln | Pro | | |
| | | | 100 | | | | | 105 | | | | | 110 | | | |

(2) INFORMATION FOR SEQ ID NO: 17 (LD2-4-VH)

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 375 Base Pairs
- (B) TYPE: Nucleic acid with corresponding protein
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

- (A) ORGANISM: Homo sapiens
- (B) STRAIN: Not applicable
- (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
- (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
- (E) HAPLOTYPE: Diploid
- (F) TISSUE TYPE: Not applicable
- (G) CELL TYPE: Peripheral lymphocyte B
- (H) CELL LINE: Not applicable
- (I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

- (A) LIBRARY: cDNA LIBRARY, LD2
- (B) CLONE: LD2-4

(viii) POSITION IN GENOME:

- (A) CHROMOSOME / SEGMENT: Chromosome 14
- (B) MAP POSITION: q32.3
- (C) UNITS: Chromosome band number

(xi) FEATURE:

- (A) NAME / KEY: CDR1, CDR2, CDR3
- (B) LOCATION: join base pair (91..105, 148..198, 295..342)
- (C) IDENTIFICATION METHOD: By similarity with a known sequence
- (D) OTHER INFORMATION:

(xi) SEQUENCE DESCRIPTION: SEQ ID NO. 17 (LD2-4-VH)

| | |
|---|-----|
| CAG GTG AAA CTG CTC GAG TCG GGG GGA GGC GTG GTC CAG CCG GGG GGG | 48 |
| Gln Val Lys Leu Leu Glu Ser Gly Gly Gly Val Val Gln Pro Gly Gly | |
| 5 10 15 | |
| TCC CTG AGA CTC TCC TGT GAA GCG TCT GGA TTC ACC CTC AGA AGT TCT | 96 |
| Ser Leu Arg Leu Ser Cys Glu Ala Ser Gly Phe Thr Leu Arg Ser Ser | |
| 20 25 30 | |
| GGC ATG CAC TGG GTC CGC CAG GCT CCT GGC AAG GGG CTG GAG TGG GTG | 144 |
| Gly Met His Typ Val Arg Gln Ala Pro Gly Lys Gly Leu Glu Typ Val | |
| 35 40 45 | |
| GCA CTT ATA TGG TTT GAT GGA AGT ATC AGA TCG TAT GCA GAA TCC GTG | 192 |
| Ala Leu Ile Typ Phe Asp Gly Ser Ile Arg Ser Tyr Ala Glu Ser Val | |
| 50 55 60 | |
| AAG GGC CGA TTC ACC ATC TCC AGA GAC ACT TCC AAG AAC ACC CTA TAT | 240 |
| Lys Gly Arg Phe Thr Ile Ser Arg Asp Thr Ser Lys Asn Thr Leu Tyr | |
| 65 70 75 80 | |
| CTC CAA ATG CGC AGT CTG AGT GCC GAC GAC ACG GCT GTG TAT TAC TGT | 288 |
| Leu Gln Met Arg Ser Leu Ser Ala Asp Asp Thr Ala Val Tyr Tyr Cys | |
| 85 90 95 | |
| GCG AGA GAC AAG GCG GTT CGG GGA ATT AGC AGG TAC AAC TAT TAC ATG | 336 |
| Ala Arg Asp Lys Ala Val Arg Gly Ile Ser Arg Tyr Asn Tyr Tyr Met | |
| 100 105 110 | |
| GAC GTC TGG GGC AAA GGG ACC ACG GTC ACC GTC TCC TCA | 375 |
| Asp Val Typ Gly Lys Gly Thr Thr Val Thr Val Ser Ser | |
| 115 120 125 | |

(2) INFORMATION FOR SEQ ID NO: 18 (LD2-4-VL)

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 312 Base Pairs

(B) TYPE: Nucleic acid with corresponding protein

(C) STRANDEDNESS: Single

(D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

(A) ORGANISM: Homo sapiens

(B) STRAIN: Not applicable

(C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor

(D) DEVELOPMENT STAGE: Adult, Rearranged development pattern

(E) HAPLOTYPE: Diploid

(F) TISSUE TYPE: Not applicable

(G) CELL TYPE: Peripheral lymphocyte B

(H) CELL LINE: Not applicable

(I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

(A) LIBRARY: cDNA LIBRARY, LD2

(B) CLONE: LD2-4

(viii) POSITION IN GENOME:

(A) CHROMOSOME / SEGMENT: Chromosome 2

(B) MAP POSITION: p12

(C) UNITS: Chromosome band number

(xi) FEATURE:

(A) NAME / KEY: CDR1, CDR2, CDR3

(B) LOCATION: join base pair (64..96, 142..162, 259..282)

(C) IDENTIFICATION METHOD: By similarity with a known sequence

(D) OTHER INFORMATION:

(2) INFORMATION FOR SEQ ID NO: 19 (LD2-5-VH)

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 375 Base Pairs
- (B) TYPE: Nucleic acid with corresponding protein
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

- (A) ORGANISM: Homo sapiens
- (B) STRAIN: Not applicable
- (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
- (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
- (E) HAPLOTYPE: Diploid
- (F) TISSUE TYPE: Not applicable
- (G) CELL TYPE: Peripheral lymphocyte B
- (H) CELL LINE: Not applicable
- (I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

- (A) LIBRARY: cDNA LIBRARY, LD2
- (B) CLONE: LD2-5

(viii) POSITION IN GENOME:

- (A) CHROMOSOME / SEGMENT: Chromosome 14
- (B) MAP POSITION: q32.3
- (C) UNITS: Chromosome band number

(xi) FEATURE:

- (A) NAME / KEY: CDR1, CDR2, CDR3
- (B) LOCATION: join base pair (91..105, 148..198, 295..342)
- (C) IDENTIFICATION METHOD: By similarity with a known sequence
- (D) OTHER INFORMATION:

(xi) SEQUENCE DESCRIPTION: SEQ ID NO. 19 (LD2-5-VH)

| | |
|---|-----|
| CAG GTG AAA CTG CTC GAG TCT GGG GGA GGC TTG GTC CAG CCG GGG GGG | 48 |
| Gln Val Lys Leu Leu Glu Ser Gly Gly Leu Val Gln Pro Gly Gly | |
| 5 10 15 | |
| TCC CTG AGA CTC TCC TGT GTA GCG TCT GGA TTC ACC TTC AGG AGT TAT | 96 |
| Ser Leu Arg Leu Ser Cys Val Ala Ser Gly Phe Thr Phe Arg Ser Tyr | |
| 20 25 30 | |
| GGA ATG CAC TGG GTC CGC CAG GCT CCA GGA AAG GGC CTG GAG TGG GTG | 144 |
| Gly Met His Typ Val Arg Gln Ala Pro Gly Lys Gly Leu Glu Typ Val | |
| 35 40 45 | |
| GCT TTT ATA TGG TTT GAT GGA AGT AAT AAA GGA TAT GTA GAC TCC GTG | 192 |
| Ala Phe Ile Typ Phe Asp Gly Ser Asn Lys Gly Tyr Val Asp Ser Val | |
| 50 55 60 | |
| AAG GGC CGA TTC ACC ATC TCC CGA GAC AAT TCC AAG AAC ATG CTC TAT | 240 |
| Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ser Lys Asn Met Leu Tyr | |
| 65 70 75 80 | |
| CTG CAA ATG AAT AGC CTG AGA GCC GAG GAC ACG GCT GTA TAT TAT TGT | 288 |
| Leu Gln Met Asn Ser Leu Arg Ala Glu Asp Thr Ala Val Tyr Tyr Cys | |
| 85 90 95 | |
| GCG AGA GAG AAG GCG CTT CGG GGA ATC AGT AGA TAC AAC TAT TAC CTG | 336 |
| Ala Arg Glu Lys Ala Leu Arg Gly Ile Ser Arg Tyr Asn Tyr Tyr Leu | |
| 100 105 110 | |
| GAC GTC TGG GGC AAG GGG GCC ACG GTC ACC GTC TCC TCA | 375 |
| Asp Val Typ Gly Lys Gly Ala Thr Val Thr Val Ser Ser | |
| 115 120 125 | |

22) INFORMATION FOR SEQ ID NO: 20 (LD2-5-VL)

- (i) SEQUENCE CHARACTERISTICS:
 - (A) LENGTH: 318 Base Pairs
 - (B) TYPE: Nucleic acid with corresponding protein
 - (C) STRANDEDNESS: Single
 - (D) TOPOLOGY: Linear
- (ii) MOLECULE TYPE: cDNA to mRNA
- (iii) HYPOTHETICAL: No
- (iv) ANTI-SENSE: No
- (vi) ORIGINAL SOURCE: N-terminal fragment
 - (A) ORGANISM: Homo sapiens
 - (B) STRAIN: Not applicable
 - (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
 - (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
 - (E) HAPLOTYPE: Diploid
 - (F) TISSUE TYPE: Not applicable
 - (G) CELL TYPE: Peripheral lymphocyte B
 - (H) CELL LINE: Not applicable
 - (I) ORGANELLE: Not applicable
- (vii) IMMEDIATE SOURCE: Escherichia coli
 - (A) LIBRARY: cDNA LIBRARY, LD2
 - (B) CLONE: LD2-5
- (viii) POSITION IN GENOME:
 - (A) CHROMOSOME / SEGMENT: Chromosome 2
 - (B) MAP POSITION: p12
 - (C) UNITS: Chromosome band number
- (xi) FEATURE:
 - (A) NAME / KEY: CDR1, CDR2, CDR3
 - (B) LOCATION: join base pair (64..96, 142..162, 259..288)
 - (C) IDENTIFICATION METHOD: By similarity with a known sequence
 - (D) OTHER INFORMATION:

(2) INFORMATION FOR SEQ ID NO: 21 (LD2-10-VH)

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 378 Base Pairs
- (B) TYPE: Nucleic acid with corresponding protein
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

- (A) ORGANISM: Homo sapiens
- (B) STRAIN: Not applicable
- (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
- (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
- (E) HAPLOTYPE: Diploid
- (F) TISSUE TYPE: Not applicable
- (G) CELL TYPE: Peripheral lymphocyte B
- (H) CELL LINE: Not applicable
- (I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

- (A) LIBRARY: cDNA LIBRARY, LD2
- (B) CLONE: LD2-10

(viii) POSITION IN GENOME:

- (A) CHROMOSOME / SEGMENT: Chromosome 14
- (B) MAP POSITION: q32.3
- (C) UNITS: Chromosome band number

(xi) FEATURE:

- (A) NAME / KEY: CDR1, CDR2, CDR3
- (B) LOCATION: join base pair (91..105, 148..198, 298..345)
- (C) IDENTIFICATION METHOD: By similarity with a known sequence
- (D) OTHER INFORMATION:

(xi) SEQUENCE DESCRIPTION: SEQ ID NO. 21 (LD2-10-VH)

| | |
|---|-----|
| CAG GTG AAA CTG CTC GAG TCT GGG GGA GGC GTG GTC CAG CCG GGG GGG | 48 |
| Gln Val Lys Leu Leu Glu Ser Gly Gly Gly Val Val Gln Pro Gly Gly | |
| 5 10 15 | |
| TCC CTG AGA CTC TCC TGT GTA GCG TCT GGA TTC ACC CTC AGG AGT TAT | 96 |
| Ser Leu Arg Leu Ser Cys Val Ala Ser Gly Phe Thr Leu Arg Ser Tyr | |
| 20 25 30 | |
| GGC ATG CAC TGG GTC CGC CAG GCT CCA GGC AAG GGC CTG GAG TGG GTG | 144 |
| Gly Met His Typ Val Arg Gln Ala Pro Gly Lys Gly Leu Glu Typ Val | |
| 35 40 45 | |
| GCT TTT ATA TGG TTT GAT GGA AGT AAT AAA GGA TAT GTA GAC TCC GTG | 192 |
| Ala Phe Ile Typ Phe Asp Gly Ser Asn Lys Gly Tyr Val Asp Ser Val | |
| 50 55 60 | |
| AAG GGC CGA TTC ACC ATC TCC CGA GAC AAT TCC AAG AAC ATG GTC TAT | 240 |
| Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ser Lys Asn Met Val Tyr | |
| 65 70 75 80 | |
| CTG CAA ATG AAC AGC CTG AGA GCC GAT GAC ACG GCT GTA TAT TAT TAT | 288 |
| Leu Gln Met Asn Ser Leu Arg Ala Asp Asp Thr Ala Val Tyr Tyr Tyr | |
| 85 90 95 | |
| TGT GCG AGA GAG AAG GCG CTT CGG GGA ATC AGC AGA TAC AAC TAT TAC | 336 |
| Cys Ala Arg Glu Lys Ala Leu Arg Gly Ile Ser Arg Tyr Asn Tyr Tyr | |
| 100 105 110 | |
| CTG GAC GTC TGG GGC AAG GGG ACC ACG GTC ACC GTC TCC TCA | 378 |
| Leu Asp Val Typ Gly Lys Gly Thr Thr Val Thr Val Ser Ser | |
| 115 120 125 | |

(2) INFORMATION FOR SEQ ID NO: 22 (LD2-10-VL)

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 333 Base Pairs

(B) TYPE: Nucleic acid with corresponding protein

(C) STRANDEDNESS: Single

(D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

(A) ORGANISM: Homo sapiens

(B) STRAIN: Not applicable

(C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor

(D) DEVELOPMENT STAGE: Adult, Rearranged development pattern

(E) HAPLOTYPE: Diploid

(F) TISSUE TYPE: Not applicable

(G) CELL TYPE: Peripheral lymphocyte B

(H) CELL LINE: Not applicable

(I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

(A) LIBRARY: cDNA LIBRARY, LD2

(B) CLONE: LD2-10

(viii) POSITION IN GENOME:

(A) CHROMOSOME / SEGMENT: Chromosome 22

(B) MAP POSITION: q11.2

(C) UNITS: Chromosome band number

(xi) FEATURE:

(A) NAME / KEY: CDR1, CDR2, CDR3

(B) LOCATION: join base pair (61..102, 148..168, 265..294)

(C) IDENTIFICATION METHOD: By similarity with a known sequence

(D) OTHER INFORMATION:

(xi) SEQUENCE DESCRIPTION: SEQ ID NO. 22 (LD2-10-VL)

| | | | | | | | | | | | | | | | | |
|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| GTG | GTG | ACT | CAG | GAG | CCC | TCA | CTG | ACT | GTG | TCC | CCA | GGA | GGG | ACA | GTC | 48 |
| Val | Val | Thr | Gln | Glu | Pro | Ser | Leu | Thr | Val | Ser | Pro | Gly | Gly | Thr | Val | |
| | | | | 5 | | | | | 10 | | | | | 15 | | |
| ACT | CTC | ACC | TGT | GCT | TCC | AGC | ACT | GGG | GCA | GTC | ACC | AGG | GGT | TAC | TAT | 96 |
| Thr | Leu | Thr | Cys | Ala | Ser | Ser | Thr | Gly | Ala | Val | Thr | Arg | Gly | Tyr | Tyr | |
| | | | 20 | | | | | 25 | | | | | 30 | | | |
| CCA | AAC | TGG | TTC | CAG | CAG | AAG | CCT | GGA | CAA | GCA | CCC | AGG | GCA | CTG | ATT | 144 |
| Pro | Asn | Typ | Phe | Gln | Gln | Lys | Pro | Gly | Gln | Ala | Pro | Arg | Ala | Leu | Ile | |
| | | 35 | | | | | 40 | | | | | 45 | | | | |
| TAT | AGT | ACA | AAC | AAA | AAA | CAC | TCC | TGG | ACC | CCT | GCC | CGG | TTC | TCA | GGC | 192 |
| Tyr | Ser | Thr | Asn | Lys | Lys | His | Ser | Typ | Thr | Pro | Ala | Arg | Phe | Ser | Gly | |
| | 50 | | | | | 55 | | | | 60 | | | | | | |
| TCC | CTC | CTT | GGG | GGC | AAA | GCT | GCC | CTG | ACA | CTG | TCA | GGT | GTG | CAG | CCT | 240 |
| Ser | Leu | Leu | Gly | Gly | Lys | Ala | Ala | Leu | Thr | Leu | Ser | Gly | Val | Gln | Pro | |
| 65 | | | | | 70 | | | | 75 | | | | | 80 | | |
| GAA | GAC | GAG | GCT | GAA | TAT | TAC | TGC | CTG | CTC | TAC | TAT | GGT | GGT | GCT | CAA | 288 |
| Glu | Asp | Glu | Ala | Glu | Tyr | Tyr | Cys | Leu | Leu | Tyr | Tyr | Gly | Gly | Ala | Gln | |
| | | | 85 | | | | | 90 | | | | | | 95 | | |
| CTC | GTA | TTC | GGC | GGA | GGG | ACC | AAG | CTG | ACC | GTC | CTA | CGT | CAG | CCC | | 333 |
| Leu | Val | Phe | Gly | Gly | Gly | Thr | Lys | Leu | Thr | Val | Leu | Arg | Gln | Pro | | |
| | | | 100 | | | | | 105 | | | | | 110 | | | |

(2) INFORMATION FOR SEQ ID NO: 23 (LD2-11-VH)

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 375 Base Pairs

(B) TYPE: Nucleic acid with corresponding protein

(C) STRANDEDNESS: Single

(D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

(A) ORGANISM: Homo sapiens

(B) STRAIN: Not applicable

(C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor

(D) DEVELOPMENT STAGE: Adult, Rearranged development pattern

(E) HAPLOTYPE: Diploid

(F) TISSUE TYPE: Not applicable

(G) CELL TYPE: Peripheral lymphocyte B

(H) CELL LINE: Not applicable

(I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

(A) LIBRARY: cDNA LIBRARY, LD2

(B) CLONE: LD2-11

(viii) POSITION IN GENOME:

(A) CHROMOSOME / SEGMENT: Chromosome 14

(B) MAP POSITION: q32.3

(C) UNITS: Chromosome band number

(xi) FEATURE:

(A) NAME / KEY: CDR1, CDR2, CDR3

(B) LOCATION: join base pair (91..105, 148..198, 295..342)

(C) IDENTIFICATION METHOD: By similarity with a known sequence

(D) OTHER INFORMATION:

(xi) SEQUENCE DESCRIPTION: SEQ ID NO. 23 (LD2-11-VH)

| | |
|---|-----|
| CAG GTG AAA CTG CTC GAG TCG GGG GGA GGC GTG GTC CAG CCG GGG GGG | 48 |
| Gln Val Lys Leu Leu Glu Ser Gly Gly Gly Val Val Gln Pro Gly Gly | |
| 5 10 15 | |
| TCC CTG AGA CTC TCC TGT GAA GCG TCT GGA TTC ACC CTC AGA AGT TCT | 96 |
| Ser Leu Arg Leu Ser Cys Glu Ala Ser Gly Phe Thr Leu Arg Ser Ser | |
| 20 25 30 | |
| GGC ATG CAC TGG GTC CGC CAG GCT CCT GGC AAG GGG CTG GAG TGG GTG | 144 |
| Gly Met His Typ Val Arg Gln Ala Pro Gly Lys Gly Leu Glu Typ Val | |
| 35 40 45 | |
| GCA CTT ATA TGG TTT GAT GGA AGT ATC AGA TCG TAT GCA GAA TCC GTG | 192 |
| Ala Leu Ile Typ Phe Asp Gly Ser Ile Arg Ser Tyr Ala Glu Ser Val | |
| 50 55 60 | |
| AAG GGC CGA TTC ACC ATC TCC AGA GAC ACT TCC AAG AAC ACC CTA TAT | 240 |
| Lys Gly Arg Phe Thr Ile Ser Arg Asp Thr Ser Lys Asn Thr Leu Tyr | |
| 65 70 75 80 | |
| CTC CAA ATG CGC AGT CTG AGT GCC GAC GAC ACG GCT GTG TAT TAC TGT | 288 |
| Leu Gln Met Arg Ser Leu Ser Ala Asp Asp Thr Ala Val Tyr Tyr Cys | |
| 85 90 95 | |
| GCG AGA GAC AAG GCG GTT CGG GGA ATT AGC AGG TAC AAC TAT TAC ATG | 336 |
| Ala Arg Asp Lys Ala Val Arg Gly Ile Ser Arg Tyr Asn Tyr Tyr Met | |
| 100 105 110 | |
| GAC GTC TGG GGC AAA GGG ACC ACG GTC ACC GTC TCC TCA | 375 |
| Asp Val Typ Gly Lys Gly Thr Thr Val Thr Val Ser Ser | |
| 115 120 125 | |

- (2) INFORMATION FOR SEQ ID NO: 24 (LD2-11-VL)
- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 315 Base Pairs
 - (B) TYPE: Nucleic acid with corresponding protein
 - (C) STRANDEDNESS: Single
 - (D) TOPOLOGY: Linear
- (ii) MOLECULE TYPE: cDNA to mRNA
- (iii) HYPOTHETICAL: No
- (iv) ANTI-SENSE: No
- (vi) ORIGINAL SOURCE: N-terminal fragment
- (A) ORGANISM: Homo sapiens
 - (B) STRAIN: Not applicable
 - (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
 - (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
 - (E) HAPLOTYPE: Diploid
 - (F) TISSUE TYPE: Not applicable
 - (G) CELL TYPE: Peripheral lymphocyte B
 - (H) CELL LINE: Not applicable
 - (I) ORGANELLE: Not applicable
- (vii) IMMEDIATE SOURCE: Escherichia coli
- (A) LIBRARY: cDNA LIBRARY, LD2
 - (B) CLONE: LD2-11
- (viii) POSITION IN GENOME:
- (A) CHROMOSOME / SEGMENT: CHROMOSOME 2
 - (B) MAP POSITION: p12
 - (C) UNITS: Chromosome band number
- (xi) FEATURE:
- (A) NAME / KEY: CDR1, CDR2, CDR3
 - (B) LOCATION: join base pair (64..96, 142..162, 259..285)
 - (C) IDENTIFICATION METHOD: By similarity with a known sequence
 - (D) OTHER INFORMATION:

(2) INFORMATION FOR SEQ ID NO: 25 (LD2-14-VH)

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 375 Base Pairs
- (B) TYPE: Nucleic acid with corresponding protein
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

- (A) ORGANISM: Homo sapiens
- (B) STRAIN: Not applicable
- (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
- (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
- (E) HAPLOTYPE: Diploid
- (F) TISSUE TYPE: Not applicable
- (G) CELL TYPE: Peripheral lymphocyte B
- (H) CELL LINE: Not applicable
- (I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

- (A) LIBRARY: cDNA LIBRARY, LD2
- (B) CLONE: LD2-14

(viii) POSITION IN GENOME:

- (A) CHROMOSOME / SEGMENT: CHROMOSOME 14
- (B) MAP POSITION: q32.3
- (C) UNITS: Chromosome band number

(xi) FEATURE:

- (A) NAME / KEY: CDR1, CDR2, CDR3
- (B) LOCATION: join base pair (91..105, 148..198, 295..342)
- (C) IDENTIFICATION METHOD: By similarity with a known sequence
- (D) OTHER INFORMATION:

(xi) SEQUENCE DESCRIPTION: SEQ ID NO. 25 (LD2-14-VH)

| | |
|---|-----|
| CAG GTG AAA CTG CTC GAG TCT GGG GGA GGC GTG GTC CAG CCG GGG GGG | 48 |
| Gln Val Lys Leu Leu Glu Ser Gly Gly Gly Val Val Gln Pro Gly Gly | |
| 5 10 15 | |
| TCC CTG AGA GTC GCC TGT GTA GCG TCT GGA TTC ACC TCC AGG AAT TTT | 96 |
| Ser Leu Arg Val Ala Cys Val Ala Ser Gly Phe Thr Ser Arg Asn Phe | |
| 20 25 30 | |
| GGC ATG CAC TGG GTC CGC CAG GCT CCA GGA AAG GGG CTG GAG TGG GTG | 144 |
| Gly Met His Typ Val Arg Gln Ala Pro Gly Lys Gly Leu Glu Typ Val | |
| 35 40 45 | |
| GTT TTT ATT TGG TTT GAT GCA AGT AAT AAA GGA TAT GGA GAC TCC GTT | 192 |
| Val Phe Ile Typ Phe Asp Ala Ser Asn Lys Gly Tyr Gly Asp Ser Val | |
| 50 55 60 | |
| AAG GGC CGA TTC ACC GTC TCC AGA GAC AAT TCC AAG AAC ACG CTC TAT | 240 |
| Lys Gly Arg Phe Thr Val Ser Arg Asp Asn Ser Lys Asn Thr Leu Tyr | |
| 65 70 75 80 | |
| CTG CAA ATG AAC GGC CTG AGA GCC GAA GAC ACG GCT GTA TAT TAT TGT | 288 |
| Leu Gln Met Asn Gly Leu Arg Ala Glu Asp Thr Ala Val Tyr Tyr Cys | |
| 85 90 95 | |
| GCG AGA GAG AAG GCG GTT CGG GGA ATT AGT AGA TAC AAC TAC TAC ATG | 336 |
| Ala Arg Glu Lys Ala Val Arg Gly Ile Ser Arg Tyr Asn Tyr Tyr Met | |
| 100 105 110 | |
| GAC GTC TGG GGC AAG GGG ACC ACG GTC ACC GTC TCC TCA | 375 |
| Asp Val Typ Gly Lys Gly Thr Thr Val Thr Val Ser Ser | |
| 115 120 125 | |

(2) INFORMATION FOR SEQ ID NO: 26 (LD2-14-VL)

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 315 Base Pairs

(B) TYPE: Nucleic acid with corresponding protein

(C) STRANDEDNESS: Single

(D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

(A) ORGANISM: Homo sapiens

(B) STRAIN: Not applicable

(C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor

(D) DEVELOPMENT STAGE: Adult, Rearranged development pattern

(E) HAPLOTYPE: Diploid

(F) TISSUE TYPE: Not applicable

(G) CELL TYPE: Peripheral lymphocyte B

(H) CELL LINE: Not applicable

(I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

(A) LIBRARY: cDNA LIBRARY, LD2

(B) CLONE: LD2-14

(viii) POSITION IN GENOME:

(A) CHROMOSOME / SEGMENT: CHROMOSOME 2

(B) MAP POSITION: p12

(C) UNITS: Chromosome band number

(xi) FEATURE:

(A) NAME / KEY: CDR1, CDR2, CDR3

(B) LOCATION: join base pair (64..96, 142..162, 259..285)

(C) IDENTIFICATION METHOD: By similarity with a known sequence

(D) OTHER INFORMATION:

(2) INFORMATION FOR SEQ ID NO: 27 (LD2-17-VH)

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 375 Base Pairs
- (B) TYPE: Nucleic acid with corresponding protein
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

- (A) ORGANISM: Homo sapiens
- (B) STRAIN: Not applicable
- (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
- (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
- (E) HAPLOTYPE: Diploid
- (F) TISSUE TYPE: Not applicable
- (G) CELL TYPE: Peripheral lymphocyte B
- (H) CELL LINE: Not applicable
- (I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

- (A) LIBRARY: cDNA LIBRARY, LD2
- (B) CLONE: LD2-17

(viii) POSITION IN GENOME:

- (A) CHROMOSOME / SEGMENT: CHROMOSOME 14
- (B) MAP POSITION: q32.3
- (C) UNITS: Chromosome band number

(xi) FEATURE:

- (A) NAME / KEY: CDR1, CDR2, CDR3
- (B) LOCATION: join base pair (91..105, 148..198, 295..342)
- (C) IDENTIFICATION METHOD: By similarity with a known sequence
- (D) OTHER INFORMATION:

(xi) SEQUENCE DESCRIPTION: SEQ ID NO. 27 (LD2-17-VH)

| | |
|---|-----|
| CAG GTG AAA CTG CTC GAG TCT GGG GGA GGC GTG GTC CAG CCG GGG GGG | 48 |
| Gln Val Lys Leu Leu Glu Ser Gly Gly Val Val Gln Pro Gly Gly | |
| 5 10 15 | |
| TCC CTG AGA CTC TCC TGT GTA GCG TCT GGA TTC ACC TCC AGG AGT TAT | 96 |
| Ser Leu Arg Leu Ser Cys Val Ala Ser Gly Phe Thr Ser Arg Ser Tyr | |
| 20 25 30 | |
| GGA ATG CAC TGG GTC CGC CAG GCT CCA GGC AAG GGC CTG GAG TGG GTG | 144 |
| Gly Met His Typ Val Arg Gln Ala Pro Gly Lys Gly Leu Glu Typ Val | |
| 35 40 45 | |
| GCT TTT ATA TGG TTT GAT GGA AGT AAT AAA GGA TAT GTA GAC TCC GTG | 192 |
| Ala Phe Ile Typ Phe Asp Gly Ser Asn Lys Gly Tyr Val Asp Ser Val | |
| 50 55 60 | |
| AAG GGC CGA TTC ACC ATC TCC CGA GAC AAT TCC AAG AAC ACG CTC TAT | 240 |
| Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ser Lys Asn Thr Leu Tyr | |
| 65 70 75 80 | |
| CTG CAA ATG AAG AGC CTG AGA GCC GAG GAC ACG GCT GTA TAT TAT TGT | 288 |
| Leu Gln Met Lys Ser Leu Arg Ala Glu Asp Thr Ala Val Tyr Tyr Cys | |
| 85 90 95 | |
| GCG AGA GAG AAG GCG CTT CGG GGA ATC AGT AGA TAC AAC TAT TAC CTG | 336 |
| Ala Arg Glu Lys Ala Leu Arg Gly Ile Ser Arg Tyr Asn Tyr Tyr Leu | |
| 100 105 110 | |
| GAC GTC TGG GGC AAG GGG ACC ACG GTC ACC GTC TCC TCA | 375 |
| Asp Val Typ Gly Lys Gly Thr Thr Val Thr Val Ser Ser | |
| 115 120 125 | |

(2) INFORMATION FOR SEQ ID NO: 28 (LD2-17-VL)

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 315 Base Pairs
- (B) TYPE: Nucleic acid with corresponding protein
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

- (A) ORGANISM: Homo sapiens
- (B) STRAIN: Not applicable
- (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
- (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
- (E) HAPLOTYPE: Diploid
- (F) TISSUE TYPE: Not applicable
- (G) CELL TYPE: Peripheral lymphocyte B
- (H) CELL LINE: Not applicable
- (I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

- (A) LIBRARY: cDNA LIBRARY, LD2
- (B) CLONE: LD2-17

(viii) POSITION IN GENOME:

- (A) CHROMOSOME / SEGMENT: CHROMOSOME 2
- (B) MAP POSITION: p12
- (C) UNITS: Chromosome band number

(xi) FEATURE:

- (A) NAME / KEY: CDR1, CDR2, CDR3
- (B) LOCATION: join base pair (64..96, 142..162, 259..285)
- (C) IDENTIFICATION METHOD: By similarity with a known sequence
- (D) OTHER INFORMATION:

(2) INFORMATION FOR SEQ ID NO: 29 (LD2-18-VH)
(i) SEQUENCE CHARACTERISTICS:
(A) LENGTH: 375 Base Pairs
(B) TYPE: Nucleic acid with corresponding protein
(C) STRANDEDNESS: Single
(D) TOPOLOGY: Linear
(ii) MOLECULE TYPE: cDNA to mRNA
(iii) HYPOTHETICAL: No
(iv) ANTI-SENSE: No
(vi) ORIGINAL SOURCE: N-terminal fragment
(A) ORGANISM: Homo sapiens
(B) STRAIN: Not applicable
(C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
(D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
(E) HAPLOTYPE: Diploid
(F) TISSUE TYPE: Not applicable
(G) CELL TYPE: Peripheral lymphocyte B
(H) CELL LINE: Not applicable
(I) ORGANELLE: Not applicable
(vii) IMMEDIATE SOURCE: Escherichia coli
(A) LIBRARY: cDNA LIBRARY, LD2
(B) CLONE: LD2-18
(viii) POSITION IN GENOME:
(A) CHROMOSOME / SEGMENT: CHROMOSOME 14
(B) MAP POSITION: q32.3
(C) UNITS: Chromosome band number
(xi) FEATURE:
(A) NAME / KEY: CDR1, CDR2, CDR3
(B) LOCATION: join base pair (91..105, 148..198, 295..342)
(C) IDENTIFICATION METHOD: By similarity with a known sequence
(D) OTHER INFORMATION:

(xi) SEQUENCE DESCRIPTION: SEQ ID NO. 29 (LD2-18-VH)

| | |
|---|-----|
| CAG GTG AAA CTG CTC GAG TCT GGG GGA GGC TTG GTC CAG CCG GGG GGG | 48 |
| Gln Val Lys Leu Leu Glu Ser Gly Gly Gly Leu Val Gln Pro Gly Gly | |
| 5 10 15 | |
| TCC CTG AGA CTC TCC TGT GTA GCG TCT GGA TTC ACC TTC AGG AGT TAT | 96 |
| Ser Leu Arg Leu Ser Cys Val Ala Ser Gly Phe Thr Phe Arg Ser Tyr | |
| 20 25 30 | |
| GGC ATG CAC TGG GTC CGC CAG GCT CCA GGC AAG GGC CTG GAG TGG GTG | 144 |
| Gly Met His Typ Val Arg Gln Ala Pro Gly Lys Gly Leu Glu Typ Val | |
| 35 40 45 | |
| GCT TTT ATA TGG TTT GAT GGA AGT AAT AAA GGA TAT GTA GAC TCC GTG | 192 |
| Ala Phe Ile Typ Phe Asp Gly Ser Asn Lys Gly Tyr Val Asp Ser Val | |
| 50 55 60 | |
| AAG GGC CGA TTC ACC ATC TCC CGA GAC AAT TCC AAG AAC ATG CTC TAT | 240 |
| Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ser Lys Asn Met Leu Tyr | |
| 65 70 75 80 | |
| CTG CAA ATG AAT AGC CTG AGA GCC GAG GAC ACG GCT GTA TAT TAT TGT | 288 |
| Leu Gln Met Asn Ser Leu Arg Ala Glu Asp Thr Ala Val Tyr Tyr Cys | |
| 85 90 95 | |
| GCG AGA GAG AAG GCG CTT CGG GGA ATC AGT AGA TAC AAC TAT TAC CTG | 336 |
| Ala Arg Glu Lys Ala Leu Arg Gly Ile Ser Arg Tyr Asn Tyr Tyr Leu | |
| 100 105 110 | |
| GAC GTC TGG GGC AAG GGG ACC ACG GTA ACC GTC TCC TCA | 375 |
| Asp Val Typ Gly Lys Gly Thr Thr Val Thr Val Ser Ser | |
| 115 120 125 | |

(2) INFORMATION FOR SEQ ID NO: 30 (LD2-18-VL)

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 318 Base Pairs
- (B) TYPE: Nucleic acid with corresponding protein
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

- (A) ORGANISM: Homo sapiens
- (B) STRAIN: Not applicable
- (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
- (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
- (E) HAPLOTYPE: Diploid
- (F) TISSUE TYPE: Not applicable
- (G) CELL TYPE: Peripheral lymphocyte B
- (H) CELL LINE: Not applicable
- (I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

- (A) LIBRARY: cDNA LIBRARY, LD2
- (B) CLONE: LD2-18

(viii) POSITION IN GENOME:

- (A) CHROMOSOME / SEGMENT: CHROMOSOME 2
- (B) MAP POSITION: p12
- (C) UNITS: Chromosome band number

(xi) FEATURE:

- (A) NAME / KEY: CDR1, CDR2, CDR3
- (B) LOCATION: join base pair (64..96, 142..162, 259..288)
- (C) IDENTIFICATION METHOD: By similarity with a known sequence
- (D) OTHER INFORMATION:

(2) INFORMATION FOR SEQ ID NO: 31 (LD2-20-VH)

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 375 Base Pairs
- (B) TYPE: Nucleic acid with corresponding protein
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

- (A) ORGANISM: Homo sapiens
- (B) STRAIN: Not applicable
- (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
- (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
- (E) HAPLOTYPE: Diploid
- (F) TISSUE TYPE: Not applicable
- (G) CELL TYPE: Peripheral lymphocyte B
- (H) CELL LINE: Not applicable
- (I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

- (A) LIBRARY: cDNA LIBRARY, LD2
- (B) CLONE: LD2-20

(viii) POSITION IN GENOME:

- (A) CHROMOSOME / SEGMENT: CHROMOSOME 14
- (B) MAP POSITION: q32.3
- (C) UNITS: Chromosome band number

(xi) FEATURE:

- (A) NAME / KEY: CDR1, CDR2, CDR3
- (B) LOCATION: join base pair (91..105, 148..198, 295..342)
- (C) IDENTIFICATION METHOD: By similarity with a known sequence
- (D) OTHER INFORMATION:

(xi) SEQUENCE DESCRIPTION: SEQ ID NO. 31 (LD2-20-VH)

| | |
|---|-----|
| CAG GTG AAA CTG CTC GAG TCT GGG GGA GGC GTG GTC CAG CCG GGG GGG | 48 |
| Gln Val Lys Leu Leu Glu Ser Gly Gly Gly Val Val Gln Pro Gly Gly | |
| 5 10 15 | |
| TCC CTG AGA CTC TCC TGT GTA GCG TCT GGA TTC ACC TCC AGG AGT TAT | 96 |
| Ser Leu Arg Leu Ser Cys Val Ala Ser Gly Phe Thr Ser Arg Ser Tyr | |
| 20 25 30 | |
| GGC ATG CAC TGG GTC CGC CAG GCT CCA GGA AAG GGC CTG GAG TGG GTG | 144 |
| Gly Met His Typ Val Arg Gln Ala Pro Gly Lys Gly Leu Glu Typ Val | |
| 35 40 45 | |
| GCT TTT ATT TGG TTT GAT GGA AGT AAT AAA GGA TAT GTA GAC TCC GTG | 192 |
| Ala Phe Ile Typ Phe Asp Gly Ser Asn Lys Gly Tyr Val Asp Ser Val | |
| 50 55 60 | |
| AAG GGC CGA TTC ACC ATC TCC CGA GAC AAT TCC AAG AAC ACG CTC TAT | 240 |
| Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ser Lys Asn Thr Leu Tyr | |
| 65 70 75 80 | |
| CTG CAA ATG AAG AGC CTG AGA GCC GAG GAC ACG GCT GTA TAT TAT TGT | 288 |
| Leu Gln Met Lys Ser Leu Arg Ala Glu Asp Thr Ala Val Tyr Tyr Cys | |
| 85 90 95 | |
| GCG AGA GAG AAG GCG CTT CGG GGA ATC AGT AGA TAC AAC TAT TAC CTG | 336 |
| Ala Arg Glu Lys Ala Leu Arg Gly Ile Ser Arg Tyr Asn Tyr Tyr Leu | |
| 100 105 110 | |
| GAC GTC TGG GGC AAG GGG ACC ACG GTC ACC GTC TCC TCA | 375 |
| Asp Val Typ Gly Lys Gly Thr Thr Val Thr Val Ser Ser | |
| 115 120 125 | |

(2) INFORMATION FOR SEQ ID NO: 32 (LD2-20-VL)

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 315 Base Pairs
- (B) TYPE: Nucleic acid with corresponding protein
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

- (A) ORGANISM: Homo sapiens
- (B) STRAIN: Not applicable
- (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
- (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
- (E) HAPLOTYPE: Diploid
- (F) TISSUE TYPE: Not applicable
- (G) CELL TYPE: Peripheral lymphocyte B
- (H) CELL LINE: Not applicable
- (I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

- (A) LIBRARY: cDNA LIBRARY, LD2
- (B) CLONE: LD2-20

(viii) POSITION IN GENOME:

- (A) CHROMOSOME / SEGMENT: CHROMOSOME 2
- (B) MAP POSITION: p12
- (C) UNITS: Chromosome band number

(xi) FEATURE:

- (A) NAME / KEY: CDR1, CDR2, CDR3
- (B) LOCATION: join base pair (64..96, 142..162, 259..285)
- (C) IDENTIFICATION METHOD: By similarity with a known sequence
- (D) OTHER INFORMATION:

(2) INFORMATION FOR SEQ ID NO: 33 (LD1-6-17-VH)

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 384 Base Pairs
- (B) TYPE: Nucleic acid with corresponding protein
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

- (A) ORGANISM: Homo sapiens
- (B) STRAIN: Not applicable
- (C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor
- (D) DEVELOPMENT STAGE: Adult, Rearranged development pattern
- (E) HAPLOTYPE: Diploid
- (F) TISSUE TYPE: Not applicable
- (G) CELL TYPE: Peripheral lymphocyte B
- (H) CELL LINE: Not applicable
- (I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

- (A) LIBRARY: cDNA LIBRARY, LD1
- (B) CLONE: LD1-6-17

(viii) POSITION IN GENOME:

- (A) CHROMOSOME / SEGMENT: Chromosome 14
- (B) MAP POSITION: q32.3
- (C) UNITS: Chromosome band number

(xi) FEATURE:

- (A) NAME / KEY: CDR1, CDR2, CDR3
- (B) LOCATION: join base pair (91..105, 148..198, 295..351)
- (C) IDENTIFICATION METHOD: By similarity with a known sequence
- (D) OTHER INFORMATION:

(xi) SEQUENCE DESCRIPTION: SEQ ID NO. 33 (LD1-6-17-VH)

| | |
|---|-----|
| CAG GTG AAA CTG CTC GAG TCT GGG GGA GGC GTG GTC CAG CCT GGG AGG | 48 |
| Gln Val Lys Leu Leu Glu Ser Gly Gly Gly Val Val Gln Pro Gly Arg | |
| 5 10 15 | |
| TCC CTG AGA CTT TCC TGT GCA GCG TCT GGA TTT ACC TTC AGT AGC TAT | 96 |
| Ser Leu Arg Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Ser Tyr | |
| 20 25 30 | |
| GGA ATG CAC TGG GTC CGC CAG GCT CCA GGA AAG GGG CTG GAG TGG GTG | 144 |
| Gly Met His Typ Val Arg Gln Ala Pro Gly Lys Gly Leu Glu Typ Val | |
| 35 40 45 | |
| ACA GAT ATA TGG TTT GAT GGA GGT AAT AAA CAT TAT GCA GAC TTC GTG | 192 |
| Thr Asp Ile Typ Phe Asp Gly Gly Asn Lys His Tyr Ala Asp Phe Val | |
| 50 55 60 | |
| AAG GGC CGA TTC ACC ATC TCC AGA GAC AAT TCC AAG AAC ACG GGG TTT | 240 |
| Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ser Lys Asn Thr Gly Phe | |
| 65 70 75 80 | |
| CTA CAA ATG AAC AGC CTG AGA GTC GAG GAC ACG GCT GTG TAT TAC TGT | 288 |
| Leu Gln Met Asn Ser Leu Arg Val Glu Asp Thr Ala Val Tyr Tyr Cys | |
| 85 90 95 | |
| GCG AGG GAT TAC TAT AGC GTT ACT AAG AAA CTC AGA CTC CAC TAC TAC | 336 |
| Ala Arg Asp Tyr Tyr Ser Val Thr Lys Lys Leu Arg Leu His Tyr Tyr | |
| 100 105 110 | |
| TAC TAC ATG GAC GTC TGG GGC AAA GGG ACC ACG GTC ACC GTC TCC TCA | 384 |
| Tyr Tyr Met Asp Val Typ Gly Lys Gly Thr Thr Val Thr Val Ser Ser | |
| 115 120 125 | |

(2) INFORMATION FOR SEQ ID NO: 34 (LD1-6-17-VL)

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 315 Base Pairs

(B) TYPE: Nucleic acid with corresponding protein

(C) STRANDEDNESS: Single

(D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: cDNA to mRNA

(iii) HYPOTHETICAL: No

(iv) ANTI-SENSE: No

(vi) ORIGINAL SOURCE: N-terminal fragment

(A) ORGANISM: Homo sapiens

(B) STRAIN: Not applicable

(C) INDIVIDUAL / ISOLATE: Hyperimmune Rhesus D donor

(D) DEVELOPMENT STAGE: Adult, Rearranged development pattern

(E) HAPLOTYPE: Diploid

(F) TISSUE TYPE: Not applicable

(G) CELL TYPE: Peripheral lymphocyte B

(H) CELL LINE: Not applicable

(I) ORGANELLE: Not applicable

(vii) IMMEDIATE SOURCE: Escherichia coli

(A) LIBRARY: cDNA LIBRARY, LD1

(B) CLONE: LD1-6-17

(viii) POSITION IN GENOME:

(A) CHROMOSOME / SEGMENT: Chromosome 2

(B) MAP POSITION: p12

(C) UNITS: Chromosome band number

(xi) FEATURE:

(A) NAME / KEY: CDR1, CDR2, CDR3

(B) LOCATION: join base pair (64..96, 142..162, 259..285)

(C) IDENTIFICATION METHOD: By similarity with a known sequence

(D) OTHER INFORMATION:

Fig. 1a

LD1-28-VH sequence

| | | | | | | | | | | | | | | | | | | |
|----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| 5' | CAG | GTG | AAA | CTG | CTC | GAG | TCT | GGG | GGA | GGC | GTG | GTC | CAG | CCG | GGG | GGG | TCC | CTG |
| | Q | V | K | L | L | E | S | G | G | G | V | V | Q | P | G | G | S | L |
| | AGA | CTC | TCC | TGT | GAA | GCG | TCT | GGA | TTC | GCC | CTC | AGA | AGT | TCT | GGC | ATG | CAC | TGG |
| | R | L | S | C | E | A | S | G | F | A | L | R | S | S | G | M | H | W |
| | GTC | CGC | CAG | GCT | CCT | GGC | AAG | GGG | ATG | GAG | TGG | GTG | GCA | CTT | ATA | TGG | TTT | GAT |
| | V | R | Q | A | P | G | K | G | M | E | W | V | A | L | I | W | F | D |
| | GGA | AGT | ATC | AGA | TCG | TAT | GCA | GAA | TCC | GCG | AAG | GGC | CGA | TTC | ACC | ATC | TCC | AGA |
| | G | S | I | R | S | Y | A | E | S | A | K | G | R | F | T | I | S | R |
| | GAC | ACT | TCC | AAG | AAC | ACC | CTA | TAT | CTC | CAT | ATG | CGC | AGT | CTG | AGT | GCC | GAC | GAC |
| | D | T | S | K | N | T | L | Y | L | H | M | R | S | L | S | A | D | D |
| | ACG | GAT | GTG | TTT | TAC | TGT | GCG | AGA | GAC | AAG | GCG | GTT | CGG | GGA | ATT | AAC | AGG | TAC |
| | T | D | V | F | Y | C | A | R | D | K | A | V | R | G | I | N | R | Y |
| | AAC | TAT | TAC | ATG | GAC | GTC | TGG | GTC | AAA | GGG | ACC | ACG | GTC | ACC | GTC | TCC | TCA | 3' |
| | N | Y | Y | M | D | V | W | V | K | G | T | T | V | T | V | S | S | |

Fig. 1b

LD1-28-VL sequence

```

5'  GTG  ATG  ACC  CAG  TCT  CCA  TCC  TCC  CTG  TCT  GCA  TCT  GTA  GGA  GAC  AGA  GTC  ACC
      9      18      27      36      45      54
      V  M  T  Q  S  P  S  S  L  S  A  S  V  G  D  R  V  T
-----
      63      72      81      90      99      108
ATC  ACT  TGC  CGG  GCA  AGT  CAG  AAC  ATT  ATC  CGC  TAT  TTA  AAT  TGG  TAT  CAG  CAG
      I  T  C  R  A  S  Q  N  I  I  R  Y  L  N  W  Y  Q  Q
      <----- CDR1 ----->
      117      126      135      144      153      162
AAG  CCA  GGG  AAA  GCC  CCT  AGG  CTC  CTG  ATC  TAT  GGT  GCG  TCC  ACT  TTG  CAA  AGT
      K  P  G  K  A  P  R  L  L  I  Y  G  A  S  T  L  Q  S
      <----- CDR2 ----->
      171      180      189      198      207      216
GGG  GTC  CCA  TCA  AGG  TTC  AGT  GGC  AGT  GGA  TCT  GGG  ACA  GAT  TTC  ACT  CTC  ACC
      G  V  P  S  R  F  S  G  S  G  S  G  T  D  F  T  L  T
      225      234      243      252      261      270
ATC  AGT  AGT  CTG  CAA  CCT  GAA  GAT  TTT  GCA  ACT  TAC  TAC  TGT  CAA  CAG  AGT  TAC
      I  S  S  L  Q  P  E  D  F  A  T  Y  Y  C  Q  Q  S  Y
      <----- CDR3 ----->
      279      288      297      306      315
CGT  ACC  CCT  CCA  TTC  ACT  TTC  GGC  CCT  GGG  ACC  AAA  GTG  GAG  ATC  AAA  3'
      R  T  P  P  F  T  F  G  P  G  T  K  V  E  I  K

```


Fig. 2b

LD1-40-VL sequence

[illegible]

Fig. 3a

LD1-52-VH sequence

```

      9      18      27      36      45      54
5'  CAG GTG AAA CTG CTC GAG TCT GGG GGA GGC GTG GTC CAG CCG GGG GGG TCC CTG
    Q  V  K  L  L  E  S  G  G  G  V  V  Q  P  G  G  S  L
    -----
      63      72      81      90      99      108
AGA CTC TCC TGT GAA GCG TCT GGA TTC GCC CTC AGA AGT TCT GGA ATG CAC TGG
    R  L  S  C  E  A  S  G  F  A  L  R  S  S  G  M  H  W
    -----
                                CDR1 -----
      117      126      135      144      153      162
GTC CGC CAG GCT CCT GGC AAG GGG CTG GAG TGG GTG GCA CTT ATA TGG TTT GAT
    V  R  Q  A  P  G  K  G  L  E  W  V  A  L  I  W  F  D
    -----
                                CDR2 -----
      171      180      189      198      207      216
GGA AGT ATC AGA TCG TAT GCA GAA TCC GTG AAG GGC CGA TTC ACC ATC TCC AGA
    G  S  I  R  S  Y  A  E  S  V  K  G  R  F  T  I  S  R
    -----
    CDR2 -----
      225      234      243      252      261      270
GAC ACT TCC AAG AAC ACC CTA TAT CTC CAA ATG CGC AGT CTG AGT GCC GAC GAC
    D  T  S  K  N  T  L  Y  L  Q  M  R  S  L  S  A  D  D
    -----
      279      288      297      306      315      324
ACG GCT GTG TAT TAC TGT GCG AGA GAC AAG GCG GTT CGG GGA ATT AGC AGG TAC
    T  A  V  Y  Y  C  A  R  D  K  A  V  R  G  I  S  R  Y
    -----
                                CDR3 -----
      333      342      351      360      369
AAC TAT TAC ATG GAC GTC TGG GGC AAA GGG ACC ACG GTC ACC GTC TCC TCA 3'
    N  Y  Y  M  D  V  W  G  K  G  T  T  V  T  V  S  S
    -----
    CDR3 -----

```

Fig. 3b

LD1-52-VL sequence

```

5'  GTG  ATG  ACC  CAG  TCT  CCA  TCC  TCC  CTG  TCT  GCA  TCT  GTA  GGA  GAC  AGA  GTC  ACC
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    V    M    T    Q    S    P    S    S    L    S    A    S    V    G    D    R    V    T

      9          18          27          36          45          54
    ATC  ACT  TGC  CGG  GCA  AGT  CAG  AAC  ATT  ATC  CGC  TAT  TTA  AAT  TGG  TAT  CAG  CAG
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    I    T    C    R    A    S    Q    N    I    I    R    Y    L    N    W    Y    Q    Q

      63          72          81          90          99          108
    AAG  CCA  GGG  AAA  GCC  CCT  AGG  CTC  CTG  ATC  TAT  GGT  GCG  TCC  ACT  TTG  CAA  AGT
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    K    P    G    K    A    P    R    L    L    I    Y    G    A    S    T    L    Q    S

      117          126          135          144          153          162
    GGG  GTC  CCA  TCA  AGG  TTC  AGT  GGC  AGT  GGA  TCT  GGG  ACA  GAT  TTC  ACT  CTC  ACC
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    G    V    P    S    R    F    S    G    S    G    S    G    T    D    F    T    L    T

      171          180          189          198          207          216
    ATC  AGT  AGT  CTG  CAA  CCT  GAA  GAT  TTT  GCA  ACT  TAC  TAC  TGT  CAA  CAG  AGT  TAC
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    I    S    S    L    Q    P    E    D    F    A    T    Y    Y    C    Q    Q    S    Y

      225          234          243          252          261          270
    CGT  ACC  CCT  CCA  TTC  ACT  TTC  GGC  CCT  GGG  ACC  AAA  GTG  GAG  ATC  AAA  3'
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    R    T    P    P    F    T    F    G    P    G    T    K    V    E    I    K

      279          286          297          306          315
    CDR1 ←-----→ CDR2 ←-----→
    CDR3 ←-----→
  
```

Fig. 4a

LD1-84-VH sequence

| | | | | | | | | | | | | | | | | | | |
|----|----------|-----|-----|-----|-----|-----|-----|-----|----------|-----|-----|-----|----------|-----|-----|-----|-----|-----|
| | | | 9 | | | 18 | | | 27 | | | 36 | | | 45 | | | 54 |
| 5' | CAG | GTA | AAA | CTG | CTC | GAG | TCT | GGG | GGG | GGC | GTG | GTC | CAG | CCG | GGG | GGG | TCC | CTG |
| | Q | V | K | L | L | E | S | G | G | G | V | V | Q | P | G | G | S | L |
| | | | 63 | | | 72 | | | 81 | | | 90 | | | 99 | | | 108 |
| | AGA | CTC | TCC | TGT | GAA | GCG | TCT | GGA | TTC | ACC | CTC | AGA | AGT | TCT | GGC | ATG | CAC | TGG |
| | R | L | S | C | E | A | S | G | F | T | L | R | S | S | G | M | H | W |
| | | | | | | | | | | | | | ← CDR1 → | | | | | |
| | GTC | CGC | 117 | | | 126 | | | 135 | | | 144 | | | 153 | | | 162 |
| | GTC | CGC | CAG | GCT | CCT | GGC | AAG | GGG | CTG | GAG | TGG | GTG | ACA | CTT | ATA | TGG | TTT | GAT |
| | V | R | Q | A | P | G | K | G | L | E | W | V | T | L | I | W | F | D |
| | | | | | | | | | | | | | ← CDR2 → | | | | | |
| | GGA | AGT | 171 | | | 180 | | | 189 | | | 198 | | | 207 | | | 216 |
| | GGA | AGT | ATC | AGA | TCG | TAT | GCA | GAA | TCC | GTG | AAG | GGC | CGA | TTC | ACC | ATC | TCC | AGA |
| | G | S | I | R | S | Y | A | E | S | V | K | G | R | F | T | I | S | R |
| | ← CDR2 → | | | | | | | | | | | | | | | | | |
| | GAC | ACT | 225 | | | 234 | | | 243 | | | 252 | | | 261 | | | 270 |
| | GAC | ACT | TCC | AAG | AAC | ACC | CTA | TAT | CTC | CAA | ATG | CGC | AGT | CTG | AGT | GCC | GAC | GAC |
| | D | T | S | K | N | T | L | Y | L | Q | M | R | S | L | S | A | D | D |
| | | | 279 | | | 288 | | | 297 | | | 306 | | | 315 | | | 324 |
| | ACG | GCT | GTG | TAT | TAC | TGT | GCG | AGA | GAC | AAG | GCG | GTT | CGG | GGA | ATT | AGC | AGG | TAC |
| | T | A | V | Y | Y | C | A | R | D | K | A | V | R | G | I | S | R | Y |
| | | | | | | | | | ← CDR3 → | | | | | | | | | |
| | AAC | TAT | 333 | | | 342 | | | 351 | | | 360 | | | 369 | | | 378 |
| | AAC | TAT | TAC | ATG | GAC | GTC | TGG | GGC | AAA | GGG | ACC | ACG | GTC | ACC | GTC | TCC | TCA | 3' |
| | N | Y | Y | M | D | V | W | G | K | G | T | T | V | T | V | S | S | |
| | ← CDR3 → | | | | | | | | | | | | | | | | | |

Fig. 5a

LD1-98-VH sequence

```

5'  CAG GTG AAA CTG CTC GAG TCT GGG GGA GGC GTG GTC CAG CCG GGG GGG TCC CTG
    Q  V  K  L  L  E  S  G  G  G  V  V  Q  P  G  G  S  L
-----
    63      72      81      90      99      108
AGA CTC TCC TGT GAA GCG TCT GGA TTC GCC CTC AGA AGT TCT GAC ATA CAC TGG
    R  L  S  C  E  A  S  G  F  A  L  R  S  S  D  I  H  W
-----
                                CDR1
    117      126      135      144      153      162
GTC CGC CAG GCT CCT GGC AAG GGG CTG GAG TGG GTG GCA CTT ATA TGG TTT GAT
    V  R  Q  A  P  G  K  G  L  E  W  V  A  L  I  W  F  D
-----
                                CDR2
    171      180      189      198      207      216
GGA AGT ATC AGA TCG TAT GCA GAA TCC GTG AAG GGC CGA TTC ACC ATC TCC AGA
    G  S  I  R  S  Y  A  E  S  V  K  G  R  F  T  I  S  R
-----
                                CDR2
    225      234      243      252      261      270
GAC ACT TCC AAG AAC ACC CTA TAT CTC CAA ATG CGC AGT CTG AGT GCC GAC GAC
    D  T  S  K  N  T  L  Y  L  Q  M  R  S  L  S  A  D  D
-----
    279      288      297      306      315      324
ACG CGT GTG TAT TAC TGT GCG AGA GAC AAG GCG GTT CGG GGA ATT AGC AGG TAC
    T  A  V  Y  Y  C  A  R  D  K  A  V  R  G  I  S  R  Y
-----
                                CDR3
    333      342      351      360      369
AAC TAT TAC ATG GAC GTC TGG GGC AAA GGG ACC ACG GTC ACC GTC TCC TCA 3'
    N  Y  Y  M  D  V  W  G  K  G  T  T  V  T  V  S  S
-----
    CDR3

```


Fig. 5b

LD1-98-VL sequence

```

5'  GTG  ATG  ACC  CAG  TCT  CCA  TCC  TCC  CTG  TCT  GCA  TCT  GTA  GGA  GAC  AGA  GTC  ACC
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    V   M   T   Q   S   P   S   S   L   S   A   S   V   G   D   R   V   T

      63      72      81      90      99      108
ATC  ACT  TGC  CGG  GCA  AGT  CAG  AAC  ATT  ATC  CGC  TAT  TTA  AAT  TGG  TAT  CAG  CAG
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    I   T   C   R   A   S   Q   N   I   I   R   Y   L   N   W   Y   Q   Q

      117      126      135      144      153      162
AAG  CCA  GGG  AAA  GCC  CCT  AGG  CTC  CTG  ATC  TAT  GGT  GCG  TCC  ACT  TTG  CAA  AGT
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    K   P   G   K   A   P   R   L   L   I   Y   G   A   S   T   L   Q   S

      171      180      189      198      207      216
GGG  GTC  CCA  TCA  AGG  TTC  AGT  GGC  AGT  GGA  TCT  GGG  ACA  GAT  TTC  ACT  CTC  ACC
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    G   V   P   S   R   F   S   G   S   G   S   G   T   D   F   T   L   T

      225      234      243      252      261      270
ATC  AGT  AGT  CTG  CAA  CCT  GAA  GAT  TTT  GCA  ACT  TAC  TAC  TGT  CAA  CAG  AGT  TAC
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    I   S   S   L   Q   P   E   D   F   A   T   Y   Y   C   Q   Q   S   Y

      279      288      297      306      315
CGT  ACC  CCT  CCA  TTC  ACT  TTC  GGC  CCT  GGG  ACC  AAA  GTG  GAG  ATC  AAA  3'
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    R   T   P   P   F   T   F   G   P   G   T   K   V   E   I   K

      CDR1  CDR2  CDR3

```

Fig. 6a

LD1-110-VH sequence

| | | | | | | | | | | | | | | | | | | | | | | | | |
|----|----------|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|----------|-----|-----|-----|-----|-----|--|--|-----|--|--|--|
| | 9 | | | | 18 | | | | 27 | | | | 36 | | | | 45 | | | | 54 | | | |
| 5' | CAG | GTG | AAA | CTG | CTC | GAG | TCT | GGG | GGA | GGC | GTG | GTC | CAG | CCT | GGG | AGG | TCC | CTG | | | | | | |
| | Q | V | K | L | L | E | S | G | G | G | V | V | Q | P | G | R | S | L | | | | | | |
| | 63 | | | | 72 | | | | 81 | | | | 90 | | | | 99 | | | | 108 | | | |
| | AGA | CTC | TCC | TGT | ATA | GCG | TCT | GGA | TTC | ACC | CTC | AGG | AAT | TAT | GCC | ATG | CAC | TGG | | | | | | |
| | R | L | S | C | I | A | S | G | F | T | L | R | N | Y | A | M | H | W | | | | | | |
| | | | | | | | | | | | | | ← CDR1 → | | | | | | | | | | | |
| | 117 | | | | 126 | | | | 135 | | | | 144 | | | | 153 | | | | 162 | | | |
| | GTC | CGC | CAG | GCT | CCA | GGA | AAG | GGG | CTG | GAG | TGG | GTG | GCA | GGT | ATA | TGG | TTT | GAT | | | | | | |
| | V | R | Q | A | P | G | K | G | L | E | W | V | A | G | I | W | F | D | | | | | | |
| | | | | | | | | | | | | | ← CDR2 → | | | | | | | | | | | |
| | 171 | | | | 180 | | | | 189 | | | | 198 | | | | 207 | | | | 216 | | | |
| | GGA | AGC | AAC | AAA | AAC | TAT | GCA | GAC | TCC | GTG | AAG | GGC | CGA | TTC | ACC | ATC | TCC | AGA | | | | | | |
| | G | S | N | K | N | Y | A | D | S | V | K | G | R | F | T | I | S | R | | | | | | |
| | CDR2 → | | | | | | | | | | | | | | | | | | | | | | | |
| | 225 | | | | 234 | | | | 243 | | | | 252 | | | | 261 | | | | 270 | | | |
| | GAC | AAC | TCC | AAG | AAC | ACT | CTG | TTT | CTG | CAC | ATG | AAC | AGC | CTG | AGA | GCC | GAG | GAC | | | | | | |
| | D | N | S | K | N | T | L | F | L | H | M | N | S | L | R | A | E | D | | | | | | |
| | 279 | | | | 288 | | | | 297 | | | | 306 | | | | 315 | | | | 324 | | | |
| | ACG | GCT | ACA | TAT | TAC | TGT | GCG | AGA | GAG | AGG | GCG | ATT | CGG | GGA | ATC | AGT | AGA | TAC | | | | | | |
| | T | A | T | Y | Y | C | A | R | E | R | A | I | R | G | I | S | R | Y | | | | | | |
| | ← CDR3 → | | | | | | | | | | | | | | | | | | | | | | | |
| | 333 | | | | 342 | | | | 351 | | | | 360 | | | | 369 | | | | 378 | | | |
| | AAT | TAC | TAC | ATG | GAC | GTC | TGG | GGC | AAG | GGG | ACC | ACG | GTC | ACC | GTC | TCC | TCA | 3' | | | | | | |
| | N | Y | Y | M | D | V | W | G | K | G | T | T | V | T | V | S | S | | | | | | | |
| | CDR3 → | | | | | | | | | | | | | | | | | | | | | | | |

Fig. 6b

LD1-110-VL sequence

```

5'  GTG  ATG  ACC  CAG  TCT  CCA  TCC  TCC  CTG  TCT  GCA  TCT  GTA  GGA  GAC  AGA  GTC  ACC
    V   M   T   Q   S   P   S   S   L   S   A   S   V   G   D   R   V   T

      9      18      27      36      45      54
-----
      63      72      81      90      99      108
ATC  ACT  TGC  CGG  GCA  AGT  CAG  AGC  ATT  CGA  AGC  TCT  TTA  AAT  TGG  TAT  CAG  CAG
    I   T   C   R   A   S   Q   S   I   R   S   S   L   N   W   Y   Q   Q
      117      126      135      144      153      162
AAA  CCA  GGG  AAA  GCC  CCT  AAA  GTC  CTG  ATC  TAT  GCT  GCA  TCC  AGT  TTG  CAA  AGT
    K   P   G   K   A   P   K   V   L   I   Y   A   A   S   S   L   Q   S
      171      180      189      198      207      216
GGG  GTC  CCA  TCC  AGG  TTC  AGT  GGC  AGA  GGA  TCT  GGG  ACA  GAT  TTC  ACT  CTC  ACC
    G   V   P   S   R   F   S   G   R   G   S   G   T   D   F   T   L   T
      225      234      243      252      261      270
ATC  AGC  AGT  CTG  CAG  CCT  GAA  GAT  TTT  GCG  ACT  TAT  TAT  TGT  CAA  CAG  AGT  TCC
    I   S   S   L   Q   P   E   D   F   A   T   Y   Y   C   Q   Q   S   S
      279      288      297      306      315
AGT  TCC  TCG  TGG  ACG  TTC  GGC  CAA  GGG  ACC  AAG  GTG  GAA  ATC  AAA 3'
    S   S   S   W   T   F   G   Q   G   T   K   V   E   I   K
      CDR3
  
```

CDR1: 135-144
 CDR2: 198-207
 CDR3: 279-288

Fig. 7a

LD1-117-VH sequence

```

5'  CAG GTG AAA CTG CTC GAG TCA GGA GGA GGC GTG GTC CAG CCT GGG AAG TCC CTG
    Q  V  K  L  L  E  S  G  G  G  V  V  Q  P  G  K  S  L

    63      72      81      90      99      108
AGA CTT TCC TGT GCA GCG TCT GGA TTC AGT TTC AAT AGC CAT GGN ATG CAC TGG
    R  L  S  C  A  A  S  G  F  S  F  N  S  H  G  M  H  W

                                CDR1
GTC CGC CAG GCT CCA GGC AAG GGG CTG GAG TGG GTG GCA TTT ATT TGG TTT GAT
    V  R  Q  A  P  G  K  G  L  E  W  V  A  F  I  W  F  D

                                CDR2
GGC AGT AAT AAA TAC TAT GCA GAC TCC GTG AAG GGC CGT TTC ACC ATC ACC AGA
    G  S  N  K  Y  Y  A  D  S  V  K  G  R  F  T  I  T  R

    CDR2
GAC AAC TCC AAG AAC ACG CTG TAT CTN CAA ATG AAC AGC CTG AGA GCC GAG GAC
    D  N  S  K  N  T  L  Y  L  Q  M  N  S  L  R  A  E  D

    279      288      297      306      315      324
ACG GCT GTC TAT TAC TGT GCG AGA GAG ACC TCA GTA AGG CTA GGG TAT AGC CGC
    T  A  V  Y  Y  C  A  R  E  T  S  V  R  L  G  Y  S  R

                                CDR3
TAC AAT TAC TAC ATG GAC GTC TGG GCC AAA GGG ACC ACG GTC ACC ATC TCG TCA 3'
    Y  N  Y  Y  M  D  V  W  A  K  G  T  T  V  T  I  S  S

    CDR3

```


Fig. 8a

LD2-1-VH sequence

```

5'  CAG GTG AAA CTG CTC GAG TCT GGG GGA GGC GTG GTC CAG CCG GGG GGG TCC CTG
    Q  V  K  L  L  E  S  G  G  G  V  V  Q  P  G  G  S  L
-----
      63      72      81      90      99      108
AGA CTC TCC TGT GTA GCG TCT GGA TTC ACC CTC AGG AGT TAT GGC ATG CAC TGG
    R  L  S  C  V  A  S  G  F  T  L  R  S  Y  G  M  H  W
-----
                                CDR1
      117      126      135      144      153      162
GTC CGC CAG GCT CCA GGA AAG GGC CTG GAG TGG GTG GCT TTT ATA TGG TTT GAT
    V  R  Q  A  P  G  K  G  L  E  W  V  A  F  I  W  F  D
-----
                                CDR2
      171      180      189      198      207      216
GGA AGT AAT AAA GGA TAT GTA GAC TCC GTG AAG GGC CGA TTC ACC ATC TCC CGA
    G  S  N  K  G  Y  V  D  S  V  K  G  R  F  T  I  S  R
-----
                                CDR2
      225      234      243      252      261      270
GAC AAT TCC AAG AAC ATG GTC TAT CTC CAA ATG AAC AGC CTG AGA GCC GAT GAC
    D  N  S  K  N  M  V  Y  L  Q  M  N  S  L  R  A  D  D
-----
      279      288      297      306      315      324
ACG GCT GTA TAT TAT TGT GCG AGA GAG AAG GCG CTT CGG GGA ATC AGC AGA TAC
    T  A  V  Y  Y  C  A  R  E  K  A  L  R  G  I  S  R  Y
-----
                                CDR3
      333      342      351      360      369
AAC TAT TAC CTG GAC GTC TGG GGC AAG GGG ACC ACG GTC ACC GTC TCC TCA 3'
    N  Y  Y  L  D  V  W  G  K  G  T  T  V  T  V  S  S
-----
                                CDR3

```

Fig. 8b

LD2-1-VL sequence

```

5'  GTG GTG ACT CAG CCA CCC TCA GCG TCT GGG ACC CCC GGA CAG AGG GTC ACC ATC
    V  V  T  Q  P  P  S  A  S  G  T  P  G  Q  R  V  T  I
-----
    TCT TGT TCT GGA AGC AAC TCC ATC CTT GGA AGT AAG TAT GTA TAC TGG TAC CAG
    S  C  S  G  S  N  S  I  L  G  S  K  Y  V  Y  W  Y  Q
-----
          63          72          81          90          99          108
          <----- CDR1 ----->
    AAA CTC CCA GGA ACG GCC CCC AAA CTC CTC ATC TAT AAG AAT GAT CAG CGG CCC
    K  L  P  G  T  A  P  K  L  L  I  Y  K  N  D  Q  R  P
-----
                                <----- CDR2 ----->
    TCA GGG GTC TCT GAC CGA TTC TCT GGC TCC AAG TCT GGC ACC TCG GCC TCC CTG
    S  G  V  S  D  R  F  S  G  S  K  S  G  T  S  A  S  L
-----
    >
    GCC ATC AGT GGG CTC CGG TCC GAG GAT GAG GCT GAC TAT TAC TGT GCA CCA TGG
    A  I  S  G  L  R  S  E  D  E  A  D  Y  Y  C  A  P  W
-----
    GAT GCC AAC CTG GGT GGC CCG GTG TTC GGC GGA GGG ACC AAG CTG ACC GTC CTA
    D  A  N  L  G  G  P  V  F  G  G  G  T  K  L  T  V  L
-----
          279          288          297          306          315          324
          <----- CDR3 ----->
    AGT CAG CCC 3'
    S  Q  P

```

Fig. 9a

LD2-4-VH sequence

| | | | | | | | | | | | | | | | | | | |
|----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| | | 9 | | 18 | | 27 | | 36 | | 45 | | 54 | | | | | | |
| 5' | CAG | GTG | AAA | CTG | CTC | GAG | TCG | GGG | GGA | GGC | GTG | GTC | CAG | CCG | GGG | GGG | TCC | CTG |
| | Q | V | K | L | L | E | S | G | G | G | V | V | Q | P | G | G | S | L |
| | | 63 | | 72 | | 81 | | 90 | | 99 | | 108 | | | | | | |
| | AGA | CTC | TCC | TGT | GAA | GCG | TCT | GGA | TTC | ACC | CTC | AGA | AGT | TCT | GGC | ATG | CAC | TGG |
| | R | L | S | C | E | A | S | G | F | T | L | R | S | S | G | M | H | W |
| | | 117 | | 126 | | 135 | | 144 | | 153 | | 162 | | | | | | |
| | GTC | CGC | CAG | GCT | CCT | GGC | AAG | GGG | CTG | GAG | TGG | GTG | GCA | CTT | ATA | TGG | TTT | GAT |
| | V | R | Q | A | P | G | K | G | L | E | W | V | A | L | I | W | F | D |
| | | 171 | | 180 | | 189 | | 198 | | 207 | | 216 | | | | | | |
| | GGA | AGT | ATC | AGA | TCG | TAT | GCA | GAA | TCC | GTG | AAG | GGC | CGA | TTC | ACC | ATC | TCC | AGA |
| | G | S | I | R | S | Y | A | E | S | V | K | G | R | F | T | I | S | R |
| | | 225 | | 234 | | 243 | | 252 | | 261 | | 270 | | | | | | |
| | GAC | ACT | TCC | AAG | AAC | ACC | CTA | TAT | CTC | CAA | ATG | CGC | AGT | CTG | AGT | GCC | GAC | GAC |
| | D | T | S | K | N | T | L | Y | L | Q | M | R | S | L | S | A | D | D |
| | | 279 | | 288 | | 297 | | 306 | | 315 | | 324 | | | | | | |
| | ACG | GCT | GTG | TAT | TAC | TGT | GCG | AGA | GAC | AAG | GCG | GTT | CGG | GGA | ATT | AGC | AGG | TAC |
| | T | A | V | Y | Y | C | A | R | D | K | A | V | R | G | I | S | R | Y |
| | | 333 | | 342 | | 351 | | 360 | | 369 | | | | | | | | |
| | AAC | TAT | TAC | ATG | GAC | GTC | TGG | GGC | AAA | GGG | ACC | ACG | GTC | ACC | GTC | TCC | TCA | 3' |
| | N | Y | Y | M | D | V | W | G | K | G | T | T | V | T | V | S | S | |

CDR1

CDR2

CDR3

Fig. 9b

LD2-4-VL sequence

```

5'  GTG  ATG  ACC  CAG  TCT  CCA  TCC  TCC  CTG  TCT  GCA  TCT  GTA  GGA  GAC  AGA  GTC  ACC
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    V   M   T   Q   S   P   S   S   L   S   A   S   V   G   D   R   V   T

      9      18      27      36      45      54
    ATC  ACT  TGC  CGG  ACA  AGT  CAG  ACC  ATT  AGC  AGA  AAT  TTA  CAT  TGG  TAT  CAG  CAG
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    I   T   C   R   T   S   Q   T   I   S   R   N   L   H   W   Y   Q   Q

      63      72      81      90      99      108
    AAA  CCA  GGG  AAA  GCC  CCT  AAG  CTC  CTG  ATC  TAT  GCT  ACA  TCC  AGT  TTG  CAA  AGT
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    K   P   G   K   A   P   K   L   L   I   Y   A   T   S   S   L   Q   S

      117      126      135      144      153      162
    GGG  GTC  CCA  TCA  AGG  TTC  AGT  GGC  AGT  GGA  TCT  GGG  ACA  GAT  TTC  ACT  CTC  ACC
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    G   V   P   S   R   F   S   G   S   G   S   G   T   D   F   T   L   T

      171      180      189      198      207      216
    ATC  AAT  AGT  CTA  CAA  CCT  GAA  GAT  TTT  GCA  ACT  TAC  TAC  TGT  CAA  CAG  AGT  TAC
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    I   N   S   L   Q   P   E   D   F   A   T   Y   Y   C   Q   Q   S   Y

      225      234      243      252      261      270
    ACT  ACC  CCT  TCG  TTC  GGC  CAA  GGG  ACC  AAG  GTG  GAA  ATC  AAA  3'
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    T   T   P   S   F   G   Q   G   T   K   V   E   I   K

      279      288      297      306      315
    CDR1 ←-----→
    CDR2 ←-----→
    CDR3 →-----→

```

Fig. 10a

LD2-5-VH sequence

```

      9      18      27      36      45      54
5'  CAG GTG AAA CTG CTC GAG TCT GGG GGA GGC TTG GTC CAG CCG GGG GGG TCC CTG
    ---
    Q  V  K  L  L  E  S  G  G  G  L  V  Q  P  G  G  S  L

      63      72      81      90      99      108
    AGA CTC TCC TGT GTA GCG TCT GGA TTC ACC TTC AGG AGT TAT GGA ATG CAC TGG
    ---
    R  L  S  C  V  A  S  G  F  T  F  R  S  Y  G  M  H  W

                                ← CDR1 →

      117      126      135      144      153      162
    GTC CGC CAG GCT CCA GGA AAG GGC CTG GAG TGG GTG GCT TTT ATA TGG TTT GAT
    ---
    V  R  Q  A  P  G  K  G  L  E  W  V  A  F  I  W  F  D

                                ← CDR2 →

      171      180      189      198      207      216
    GGA AGT AAT AAA GGA TAT GTA GAC TCC GTG AAG GGC CGA TTC ACC ATC TCC CGA
    ---
    G  S  N  K  G  Y  V  D  S  V  K  G  R  F  T  I  S  R

    ----- CDR2 ----->

      225      234      243      252      261      270
    GAC AAT TCC AAG AAC ATG CTC TAT CTG CAA ATG AAT AGC CTG AGA GCC GAG GAC
    ---
    D  N  S  K  N  M  L  Y  L  Q  M  N  S  L  R  A  E  D

      279      288      297      306      315      324
    ACG GCT GTA TAT TAT TGT GCG AGA GAG AAG GCG CTT CGG GGA ATC AGT AGA TAC
    ---
    T  A  V  Y  Y  C  A  R  E  K  A  L  R  G  I  S  R  Y

                                ← CDR3 →

      333      342      351      360      369
    AAC TAT TAC CTG GAC GTC TGG GGC AAG GGG GCC ACG GTC ACC GTC TCC TCA 3'
    ---
    N  Y  Y  L  D  V  W  G  K  G  A  T  V  T  V  S  S

    ----- CDR3 ----->

```

Fig. 10b

LD2-5-VL sequence

```

      9      18      27      36      45      54
5'  GTG ATG ACC CAG TCT CCA TCC TCC CTG TCT GTA TCT ATA GGC GAC AGA GTC ACC
    ---
      V  M  T  Q  S  P  S  S  L  S  V  S  I  G  D  R  V  T

      63      72      81      90      99      108
ATC ACT TGC CGG GCA AGT CAG AGC GTT ACC AGG TCT TTA AAT TGG TAT CAG CAG
    ---
      I  T  C  R  A  S  Q  S  V  T  R  S  L  N  W  Y  Q  Q

      CDR1
      117      126      135      144      153      162
AAA CCA GGG AAA GCC CCT AGG CTC CTA ATC TTT GGT GCG TCC ACT TTG CAA AGT
    ---
      K  P  G  K  A  P  R  L  L  I  F  G  A  S  T  L  Q  S

      CDR2
      171      180      189      198      207      216
GGG GTC CCA TCA AGG TTC AGT GGC AGT GGA TCT GGG ACA GAT TTC ACC CTC ACC
    ---
      G  V  P  S  R  F  S  G  S  G  S  G  T  D  F  T  L  T

      225      234      243      252      261      270
ATC AGC AGT CTG CAA CCT GAG GAT TTT GGA ACT TAC TAC TGT CAA CAG AAT TAC
    ---
      I  S  S  L  Q  P  E  D  F  G  T  Y  Y  C  Q  Q  N  Y

      279      288      297      306      315
AGG ACC CCT CAG TGG ACG TTC GGC CAA GGG ACC AAG GTA GAA ATC AAA 3'
    ---
      R  T  P  Q  W  T  F  G  Q  G  T  K  V  E  I  K

      CDR3

```

Fig. 11a

LD2-10-VH sequence

[illegible]

Fig. 11b

LD2-10-VL sequence

```

5'  GTG GTG ACT CAG GAG CCC TCA CTG ACT GTG TCC CCA GGA GGG ACA GTC ACT CTC
    ---
    V  V  T  Q  E  P  S  L  T  V  S  P  G  G  T  V  T  L

      9      18      27      36      45      54
    ACC TGT GCT TCC AGC ACT GGG GCA GTC ACC AGG GGT TAC TAT CCA AAC TGG TTC
    ---
    T  C  A  S  S  T  G  A  V  T  R  G  Y  Y  P  N  W  F

      63      72      81      90      99      108
    CAG CAG AAG CCT GGA CAA GCA CCC AGG GCA CTG ATT TAT AGT ACA AAC AAA AAA
    ---
    Q  Q  K  P  G  Q  A  P  R  A  L  I  Y  S  T  N  K  K

      117      126      135      144      153      162
    CAC TCC TGG ACC CCT GCC CGG TTC TCA GGC TCC CTC CTT GGG GGC AAA GCT GCC
    ---
    H  S  W  T  P  A  R  F  S  G  S  L  L  G  G  K  A  A

      171      180      189      198      207      216
    CTG ACA CTG TCA GGT GTG CAG CCT GAA GAC GAG GCT GAA TAT TAC TGC CTG CTC
    ---
    L  T  L  S  G  V  Q  P  E  D  E  A  E  Y  Y  C  L  L

      225      234      243      252      261      270
    TAC TAT GGT GGT GCT CAA CTC GTA TTC GGC GGA GGG ACC AAG CTG ACC GTC CTA
    ---
    Y  Y  G  G  A  Q  L  V  F  G  G  G  T  K  L  T  V  L

      279      288      297      306      315      324
    CGT CAG CCC 3'
    ---
    R  Q  P
  
```

CDR1: 117-153
 CDR2: 144-162
 CDR3: 279-333

Fig. 12a

LD2-11-VH sequence

| | | | | | | | | | | | | | | | | | | | | | | | | |
|----|------|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|------|-----|-----|-----|-----|-----|--|--|-----|--|--|--|
| 5' | 9 | | | | 18 | | | | 27 | | | | 36 | | | | 45 | | | | 54 | | | |
| | CAG | GTG | AAA | CTG | CTC | GAG | TCG | GGG | GGA | GGC | GTG | GTC | CAG | CCG | GGG | GGG | TCC | CTG | | | | | | |
| | Q | V | K | L | L | E | S | G | G | G | V | V | Q | P | G | G | S | L | | | | | | |
| | 63 | | | | 72 | | | | 81 | | | | 90 | | | | 99 | | | | 108 | | | |
| | AGA | CTC | TCC | TGT | GAA | GCG | TCT | GGA | TTC | ACC | CTC | AGA | AGT | TCT | GGC | ATG | CAC | TGG | | | | | | |
| | R | L | S | C | E | A | S | G | F | T | L | R | S | S | G | M | H | W | | | | | | |
| | | | | | | | | | | | | | CDR1 | | | | | | | | | | | |
| | 117 | | | | 126 | | | | 135 | | | | 144 | | | | 153 | | | | 162 | | | |
| | GTC | CGC | CAG | GCT | CCT | GGC | AAG | GGG | CTG | GAG | TGG | GTG | GCA | CTT | ATA | TGG | TTT | GAT | | | | | | |
| | V | R | Q | A | P | G | K | G | L | E | W | V | A | L | I | W | F | D | | | | | | |
| | | | | | | | | | | | | | CDR2 | | | | | | | | | | | |
| | 171 | | | | 180 | | | | 189 | | | | 198 | | | | 207 | | | | 216 | | | |
| | GGA | AGT | ATC | AGA | TCG | TAT | GCA | GAA | TCC | GTG | AAG | GGC | CGA | TTC | ACC | ATC | TCC | AGA | | | | | | |
| | G | S | I | R | S | Y | A | E | S | V | K | G | R | F | T | I | S | R | | | | | | |
| | CDR2 | | | | | | | | | | | | | | | | | | | | | | | |
| | 225 | | | | 234 | | | | 243 | | | | 252 | | | | 261 | | | | 270 | | | |
| | GAC | ACT | TCC | AAG | AAC | ACC | CTA | TAT | CTC | CAA | ATG | CGC | AGT | CTG | AGT | GCC | GAC | GAC | | | | | | |
| | D | T | S | K | N | T | L | Y | L | Q | M | R | S | L | S | A | D | D | | | | | | |
| | 279 | | | | 288 | | | | 297 | | | | 306 | | | | 315 | | | | 324 | | | |
| | ACG | GCT | GTG | TAT | TAC | TGT | GCG | AGA | GAC | AAG | GCG | GTT | CGG | GGA | ATT | AGC | AGG | TAC | | | | | | |
| | T | A | V | Y | Y | C | A | R | D | K | A | V | R | G | I | S | R | Y | | | | | | |
| | CDR3 | | | | | | | | | | | | | | | | | | | | | | | |
| | 333 | | | | 342 | | | | 351 | | | | 360 | | | | 369 | | | | | | | |
| | AAC | TAT | TAC | ATG | GAC | GTC | TGG | GGC | AAA | GGG | ACC | ACG | GTC | ACC | GTC | TCC | TCA | 3' | | | | | | |
| | N | Y | Y | M | D | V | W | G | K | G | T | T | V | T | V | S | S | | | | | | | |
| | CDR3 | | | | | | | | | | | | | | | | | | | | | | | |

Fig. 12b

LD2-11-VL sequence

```

5'  GTG TTG ACC CAG TCT CCA TCC TCC CTG TCT GCA TCT ATA CGA GAC AGA GTC ACC
    V  L  T  Q  S  P  S  S  L  S  A  S  I  R  D  R  V  T
-----
    63      72      81      90      99      108
ATC ACT TGC CGG GCA AGT CAG AAC ATT GGC AGT TAT TTA AAT TGG TAT CAG CAC
    I  T  C  R  A  S  Q  N  I  G  S  Y  L  N  W  Y  Q  H
-----
                CDR1 -----
    117      126      135      144      153      162
AAA CCA GGG ACA GCC CCT AAA CTC CTG ATC TAT GCT GTA TCC GCT TTG CAA AGT
    K  P  G  T  A  P  K  L  L  I  Y  A  V  S  A  L  Q  S
-----
                                CDR2 -----
    171      180      189      198      207      216
GGG GTC CCA TCG AGG TTC AGT GGC AGT AGA TCT GGG ACA GAT TTC ACT CTC ACC
    G  V  P  S  R  F  S  G  S  R  S  G  T  D  F  T  L  T
-----
    225      234      243      252      261      270
ATC AGC AGT CTG CAA CCT GAA GAT TTT GCA ACT TAC TAC TGT CAA CAG AGT TAC
    I  S  S  L  Q  P  E  D  F  A  T  Y  Y  C  Q  Q  S  Y
-----
    279      288      297      306      315
AGT CCC CCG TAC ACT TTC GGG CAG GGG ACC AAC CTG CAG ATC AAA 3'
    S  P  P  Y  T  F  G  Q  G  T  N  L  Q  I  K
-----
    CDR3 -----

```

Fig. 13a

LD2-14-VH sequence

| | | | | | | | | | | | | | | | | | | |
|----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|----------|-----|-----|-----|-----|-----|
| | | 9 | | 18 | | 27 | | 36 | | 45 | | 54 | | | | | | |
| 5' | CAG | GTG | AAA | CTG | CTC | GAG | TCT | GGG | GGA | GGC | GTG | GTC | CAG | CCG | GGG | GGG | TCC | CTG |
| | Q | V | K | L | L | E | S | G | G | G | V | V | Q | P | G | G | S | L |
| | | 63 | | 72 | | 81 | | 90 | | 99 | | 108 | | | | | | |
| | AGA | GTC | GCC | TGT | GTA | GCG | TCT | GGA | TTC | ACC | TCC | AGG | AAT | TTT | GGC | ATG | CAC | TGG |
| | R | V | A | C | V | A | S | G | F | T | S | R | N | F | G | M | H | W |
| | | | | | | | | | | | | | ← CDR1 → | | | | | |
| | GTC | CGC | CAG | GCT | CCA | GGA | AAG | GGG | CTG | GAG | TGG | GTG | GTT | TTT | ATT | TGG | TTT | GAT |
| | V | R | Q | A | P | G | K | G | L | E | W | V | V | F | I | W | F | D |
| | | | | | | | | | | | | | ← CDR2 → | | | | | |
| | GCA | AGT | AAT | AAA | GGA | TAT | GGA | GAC | TCC | GTT | AAG | GGC | CGA | TTC | ACC | GTC | TCC | AGA |
| | A | S | N | K | G | Y | G | D | S | V | K | G | R | F | T | V | S | R |
| | | | | | | | | | | | | | ← CDR2 → | | | | | |
| | GAC | AAT | TCC | AAG | AAC | ACG | CTC | TAT | CTG | CAA | ATG | AAC | GGC | CTG | AGA | GCC | GAA | GAC |
| | D | N | S | K | N | T | L | Y | L | Q | M | N | G | L | R | A | E | D |
| | | 225 | | 234 | | 243 | | 252 | | 261 | | 270 | | | | | | |
| | ACG | GCT | GTA | TAT | TAT | TGT | GCG | AGA | GAG | AAG | GCG | GTT | CGG | GGA | ATT | AGT | AGA | TAC |
| | T | A | V | Y | Y | C | A | R | E | K | A | V | R | G | I | S | R | Y |
| | | | | | | | | | | | | | ← CDR3 → | | | | | |
| | AAC | TAC | TAC | ATG | GAC | GTC | TGG | GGC | AAG | GGG | ACC | ACG | GTC | ACC | GTC | TCC | TCA | 3' |
| | N | Y | Y | M | D | V | W | G | K | G | T | T | V | T | V | S | S | |
| | | | | | | | | | | | | | ← CDR3 → | | | | | |

Fig. 13b

LD2-14-VL sequence

| | | | | | | | | | | | | | | | | | | | |
|----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|--|
| 5' | GTG | ATG | ACC | CAG | TCT | CCA | TCC | TCC | CTG | TCT | GCA | TCT | GTG | GGA | GAC | AGA | GTC | ACC | |
| | V | M | T | Q | S | P | S | S | L | S | A | S | V | G | D | R | V | T | |
| | | | | | | | | | | | | | | | | | | | |
| | ATC | ACT | TGC | CGG | GCA | AGT | CAG | AGC | ATT | ATC | AAC | AAT | TTA | AAT | TGG | TAT | CAG | CAG | |
| | I | T | C | R | A | S | Q | S | I | I | N | N | L | N | W | Y | Q | Q | |
| | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | |
| | AAA | CCA | GGC | AAA | GCC | CCT | GAA | CTC | CTG | ATC | TAT | GCT | GCA | TCC | AGT | TTG | CAA | AGT | |
| | K | P | G | K | A | P | E | L | L | I | Y | A | A | S | S | L | Q | S | |
| | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | |
| | GGG | GTC | CCT | TCA | AGG | TTC | CGT | GGC | AGT | GGA | TCT | GGG | AGA | GAT | TTC | ACT | CTC | ACC | |
| | G | V | P | S | R | F | R | G | S | G | S | G | R | D | F | T | L | T | |
| | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | |
| | GTC | ACC | AGT | CTG | CAA | CCT | GAA | GAT | TTT | GCA | ACT | TAC | TAC | TGT | CAA | CAG | AGT | TAC | |
| | V | T | S | L | Q | P | E | D | F | A | T | Y | Y | C | Q | Q | S | Y | |
| | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | |
| | AGT | AAC | CCT | GTG | GAC | GTT | CGG | CAA | GGG | ACC | AAG | GTG | GAA | ATC | AAA | | | | |
| | S | N | P | V | D | V | R | Q | G | T | K | V | E | I | K | | | | |
| | | | | | | | | | | | | | | | | | | | |

CDR3 →

Fig. 14a

LD2-17-VH sequence

```

5'  CAG GTG AAA CTG CTC GAG TCT GGG GGA GGC GTG GTC CAG CCG GGG GGG TCC CTG
    Q  V  K  L  L  E  S  G  G  G  V  V  Q  P  G  G  S  L
-----
    AGA CTC TCC TGT GTA GCG TCT GGA TTC ACC TCC AGG AGT TAT GGA ATG CAC TGG
    R  L  S  C  V  A  S  G  F  T  S  R  S  Y  G  M  H  W
-----
                                CDR1
GTC  CGC  CAG  GCT  CCA  GGC  AAG  GGC  CTG  GAG  TGG  GTG  GCT  TTT  ATA  TGG  TTT  GAT
    V  R  Q  A  P  G  K  G  L  E  W  V  A  F  I  W  F  D
-----
                                CDR2
GGA  AGT  AAT  AAA  GGA  TAT  GTA  GAC  TCC  GTG  AAG  GGC  CGA  TTC  ACC  ATC  TCC  CGA
    G  S  N  K  G  Y  V  D  S  V  K  G  R  F  T  I  S  R
-----
                                CDR2
GAC  AAT  TCC  AAG  AAC  ACG  CTC  TAT  CTG  CAA  ATG  AAG  AGC  CTG  AGA  GCC  GAG  GAC
    D  N  S  K  N  T  L  Y  L  Q  M  K  S  L  R  A  E  D
-----
ACG  GCT  GTA  TAT  TAT  TGT  GCG  AGA  GAG  AAG  GCG  CTT  CGG  GGA  ATC  AGT  AGA  TAC
    T  A  V  Y  Y  C  A  R  E  K  A  L  R  G  I  S  R  Y
-----
                                CDR3
AAC  TAT  TAC  CTG  GAC  GTC  TGG  GGC  AAG  GGG  ACC  ACG  GTC  ACC  GTC  TCC  TCA  3'
    N  Y  Y  L  D  V  W  G  K  G  T  T  V  T  V  S  S
-----
                                CDR3

```

Fig. 14b

LD2-17-VL sequence

```

5'  GTG  ATG  ACC  CAG  TCT  CCA  TTC  TCC  CTG  TCT  GCA  TCT  GTA  GGA  GAC  AGA  GTC  ACC
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    V   M   T   Q   S   P   F   S   L   S   A   S   V   G   D   R   V   T

      63      72      81      90      99      108
ATC  ACT  TGC  CGG  GCA  AGT  CAG  AAC  ATT  AGG  AGT  TTT  TTA  AGT  TGG  TAT  CAG  CAG
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    I   T   C   R   A   S   Q   N   I   R   S   F   L   S   W   Y   Q   Q

      117      126      135      144      153      162
AAA  CCA  GGG  ACA  GCC  CCT  AAG  CTC  CTG  ATC  TAT  GCT  GCA  TCC  AGG  TTG  CAA  AGT
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    K   P   G   T   A   P   K   L   L   I   Y   A   A   S   R   L   Q   S

      171      180      189      198      207      216
GGG  GTC  CCA  TCA  AGG  TTC  AGT  GGC  AGT  GGG  TCT  GGG  ACA  GAT  TTC  ACT  CTC  ACC
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    G   V   P   S   R   F   S   G   S   G   S   G   T   D   F   T   L   T

      225      234      243      252      261      270
ATC  AGC  ACT  CTG  CAA  CCT  GAA  GAT  TTT  GCG  ACT  TAC  TAC  TGT  CAA  CAG  AGT  TAC
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    I   S   T   L   Q   P   E   D   F   A   T   Y   Y   C   Q   Q   S   Y

      279      288      297      306      315
AGT  GCC  CCT  TGG  ACG  TTC  GGC  CAA  GGG  ACC  AAG  CTG  GAA  ATC  AAA  3'
    ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---  ---
    S   A   P   W   T   F   G   Q   G   T   K   L   E   I   K

      CDR1  CDR2  CDR3

```

Fig. 15a

LD2-18-VH sequence

```

5'      9      18      27      36      45      54
CAG GTG AAA CTG CTC GAG TCT GGG GGA GGC TTG GTC CAG CCG GGG GGG TCC CTG
-----
Q   V   K   L   L   E   S   G   G   G   L   V   Q   P   G   G   S   L

      63      72      81      90      99      108
AGA CTC TCC TGT GTA GCG TCT GGA TTC ACC TTC AGG AGT TAT GGC ATG CAC TGG
-----
R   L   S   C   V   A   S   G   F   T   F   R   S   Y   G   M   H   W

                                  CDR1
GTC CGC CAG GCT CCA GGC AAG GGC CTG GAG TGG GTG GCT TTT ATA TGG TTT GAT
-----
V   R   Q   A   P   G   K   G   L   E   W   V   A   F   I   W   F   D

                                  CDR2
GGA AGT AAT AAA GGA TAT GTA GAC TCC GTG AAG GGC CGA TTC ACC ATC TCC CGA
-----
G   S   N   K   G   Y   V   D   S   V   K   G   R   F   T   I   S   R

      CDR2
GAC AAT TCC AAG AAC ATG CTC TAT CTG CAA ATG AAT AGC CTG AGA GCC GAG GAC
-----
D   N   S   K   N   M   L   Y   L   Q   M   N   S   L   R   A   E   D

      279      288      297      306      315      324
ACG GCT GTA TAT TAT TGT GCG AGA GAG AAG GCG CTT CGG GGA ATC AGT AGA TAC
-----
T   A   V   Y   Y   C   A   R   E   K   A   L   R   G   I   S   R   Y

                                  CDR3
AAC TAT TAC CTG GAC GTC TGG GGC AAG GGG ACC ACG GTA ACC GTC TCC TCA 3'
-----
N   Y   Y   L   D   V   W   G   K   G   T   T   V   T   V   S   S

      CDR3

```

Fig. 15b

LD2-18-VL sequence

```

5'  GTG  ATG  ACC  CAG  TCT  CCA  TCC  TCC  CTG  TCT  GTA  TCT  ATA  GGG  GAA  AGA  GTC  ACC
      9      18      27      36      45      54
      V   M   T   Q   S   P   S   S   L   S   V   S   I   G   E   R   V   T
-----
      63      72      81      90      99      108
ATC  ACT  TGC  CGG  GAA  AGT  CAG  AGC  GTT  ACC  AGG  TCT  TTA  ATT  TGG  TTT  CAG  AAG
      I   T   C   R   E   S   Q   S   V   T   R   S   L   I   W   F   Q   K
-----
      117      126      135      144      153      162
AAA  CCA  GGG  AAA  GCC  CCT  AGG  CTC  CTA  ATC  TTT  GTT  GCG  TCC  ACT  TGG  AAA  AGT
      K   P   G   K   A   P   R   L   L   I   F   V   A   S   T   W   K   S
-----
      171      180      189      198      207      216
GGG  GTC  CCA  TCA  AGG  TTC  AGT  GGC  AGT  GGA  TCT  GGG  ACA  GAT  TTC  ACC  CTC  ACC
      G   V   P   S   R   F   S   G   S   G   S   G   T   D   F   T   L   T
-----
      225      234      243      252      261      270
ATC  AGC  AGT  CTG  CAA  CCT  GAG  GAT  TTT  GGA  ACT  TAC  TAC  TGT  CAA  CAG  AAT  TAC
      I   S   S   L   Q   P   E   D   F   G   T   Y   Y   C   Q   Q   N   Y
-----
      279      288      297      306      315      324
AGG  ACC  CCT  CAG  TGG  ACG  TTC  GGC  CAA  GGG  ACC  AAG  GTA  GAA  ATC  AAA  3'
      R   T   P   Q   W   T   F   G   Q   G   T   K   V   E   I   K
-----
      CDR3 ----->
      CDR1 ----->
      CDR2 ----->

```

Fig. 16a

LD2-20-VH sequence

```

5'  CAG GTG AAA CTG CTC GAG TCT GGG GGA GGC GTG GTC CAG CCG GGG GGG TCC CTG
    Q  V  K  L  L  E  S  G  G  G  V  V  Q  P  G  G  S  L

    63      72      81      90      99      108
AGA CTC TCC TGT GTA GCG TCT GGA TTC ACC TCC AGG AGT TAT GGC ATG CAC TGG
R  L  S  C  V  A  S  G  F  T  S  R  S  Y  G  M  H  W

    117      126      135      144      153      162
GTC CGC CAG GCT CCA GGA AAG GGC CTG GAG TGG GTG GCT TTT ATT TGG TTT GAT
V  R  Q  A  P  G  K  G  L  E  W  V  A  F  I  W  F  D

    171      180      189      198      207      216
GGA AGT AAT AAA GGA TAT GTA GAC TCC GTG AAG GGC CGA TTC ACC ATC TCC CGA
G  S  N  K  G  Y  V  D  S  V  K  G  R  F  T  I  S  R

    225      234      243      252      261      270
GAC AAT TCC AAG AAC ACG CTC TAT CTG CAA ATG AAG AGC CTG AGA GCC GAG GAC
D  N  S  K  N  T  L  Y  L  Q  M  K  S  L  R  A  E  D

    279      288      297      306      315      324
ACG GCT GTA TAT TAT TGT GCG AGA GAG AAG GCG CTT CGG GGA ATC AGT AGA TAC
T  A  V  Y  Y  C  A  R  E  K  A  L  R  G  I  S  R  Y

    333      342      351      360      369
AAC TAT TAC CTG GAC GTC TGG GGC AAG GGG ACC ACG GTC ACC GTC TCC TCA 3'
N  Y  Y  L  D  V  W  G  K  G  T  T  V  T  V  S  S

    CDR1  ←-----→
    CDR2  ←-----→
    CDR3  ←-----→
  
```

Fig. 16b

LD2-20-VL sequence

```

5'  GTG  ATG  ACC  CAG  TCT  CCA  TCC  TCC  CTG  TCT  GCA  TCT  GTA  GGA  GAC  AGA  GTC  ACC
      9      18      27      36      45      54
      V  M  T  Q  S  P  S  S  L  S  A  S  V  G  D  R  V  T
-----
      63      72      81      90      99      108
ATC  ACT  TGC  CGG  GCA  AGT  CAG  AGC  ATT  AGC  AGC  TAT  TTA  AAT  TGG  TAT  CAG  CAG
-----
      I  T  C  R  A  S  Q  S  I  S  S  Y  L  N  W  Y  Q  Q
      <----- CDR1 ----->
      117      126      135      144      153      162
AAA  CCA  GGG  AAA  GCC  CCT  AAG  CTC  CTG  ATC  TAT  GCT  GCA  TCC  AGT  TTG  CAA  AGT
-----
      K  P  G  K  A  P  K  L  L  I  Y  A  A  S  S  L  Q  S
      <----- CDR2 ----->
      171      180      189      198      207      216
GGG  GTC  CCA  TCA  AGG  TTC  AGT  GGC  AGT  GGA  TCT  GGG  ACA  GAT  TTC  ACT  CTC  ACC
-----
      G  V  P  S  R  F  S  G  S  G  S  G  T  D  F  T  L  T
      225      234      243      252      261      270
ATC  AGC  AGT  CTG  CAA  CCT  GAA  GAT  TTT  GCA  ACT  TAC  TAC  TGT  CAA  CAG  AGT  TAC
-----
      I  S  S  L  Q  P  E  D  F  A  T  Y  Y  C  Q  Q  S  Y
      <----- CDR3 ----->
      279      288      297      306      315
AGT  ACC  CGA  TTC  ACT  TTC  GGC  CCT  GGG  ACC  AAA  GTG  GAT  ATC  AAA  3'
-----
      S  T  R  F  T  F  G  P  G  T  K  V  D  I  K

```

Fig. 17a

LD1-6-17-VH sequence

```

5'  CAG GTG AAA CTG CTC GAG TCT GGG GGA GGC GTG GTC CAG CCT GGG AGG TCC CTG
    Q  V  K  L  L  E  S  G  G  G  V  V  Q  P  G  R  S  L
-----
    AGA CTT TCC TGT GCA GCG TCT GGA TTT ACC TTC AGT AGC TAT GGA ATG CAC TGG
    R  L  S  C  A  A  S  G  F  T  F  S  S  Y  G  M  H  W
-----
                                CDR1
    GTC CGC CAG GCT CCA GGA AAG GGG CTG GAG TGG GTG ACA GAT ATA TGG TTT GAT
    V  R  Q  A  P  G  K  G  L  E  W  V  T  D  I  W  F  D
-----
                                CDR1
    GGA GGT AAT AAA CAT TAT GCA GAC TTC GTG AAG GGC CGA TTC ACC ATC TCC AGA
    G  G  N  K  H  Y  A  D  F  V  K  G  R  F  T  I  S  R
-----
                                CDR2
    GAC AAT TCC AAG AAC ACG GGG TTT CTA CAA ATG AAC AGC CTG AGA GTC GAG GAC
    D  N  S  K  N  T  G  F  L  Q  M  N  S  L  R  V  E  D
-----
    ACG GCT GTG TAT TAC TGT GCG AGG GAT TAC TAT AGC GTT ACT AAG AAA CTC AGA
    T  A  V  Y  Y  C  A  R  D  Y  Y  S  V  T  K  K  L  R
-----
                                CDR3
    CTC CAC TAC TAC TAC TAC ATG GAC GTC TGG GGC AAA GGG ACC ACG GTC ACC GTC
    L  H  Y  Y  Y  Y  M  D  V  W  G  K  G  T  T  V  T  V
-----
                                CDR3
TCC TCA 3'
S  S

```


Fig. 17b

LD1-6-17-VL sequence





Fig. 18

The pComb3 Expression System






Fig. 19




Fig. 20

